

GSA prices effective April 2, 2018 (based on February 6, 2017, commercial price list) Published August 2018

Laboratory Furniture

Co/Struc® System

General Services Administration Federal Supply Service Authorized Federal Supply Schedule Price List

Online access to contract ordering information, terms and conditions, up-to-date pricing, and the option to create an electronic delivery order is available through GSA Advantage!*, a menu-driven database system. The Internet address for GSA Advantage! is www.gsaadvantage.gov. For more information on ordering from Federal Supply Schedules, click on the FSS Schedules button at www.fss.gsa.gov.

FSC Group 66, Part II, Section N FSC Classes 6640 Contract Number: GS-07F-5623R

Contract Period: February 1, 2015, through January 31, 2020

Herman Miller, Inc. 855 East Main Avenue Zeeland, Michigan 49464 (616) 654 3000 Phone (616) 654 8278 Fax HermanMiller.com/government

Business Size: Large

1a. Special Item Numbers (SIN)

- 566-1: Modular Laboratory Furniture Systems, all product except non-modular Laboratory Tables, see 566-2 below.
- 566-2:Individual non-modular Laboratory Tables (C4115. and C4215.).
- 566-5: Workstation Design/Installation option. Including:

Workstation Design/Layout Option: Includes the design/redesign of lab furniture. Design service will be provided by Herman Miller, Inc., or its participating dealer based on a fee of \$65 per hour.

Workstation Furniture and Equipment Installation: Includes delivery coordination of furniture and related equipment, unpacking, inspection, inventory and assembly of furniture components. Pricing shall be negotiated on a project-by-project basis.

Workstation Disassembly and Reinstallation Option: Includes the taking down of installed furniture and related hardware for the purpose of packing and moving or storage. Pricing shall be negotiated on a project-by-project basis.

Workstation Repair: Includes the repair and/or replacement of damaged furniture and related hardware. Pricing shall be negotiated on a project-by-project basis.

1b. Lowest Price Model:

566-1 = C1315.06 \$19.00 566-2 = C4700.2436L \$229.00

- 1c. Not applicable.
- 2. Maximum Order. \$400,000 net product value.
- 3. Minimum Order. \$50 net.
- Geographic Coverage. 48 contiguous states and District of Columbia. Other geographic locations see box below.
- 5. Production Points. Zeeland, Michigan, County of Ottawa.
- Pricing. Prices shown are list. To calculate contract prices, deduct 58.5%.
- **7. Quantity Discounts.** Reference item #6.
- 8. Payment Terms. Net 30 days.
- 9a. Government Purchase Cards. ARE accepted below the micropurchase threshold.
- 9b. Government Purchase Cards. ARE accepted above the micropurchase threshold.
- 10. Foreign Items. Not applicable.
- 11a. Time of Delivery. 56 days ARO.
- 11b. Expedited Delivery. Select items may be available for expedited delivery. Contact Herman Miller, Inc., or local authorized dealer for availability.
- 11c. Overnight and 2-Day Delivery. Select items may be available for overnight or 2-day delivery. Contact Herman Miller, Inc., or local authorized dealer for availability.

- 11d. Urgent Requirements. Select items may be available for urgent delivery. Contact Herman Miller, Inc., or local authorized dealer for availability.
- **12. F.O.B. Point.** Destination within the 48 contiguous states and Dist. of Columbia.
- **13a. Order Address.** Herman Miller, Inc., Government Customer Care o161, 855 East Main Avenue, Zeeland, MI 49464.
- **13b. Order Procedures.** For supplies and services, the ordering procedures, information on Blanket Purchase Agreements (BPA's), and a sample BPA can be found at the GSA/FSS Schedule home page (www.fss.gsa.gov/schedules).
- Payment Address. Herman Miller, Inc., 22764 Network Place, Chicago, IL 60673-1227.
- Warranty Provisions. Herman Miller, Inc., commercial warranty applies.
- 16. Export Packing Charges. Prices supplied on request.
- Terms and Conditions of Government Commercial Credit Card Acceptance. None.
- Terms and Conditions of Rental, Maintenance, and Repair. Not applicable.
- 19. Terms and Conditions of Installation. Installation services shall be negotiated between ordering agency and Herman Miller local authorized dealer.
- 20. Terms and Conditions of Repair Parts. Not applicable.
- 20a. Terms and Conditions for Any Other Services. Services shall be negotiated between ordering agency and Herman Miller local authorized dealer.
- List of Service and Distribution Points. Access Herman Miller web site for current list of participating dealers: www.HermanMiller.com.
- **22. List of Participating Dealers.** Access Herman Miller web site for current list of participating dealers: www.HermanMiller.com.
- 23. Preventive Maintenance. Not applicable.
- 24. Special Attributes. Reconfigurable, renewable.
- **25. DUNS Number.** 006012801
- 26. Registered in Central Contractor Registration (CCR) Database. "Herman Miller, Inc."

Herman Miller Miscellaneous Information:

Federal ID Number: 38-0837640

Commercial and Government Entity Code (CAGE): 40636

Co/Struc® System Structure Work Surfaces Storage Transport Indices By Name By Number Appendices	page 2
Co/Struc® System	3
Structure	5
Work Surfaces	7.6
Storage	110
	222
Indices	229
By Name	229
By Number	231
Appendices	
20-Day or Less/Assigned Lead-Time Order Inform	ation - Healthcare
Packaging/Keyed-Alike Information-Co/Struc Sys	tem
For the latest materials information, refer to the Ma	aterials pages on HermanMiller.com.

A Note on the Organization of This Book

There are a few signals that will help you understand the organization of this book. Once you know them, you should be able to find your way around easily.

Like a newspaper, this book is formatted with columns of text and illustrations that run vertically.

Information wraps from one column to the next and continues for as many pages as it needs to.

Black bars are clues.

A black bar at the top of a page signals the start of information about a product.

This information is divided into two sections, each signaled by a black bar. Product Information includes a description, additional notes, and dimensional drawings. Specification Information, which always begins at the top of a column, includes a series of steps. Each step represents a decision you need to make in order to complete your product number. Steps must be followed sequentially. The number of pages needed to complete information on a product varies; some will be complete in one page; others may fill three or four pages. Just continue going through the steps until you run into the next product, signaled by the black bar at the top of the page.

Black tabs running along the outside edge of the page help you locate what section you are in or what section you want to go to. You can thumb forwards or backwards to find a section; tabs and text appear on both sides of the page.

The additional pages in the back of the book include the appendices and indices. They are signaled by a long black bar running along the outside edge of the page.

The index is presented two ways.

The first index is in alphabetical order by product name. It can be helpful if you are new to the products and numbering scheme or if you are looking for a type of product, such as work chair or round table. The second index is by product number. It lists, in alphanumeric order, the base product number, typically a five-digit number.

General Information

This book is effective April 2, 2018, subject to change without notice. Products may be purchased from authorized Herman Miller dealers who will quote prices upon request. For more information about our products and services or to see a list of dealers, visit www. HermanMiller.com.

All prices are list prices. Additional services, such as planning services, design, storage, and installation, are not included and must be added to these prices when the additional services are requested.

Illustrations, specifications, and prices are based on the latest information at the time of publication. The right is reserved to make changes at any time without notice in prices, fabrics, finishes, materials, specifications, and models, and also to discontinue models and finishes.

Dimensions shown in the drawings are engineering drawing dimensions rounded up to the nearest $^{1}/_{8}$ ". Product is shipped set up or knocked down. Products that include "Shipped knocked down" in the description require some assembly. If this statement does not appear in the description, the product is shipped assembled.

Service problems are normally handled by a Herman Miller authorized dealer. Where this cannot be accomplished on a local level, service problems should be referred to the Customer Care Representative for your region, Herman Miller, Inc., Zeeland, Michigan 49464.


20-Day or Less/Assigned Lead-Time Shipments

All products and options in this price book not designated by an Assigned Program icon <a> a will ship in **20 business days or less** after being acknowledged by Herman Miller.

Products and options designated by an A are on the assigned leadtime program and may ship in 20 days, less than 20 days, or more than 20 days from order acknowledgement.

Structure Work Surfaces Storage

Transport



Description


This assembly attaches to a wall and supports Co/Struc hanging components. It includes a 12' rail, a 12' cover strip, a 12' back plate, 25 stand-offs, and 5 spacers.

Notes


To attach assembly to wall, order rail hardware (CO28441, or CO28541) separately.

To finish end, order end cap for standard or seismic rail (CO₃₃₈) separately.

When seismic components are required, order seismic wall rail (CO_{479}) instead of standard rail (CO_{282}).



Spec	cification Information				
Step	1.				
CO28	CO282FF				
Step	2. Surface Finish				
HF	inner tone light	+\$o			
LT	light tone	+ \$o			
LU	soft white	+\$ o			
WN	warm grey neutral	+\$0			



Description

This package contains enough fasteners to install 12' of standard rail. Package contains 300 fastener sets.

The 3" sheet metal screws with plastic anchors are #14, flat-head, Phillips-drive screws with $1^3/8$ "-long plastic anchors. The 4" toggle bolts with wall bushings are 1/4-20, flat-head, Phillips-drive toggle bolts.

Specification Information

Step 1.


CO28

Step 2. Hardware

441 3" sheet metal screws with plastic anchors

541 4" toggle bolts with wall bushings


Prices for Steps 1-2.	
CO28441	\$222
CO28541	\$650



Description

This rail attaches to a panel, a module, or wall strips to support Co/Struc hanging components or a storage case in a straight-line configuration. It does not hang directly on a wall. 2 rails are required to hang a locker from an open panel frame or glazed panel.

Dimensions



Specification Information


Step 1.

C1610.

Step 2	Step 2. Width					
24	24" wide					
30	30" wide					
36	36" wide					
42	42" wide					
48	48" wide					
60	6o" wide					

Prices for Steps 1-2.	
C1610. 24	\$184
30	\$191
36	\$199
42	\$219
48	\$233
60	\$250

Step 3	. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0



Description

This aluminum rail attaches to an architectural wall. It holds work tools or accessories and can be field cut to the appropriate width. The rail includes plastic end caps to accommodate 2 field-cut rail sections.

Notes

Customer must supply attachment hardware for wall conditions.

Order rail-mounted accessories separately:


- Chart holder (CT730.)
- Chest tube rack (CT737.)
- Glove box holder (CT726.)
- Sharps container holder (CT734.)
- Supply basket (CT720. or CT722.)
- Suture basket (CT723.)
- Tilt-out storage bin (CT9303., CT9304., CT9305., CT9306.)
- Utility hook (CT776.)
- Wastebasket and holder (CT728.)

Order work tools separately:

- Arch (Y7150.)
- Diagonal tray (Y7218.)
- Mini-shelf (Y7310.)
- Organizer tray (Y7210.)
- Paper tray (Y7216.)
- Rail divider (Y7320.)
- Vertical tray (Y7217.)

To mount customer-supplied accessories, order rail clamp (CT729.) separately.

Dimensions



Specification Information


Step 1.

C1170.

Step :	2. Width			
24	24" wide			
30	30" wide			
36	36" wide			
42	42" wide			
48	48" wide			
60	6o" wide			

Prices for Steps 1-2.	
C1170. 24	\$91
30	\$110
36	\$126
42	\$136
48	\$148
60	\$169

Step 3.	End Cap Finish	
LU	soft white	+\$0
WN	warm grev neutral	+\$c



Description

This aluminum rail attaches to a panel, freestanding module, or wall strips and holds work tools or accessories. It includes steel end caps. Attachment hardware is included.

Notes

Order rail-mounted accessories separately:


- Chart holder (CT730.)
- Chest tube rack (CT737.)
- Glove box holder (CT726.)
- Sharps container holder (CT734.)
- Supply basket (CT720. or CT722.)
- Suture basket (CT723.)
- Tilt-out storage bin (CT9303., CT9304., CT9305., CT9306.)
- Utility hook (CT776.)
- Wastebasket and holder (CT728.)

Order work tools separately:

- Arch (Y7150.)
- Diagonal tray (Y7218.)
- Mini-shelf (Y7310.)
- Organizer tray (Y7210.)
- Paper tray (Y7216.)
- Rail divider (Y7320.)
- Vertical tray (Y7217.)

To mount customer-supplied accessories, order rail clamp (CT729.) separately.

Dimensions



Specification Information


Step 1.

C1172.

Step :	Step 2. Width							
24	24" wide							
30	30" wide							
36	36" wide							
48	48" wide							

Prices for Steps 1-2.	
C1172. 24	\$129
30	\$136
36	\$148
48	\$185

Step 3.	End Cap Finish	
LU	soft white	+\$o
WN	warm grey neutral	+ \$o




Description

This end cap provides a finished end for a standard or seismic wall rail and a positive stop for a locker.


Notes

For areas requiring seismic components, order 1 end cap for each end of rail, and order seismic wall rail (CO479) separately.

For areas not requiring seismic components, order 1 end cap for each end of rail, and order standard rail (CO282) separately.



Spec	ification Information	
Step :	1.	
CO33	CO338FF	
Step :	2. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+ \$o
LU	soft white	+ \$o
WN	warm grey neutral	+ \$o



Description

This aluminum rail attaches to a wall and supports Co/Struc hanging components. It includes a 12' rail, a 12' cover strip, a 12' back plate, 26 stand-offs, and 5 spacers.

Notes


Order rail for areas requiring seismic components, when attaching to studs 16" on center, or when attaching to plaster in accordance with OPA-0117 in California.

To meet the conditions of appropriate OPA, refer to the OPA document for specific hardware. For other applications, order rail hardware (CO28441, or CO28541) separately.


Order additional seismic components separately:

• End cap for standard or seismic rail (CO338)

Field drilling is required for proper installation.



Spec	ification Information		
Step 1			
CO479	CO479FF		
Step 2	2. Surface Finish		
HF	inner tone light	+\$0	
LT	light tone	+\$o +\$o	
LU	soft white	+\$o +\$o	
WN	warm grey neutral	+\$o	




Description

This bracket secures Co/Struc rail-hanging components to an architectural wall. It is designed to be used in areas requiring seismic components. Finish is soft white. Package contains 10.

Notes

Customer must supply required fasteners.


Dimensions



Specification Information

Step 1.

CO569.



Description

This draw rod connects Action Office® Series 1 or 2 panels or Co/Struc® support panels in a straight line. It has a threaded hole in the draw block for attachment to a seismic floor anchor bracket.

Notes

For Action Office Series 1 panels or Co/Struc support panels, specify height of draw rod to match height of panel.

For Action Office Series 2 panels, specify draw rod as follows:

- 32"-high panel, CO215.27
- 39"-high panel, CO215.34
- 47"-high panel, CO215.42
- 53"-high panel, CO215.48
- 62"-high panel, CO215.57
- 67"-high panel, CO215.62
- 85"-high panel, CO215.80

Order floor anchor bracket, Action Office Series 1 (CO481.) or floor anchor bracket, Action Office Series 2 (CO485.) separately.

Dimensions




Specification Information

Step 1.

CO215.

Step 2	. Height			
27	27" high			
34	34" high			
42	42" high			
48	48" high			
57	57" high			
62	62" high			
80	8o" high			

Prices for Steps 1-2.	
CO215. 27	\$25
34	\$26
42	\$26
48	\$27
57	\$27
62	\$27
80	\$31



Description


This bracket fastens both sides of a module frame to the floor. Each anchor location requires 1 bracket. Finish is light tone. Package contains 10. Hardware to attach module is included. Floor anchor and bolt must be customer supplied.

Notes

Order bracket for areas requiring seismic components.

Lab modules need to be mounted flush for anchor to install.


Dimensions



Specification Information

Step 1.

CO568FF



Description

This bracket fastens both sides of an Action Office® Series 1 panel to the floor and is designed for use in areas requiring seismic components. Finish is light tone. Package contains 10.

Notes

Order the following products separately:


- Draw rod (CO215.)
- Floor anchor adapter kit (CO552.)
- Seismic finished end (A1250.) when using only Action Office Series 2 components
- Seismic hingeable end (AO258.) when using Action Office Series 1 or Co/Struc® components


Customer must supply required bolts and concrete anchors.

Dimensions










Specification Information

Step 1.

CO481



Description

This bracket fastens both sides of an Action Office® Series 2 panel to the floor and is designed for use in areas requiring seismic components. Package contains 10.


Notes

Order the following products separately:

- Draw rod (CO215.)
- Floor anchor adapter kit (CO552.)
- Seismic finished end (A1250.)

Customer must supply required concrete anchors.


Dimensions



Specification Information

Step 1.

CO485




Description

This kit replaces the bottom block of an Action Office® Series 1 or 2 spacer or 2-way, 3-way, or 4-way connector. It is designed to be used in areas requiring seismic components and has a threaded hole in the draw block for attachment to a seismic floor anchor bracket. Package contains 10.

Notes

Order floor anchor bracket, Action Office Series 1 (CO481.) or floor anchor bracket, Action Office Series 2 (CO485.) separately.


Dimensions



Specification Information


Step 1.

CO552




Description

This top cap fits over a standard top cap and spans multiple Action Office® Series 1 or 2 panels. It provides extra rigidity along the panel run. The 12'-long top cap can be field cut to the appropriate length. Attachment hardware is included.



Spe	ecification Information	
Step	1.	
CO55	53FF	\$520
Step	2. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+ \$o
LU	soft white	+ \$o
MT	medium tone	+ \$o
WN	warm grey neutral	+\$0


Frame Stiffener, Ethospace® Top CO554 Cap




Product Information

Description

This top cap replaces a standard top cap and spans multiple Ethospace® frames. It provides extra rigidity along the frame run. The 12'-long top cap can be field cut to the appropriate length. Attachment hardware is included.



Spec	ification Information	
Step 1	ı.	
CO55	4FF	\$1024
	2. Surface Finish	
BU	black umber	+\$0
HF	inner tone light	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0




Description


This plastic cover snaps onto a wall strip to cover the slots. It helps with infection control in healthcare environments by preventing bacteria and germs from growing in these slots. It is 84" and must be field cut to match the wall strip height. Package contains 6 covers.

Notes

Order wall strips (AO213. or E1130.) separately.



Spe	cification Information	
Step	1.	
CO48	87.	\$439
Step	2. Finish	
BU	black umber	+\$ o
HF	inner tone light	+ \$o
LU	soft white	+ \$o
ΜT	medium tone	+\$o



Description

This open panel has hanger frames at 1" intervals to support components. It can be used as a return panel to stabilize a panel run. The panel cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.


Notes

Order draw rods (AO215.) or seismic draw rods (CO215.) separately. To attach panels of unequal height, order change of height draw rods (AO219.) and change of height panel/panel finished ends (AO272.) separately.

To join 2 support panels of equal or unequal heights at a 90° angle, order 2-way 90° corner connector (C1221.) separately.

To join 3 support panels of equal or unequal heights at a 90° angle, order 3-way 90° corner connector (C1231.) separately.

Dimensions




FRONT

Specification Information

Co/Struc® Structure

Step 1.	_					
C1150.						
Stona	Height					
48						
	48" high					
80	8o" high					
Step 3.	Width					
12	12" wide					
24	24" wide					
30	30" wide					
36	36" wide					
48	48" wide					
Prices 1	for Steps 1-3.					
		12	24	30	36	48
C1150.	48	\$455	473	507	546	591
	80	\$628	667	702	746	813
Step 4.	Surface Finish					
HF	inner tone light					+\$0
LT	light tone					+\$o
LU	soft white					+\$0
WN	warm grey neutral					+\$0
	5. c, neatrat					. 40



Description


This panel has a 48"-high recessed filler and hanger frames at 1" intervals to support components. It can be used as a return panel to stabilize a panel run. The panel has an area between the filler and the rear edge of an attached work surface for passage of wires or tubes. It cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

Notes

Order draw rods (AO215.) or seismic draw rods (CO215.) separately. To attach panels of unequal height, order change of height draw rods (AO219.) and change of height panel/panel finished ends (AO272.) separately.


To join 2 support panels of equal or unequal heights at a 90° angle, order 2-way 90° corner connector (C1221.) separately.

To join 3 support panels of equal or unequal heights at a 90° angle, order 3-way 90° corner connector (C1231.) separately.



Specification Information

Step:	1.					
C115	1.					
Step:	2. Height					
48	48" high					
80	8o" high					
Step	3. Width					
12	12" wide					
24	24" wide					
30	30" wide					
36	36" wide					
48	48" wide					
Prices	s for Steps 1-3.					
		12	24	30	36	48
C115	1. 48	\$606	751	833	904	1017
	80	\$811	915	1010	1135	1254
Step	4. Surface Finish					
HF	inner tone light					+\$o
LT	light tone					+\$o
LU	soft white					+\$o
WN	warm grey neutral					+\$o
	· ,					



Description

This panel has a vertical raceway for routing and concealing electrical wires and communication cables between the floor and a wire chase rail. It has a 48"-high recessed filler, hanger frames at 1" intervals to support components, and 6 knockouts: 3 at the floor level and 3 at 47" above floor level. It cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

Notes


To supply power from both sides of panel, order wire chase rail (C1313. or C1314.) separately and mount 48" above floor.

For power feed to wire chase rail mounted at heights other than 48", field cut removable cover on 1 side of vertical raceway.


Order draw rods (AO215.) or seismic draw rods (CO215.) separately. To attach panels of unequal height, order change of height draw rods (AO219.) and change of height panel/panel finished ends (AO272.) separately.

To join 2 support panels of equal or unequal heights at a 90° angle, order 2-way 90° corner connector (C1221.) separately.

To join 3 support panels of equal or unequal heights at a 90° angle, order 3-way 90° corner connector (C1231.) separately.



Spec	cification Information			
Step				
C115	5.			
Step:	2. Height			
48	48" high			
80	8o" high			
Step	3. Width			
24	24" wide			
30	30" wide			
48	48" wide			
Prices	s for Steps 1-3.			
		24	30	48
C115	5. 48	\$1103	1220	145
	80	\$1300	1378	162
Step	4. Surface Finish			
HF	inner tone light			+\$0
LT	light tone			+\$0
LU	soft white			+\$
WN	warm grey neutral			+\$0



Description

This panel has a vertical raceway for routing and concealing electrical wires and communication cables between the ceiling and a wire chase rail. It has a 48"-high recessed filler, hanger frames at 1" intervals to support components, and 6 knockouts: 3 at the floor level and 3 at 47" above floor level. It cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

Notes

To carry different types of electrical power (e.g. regular and emergency power) or to carry electrical power and low-voltage cables, vertical raceway can be divided into separate internal raceways; contact local contractor.


To supply power from both sides of panel, order wire chase rail (C1313. or C1314.) separately and mount 48" above floor.

For power feed to wire chase rail mounted at heights other than 48", field cut removable cover on 1 side of vertical raceway.

Order draw rods (AO215.) or seismic draw rods (CO215.) separately. To attach panels of unequal height, order change of height draw rods (AO219.) and change of height panel/panel finished ends (AO272.) separately.


To join 2 support panels of equal or unequal heights at a 90° angle, order 2-way 90° corner connector (C1221.) separately.

To join 3 support panels of equal or unequal heights at a 90° angle, order 3-way 90° corner connector (C1231.) separately.



Co/Struc® Structure

Spec	ification Information			
Step 1				
C1156				
Step 2	. Height			
48	48" high			
80	8o" high			
Step 3	g. Width			
24	24" wide			
30	30" wide			
48	48" wide			
	•			
Prices	for Steps 1-3.			
Prices		24	30	48
Prices	for Steps 1-3.	24 \$2674		48 2911
	for Steps 1-3.		2895	
	for Steps 1-3.	\$2674	2895	2911
C1156	for Steps 1-3.	\$2674	2895	2911 3166
C1156	for Steps 1-3. 5. 48 80	\$2674	2895	2911
C1156 Step 4	for Steps 1-3. 5. 48 80 6. Surface Finish	\$2674	2895	2911 3166
Step 4	for Steps 1-3. 6. 48 80 4. Surface Finish inner tone light	\$2674	2895	2911 3166 +\$0



Description


The center point of this 60"-wide panel attaches at a 90° angle to other panels and provides support to a panel run. Each panel has adjustable glides, finished ends, and a 48"-high recessed filler and hanger frames at 1" intervals to support components. The panel has an area between the filler and the rear edge of an attached work surface for passage of wires or tubes.


The panel is available in 2 heights. The 80"-high panel supports a locker; 48"- and 80"-high panels support work surfaces and other 60"-wide components.


Notes

To support locker, order 2 adapter rails (C1610.) separately.


To connect to panel at 90° angle, order draw rod (AO215.) separately. Order draw rods (AO215.) or seismic draw rods (CO215.) separately.







Specification Information		
Step 1.		
C1161		
Step 2	Height	
4860	48" high	
8060	8o" high	
Prices	for Steps 1-2.	
C1161	4860	\$1699
	8060	\$2090
Step 3	Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o



Description


This connector attaches a support panel to a wall or column for a peninsula configuration.


Notes

Order draw rod (AO215.) separately.

Package does not contain attachment hardware; designer must select fasteners appropriate for wall conditions and installer must field install.

Dimensions





SIDE FRONT


Specification Information Step 1. C1211.

Step 2. Height

48 48" high80 80" high

Prices for Steps 1-2.	
C1211. 48	\$301
80	\$360

Step 3	3. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+ \$o
LU	soft white	+ \$o
WN	warm grey neutral	+\$o



Description

This connector joins 2 support panels of equal or unequal heights at a 90° angle. It cannot be used with a panel run that supports Action Office® Series 2 work surfaces.


Notes


Order support panels (C1150., C1151., C1155., C1156.) separately.

Order 2 draw rods (AO215.) or 2 seismic draw rods (CO215.) separately.

When connecting panels of unequal heights, specify height of connector to match height of higher panel.

Dimensions





Specification Information

Step 1.

C1221.

Step 2.	Height
40	. 0// - : -

48 48" high


80 80" high

Prices for Steps 1-2.

C1221. 48 \$288 **80** \$346

Step 3. Surface Finish

HF	inner tone light	+\$o
LT	light tone	+ \$0
LU	soft white	+ \$0
WN	warm grev neutral	+\$o



Description

This connector joins 3 support panels of equal or unequal heights at a 90° angle. It cannot be used with a panel run that supports Action Office® Series 2 work surfaces.


Notes


Order support panels (C1150., C1151., C1155., C1156.) separately.

Order 3 draw rods (AO215.) or 3 seismic draw rods (CO215.) separately.

When connecting panels of unequal heights, specify height of connector to match height of higher panel.

Dimensions





Specification Information

Step 1.

C1231.


Step 2. Height

48 48" high80 80" high

Drices for Stens 1-

Prices for Steps 1-2.	
C1231. 48	\$316
80	\$389

Step	3. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$o



Description

This 11"-wide connector joins a module to a support panel. It has a removable cover strip that conceals predrilled holes used for attaching connectors.

Notes


To complete shelf surface at end of module run, order shelf end cap (C1510.E or C1511.E) separately.


Order panel connector (C1281.) and 3 draw rods (AO215. or CO215.) separately.

To attach support panel to interface connector (C1261.) of unequal height, order the following products separately:

- Panel connector (C1281.)
- Change of height draw rod (AO219.)
- Change of height panel/panel finished end (AO272.)

Dimensions





Specification Information

Step 1.

WN

C1261.

Cto	n a	Цаі	aht
Ste	μ Ζ.	1161	giii

4812 48" high

8012 80" high

Prices for Steps 1-2.

C1261. 4812	\$769
8012	\$938

Step 3. Surface Finish


warm grey neutral

1	2	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o

FRONT

+\$0

+\$o



Product Information

Description


This connector attaches a support panel to an interface connector.

Notes

Specify connector height to match panel height.


Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

Dimensions



Spe	cification Information	
Step	1.	
C128	1.	
Step	2. Height	
48	48" high	
80	8o" high	
Price	s for Steps 1-2.	
C128	1. 48	\$67
	80	\$87
Step	3. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$o

warm grey neutral



Description


This brace provides additional support to a straight panel run that has components hanging on 1 side. It includes a leveling glide.


Notes

To support 24"-deep work surface, specify 21"-wide mid-run support; to support 30"-deep work surface, specify 27"-wide mid-run support.

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

Dimensions





Specification Information

Step 1.

C1272.


Step :	2. Height			
48	48" high			
80	8o" high			

Step 3.	Width
21	21" wide

27 27" wide

Prices for Steps 1-3.		
	21	27
C1272. 48	\$524	580
80	\$656	683

Step 4	. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$ 0



Description

This brace provides additional support to a straight panel run that has components hanging on both sides. It includes leveling glides.

Notes


To support 24"-deep work surface, specify 41"-wide mid-run support; to support 30"-deep work surface, specify 53"-wide mid-run support. To support 24"-deep work surface on 1 side and 30"-deep work surface on other side, specify 53"-wide mid-run support.

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

Dimensions



TOP



23" (584) 23" (584) 23" (584) 23" (584) 23" (584) 41" WIDE 52-5/8" (1337) 41" WIDE SIDE SIDE

Specification Information

Step 1.

C1273.

Step 2. Height

48 48" high80 80" high


Step 3. Width

41 41" wide

53 53" wide


Prices for Steps 1-3.		
	41	53
C1273. 48	\$593	689
80	\$808	901

Step 4.	. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o



Description

This 20"-wide aluminum top cap attaches to 2 support panels in a straight connection. It provides additional support for a panel run.



Specification Information		
Step 1.		
C1274.20		\$127
Step 2.	Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o



Description

This 48"-high module attaches to a wall or connects to other freestanding modules. It encloses and organizes mechanical services. The module includes a top shelf and 4 adjustable glides. It has hanger frames at 1" intervals to support components.

Notes

Order access panels separately.


For L, T, or X configuration, order module corner connector separately:


- 2-way connector (C1220.)
- 3-way connector (C1230.)

Order draw rods (AO215.) or seismic draw rods (CO215.) separately. To attach modules of unequal height, order 2 change of height draw rods (AO219.) and 2 change of height panel/panel finished ends (AO272.) separately.

To enclose the end of a module when a terminal panel (C1122.) is not used, order the end closure assembly (C1451.) separately.

Dimensions





Specification Information


Step 1.

C1112.48

Step 2	2. Width	
24	24" wide	
30	30" wide	
36	36" wide	
48	48" wide	

Prices for S	Prices for Steps 1-2.				
C1112.48	24	\$1185			
	30	\$1394			
	36	\$1422			
	48	\$1552			

Step 3. Surface Finish		
JD	warm grey neutral Chemsurf®	+ \$o
KH	light tone Chemsurf®	+\$ o
KR	inner tone light Chemsurf®	+\$ o
LE	soft white Chemsurf®	+\$ o



Description

This 80"-high module attaches to a wall or connects to other freestanding modules. It encloses and organizes mechanical services and provides pass-through access. The module includes an interior shelf that is field adjustable to 30", 36", or 48" high. It has hanger frames at 1" intervals to support components and 4 adjustable glides.

Notes

Order access panels separately.

For L, T, or X configuration, order module corner connector separately:


- 2-way connector (C1220.)
- 3-way connector (C1230.)

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

To attach modules of unequal height, order 2 change of height draw rods (AO219.) and 2 change of height panel/panel finished ends (AO272.) separately.


To enclose the top section of a pass-through module, order the end closure panel (C1450.) and 2 change of height panel/panel finished ends (A0272.) separately.

To enclose the end of a module when a terminal panel (C1122., C1124.) is not used, order the end closure assembly (C1451. or C1452.) separately.



Spec	cification Information	
Step	1.	
C111	4.80	
Step	2. Width	
24	24" wide	
30	30" wide	
36	36" wide	
48	48" wide	
Prices	s for Steps 1-2.	
C111	4.80 24	\$1533
	30	\$1696
	36	\$1800
	48	\$1910
Step	3. Surface Finish	
JD	warm grey neutral Chemsurf®	+\$0
KH	light tone Chemsurf®	+\$0
KR	inner tone light Chemsurf®	+\$0
LE	soft white Chemsurf®	+\$0

+\$0



Description


This 48"-high, hard-surfaced end panel attaches to a freestanding module. It supports a module run or divides space. The panel includes a top shelf end cap, adjustable glides, and finished ends, and comes in 3 widths.

The 60"- and 72"-wide panel supports hanging components.

Notes

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

Dimensions



Specification Information

Step 1.

C1122.48

Sten 2	Width	/Position

36L 36" wide and left-hand panel36R 36" wide and right-hand panel

60 6o" wide72 72" wide

Prices for Steps 1-2.


C1122.48	36L	\$1573
	36R	\$1573
	60	\$1926
	72	\$2844

Step 3. Surface Finish

Otop	J. 24.1432	
JD	warm grey neutral Chemsurf®	+\$0
KH	light tone Chemsurf®	+\$o
KR	inner tone light Chemsurf®	+\$o
LE	soft white Chemsurf®	+\$o

Pass-Through Module Terminal Panel

C1122. C1124.




Product Information


Description


This 80"-high, hard-surfaced end panel attaches to a freestanding module. It supports a module run or divides space. The panel is available in 3 widths. The 60"- and 72"-wide, full-height panel supports hanging components. Each panel has an interior shelf end cap, adjustable glides, and finished ends.

Notes

To attach terminal panel to module, order 2 draw rods (AO215.) or 2 seismic draw rods (CO215.) separately.








FULL-HEIGHT INSIDE

Q)
=
3
*
\approx
2
•
S
®
Ü
\boldsymbol{z}
=
is.
₹′
0
ü
_

Speci	ification Information	
Step 1		
C112		
Step 2	2. Height	
2.80	full-height	
4.80	half-height	
Step 3	B. Width/Position	
For ful	ll-height (2.80)	
60	6o" wide	
72	72" wide	
For ha	ılf-height (4.80)	
36L	36" wide and left-hand panel	
36R	36" wide and right-hand panel	
60	6o" wide	
72	72" wide	
Prices	for Steps 1-3.	
C1122	2.80 60	\$2461
	72	\$2948
C1124	i.80 36L	\$1723
	36R	\$1723
	60	\$2385
	72	\$2940
Step 4	4. Surface Finish	
JD	warm grey neutral Chemsurf®	+\$0
KH	light tone Chemsurf®	+\$0
KR	inner tone light Chemsurf®	+\$0
LE	soft white Chemsurf®	+\$0



Description


This 48"-high, 48"-wide laminate panel can be field cut to fit a non-modular space in a run of freestanding modules. It includes 3 top shelf brackets, 2 floor base brackets, and 2 wall-attached wall start connectors. The panel cannot support components.

Notes


To cover both sides of a freestanding module run, specify 2 diplomat panels.

For diplomat panel used with low or pass-through module, order top shelf (C1510.) separately; shelf must be field modified.





Specification Information		
Step 1.		
C1146.		\$715
Step 2.	Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+\$o
LU	soft white	+\$o +\$o
WN	warm grey neutral	+ \$o



Description


This 72"-high tower encloses services between the ceiling and the top shelf of a low module, the interior shelf of a pass-through module, or a 2- or 3-way module corner connector. It includes a ceiling trim plate. The cover cannot be used with a 4-way module corner connector.

Notes


Order module or connector separately:

- 2-way module corner connector (C1220.)
- 3-way module corner connector (C1230.)
- Low module (C1112.)
- Pass-through module (C1114.)


Dimensions



TOP



Specification Information			
Step 1.			
C1460		\$700	
Step 2.	Surface Finish		
HF	inner tone light	+ \$o	
LT	light tone	+\$0 +\$0	
LU	soft white	+\$o	
WN	warm grey neutral	+\$o	



Description

This connector starts a peninsula configuration by attaching a low or pass-through module to a wall or column. It includes a shelf end cap.


Notes


For low module, specify 48"-high wall start; for pass-through module, specify 80"-high wall start.

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

Package does not contain attachment hardware; designer must select fasteners appropriate for wall conditions and installer must field install.

Dimensions





Specification Information

Step 1.

C1212.

Step 2. Height

48 48" high80 80" high

Prices for Steps 1-2. C1212. 48 \$535


\$576 Step 3. Surface Finish JD warm grey neutral Chemsurf® +\$0

 JD
 warm grey neutral Chemsurf®
 +\$0

 KH
 light tone Chemsurf®
 +\$0

 KR
 inner tone light Chemsurf®
 +\$0

 LE
 soft white Chemsurf®
 +\$0



Description


This connector joins 2 modules of equal or unequal heights at a 90° angle. It includes a top shelf and filler panels to enclose the utility chase.

Notes

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

To attach module of unequal height, order 2 change of height draw rods (AO219.) and 2 change of height panel/panel finished ends (AO272.) separately.

Dimensions



Specification Information

Step 1.

C1220.

Step 2. Height

48 48" high


80 80" high

Prices for Steps 1-2.

C1220. 48	\$1254
80	\$1504

Step 3. Surface Finish

JD	warm grey neutral Chemsurf®	+\$o
KH	light tone Chemsurf®	+\$o
KR	inner tone light Chemsurf®	+\$o
LE	soft white Chemsurf®	+\$o



Description


This connector joins 3 modules of equal or unequal heights at 90° angles. It includes a top shelf and filler panels to enclose the utility chase.

Notes

Order draw rods (AO215.) or seismic draw rods (CO215.) separately.

To attach module of unequal height, order 2 change of height draw rods (AO219.) and 2 change of height panel/panel finished ends (AO272.) separately.

Dimensions



Specification Information

Step 1.

C1230.


Step 2. Height

48 48" high

80 80" high


Prices for Steps 1-2.	
C1230. 48	\$1254
80	\$1507


Step	3. Surface Finish	
JD	warm grey neutral Chemsurf®	+\$0
KH	light tone Chemsurf®	+\$0
KR	inner tone light Chemsurf®	+\$0
LE	soft white Chemsurf®	+\$0




Description

This leg provides additional support for normal loading on a low or pass-through module. It has a leveling glide.





Specif	ication Information	
Step 1.		
C1275.		\$499
Step 2.	Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o




Description

This kit fastens a module to the floor. It includes a bolt, spacer, and plate. Package contains 10.


Notes

Anchor support kit is not acceptable for seismic anchoring. Source anchor component separately.

Dimensions








Specification Information

Step 1.

C1276.

\$285



Description

This 11"-deep shelf encloses the top of a low module service chase and provides a storage shelf. It has a Chemsurf® surface. The shelf can be drilled for fixtures or piping. It cannot attach to a pass-through module.

Dimensions



Specification Information


Step 1.

C1510.12

Step	2. Width	
12	12" wide	
24	24" wide	
30	30" wide	
36	36" wide	
48	48" wide	
Prices	s for Steps 1-2.	
C151	0.12 12	\$225
	24	\$301

1 11003 101 3	лерэ 1-2.	
C1510.12	12	\$225
	24	\$301
	30	\$377
	36	\$396
	48	\$407

Step 3.	. Surface Finish	
JD	warm grey neutral Chemsurf®	+\$0
KH	light tone Chemsurf®	+\$0
KR	inner tone light Chemsurf®	+\$0
LE	soft white Chemsurf®	+\$0




Description

This 11"-deep shelf fits inside a pass-through module. It encloses the top of a 48"-high service chase and provides a storage shelf. The shelf has a Chemsurf® surface. It can be drilled for fixtures or piping. Attachment hardware is included.

Notes

To install shelf higher than 48", order interior shelf and support kit (C1515.) separately.

Dimensions



Specification Information


Step 1.

C1511.12

					_
Step 2.	. Width				
24	24" wide				
30	30" wide				
36	36" wide				
48	48" wide				

Prices for S	Steps 1-2.	
C1511.12	24	\$358
	30	\$399
	36	\$424
	48	\$448


Step 3. Surface Finish		
JD	warm grey neutral Chemsurf®	+\$o
KH	light tone Chemsurf®	+\$ 0
KR	inner tone light Chemsurf®	+ \$o
LE	soft white Chemsurf®	+\$ 0



Description

This cap finishes the end of an interior or top shelf at the connection of a low or pass-through module and an interface connector. It has a Chemsurf® surface. Attachment hardware is included.

Dimensions



Specification Information

Step 1.

C151


Step 2. Depth/Usage

1.E 7" deep for interior shelf

0.E 12" deep for top shelf

Prices for Steps 1-2.	
C1511.E	\$197
C1510.E	\$193


Step	3. Surface Finish	
JD	warm grey neutral Chemsurf®	+ \$o
KH	light tone Chemsurf®	+\$o
KR	inner tone light Chemsurf®	+\$o
LE	soft white Chemsurf®	+\$0



Description

This 11"-deep interior shelf is installed within a pass-through module. Attachment hardware is included.

Dimensions



Specification Information


Step 1.

C1515.12

Step 2	. Width
24	24" wide
30	30" wide
36	36" wide
48	48" wide

Prices	for Steps 1-2.	
C1515	.12 24	\$706
	30	\$762
	36	\$847
	48	\$967

Step 3. Surface Finish			
JD	warm grey neutral Chemsurf®	+\$o	
KH	light tone Chemsurf®	+\$o	
KR	inner tone light Chemsurf®	+\$o	
LE	soft white Chemsurf®	+\$o	




Description


This plastic barrier attaches to the back of a top shelf and creates a backsplash. It is used with a module that sits against a wall. The barrier can also compensate for variations in a wall surface. Package contains 4.

Notes

For use with interior shelf, field modify shelf backsplash.

Dimensions





Specification Information


Step 1.

C1520.

Step 2.	Width			
24	24" wide			
30	30" wide			
36	36" wide			
48	48" wide			

Prices for Steps 1-2.	
C1520. 24	\$83
30	\$103
36	\$112
48	\$147


Step 3. Surface Finish				
HF	inner tone light	+\$o		
.T	light tone	+\$o		
.U	soft white	+\$o		
٧N	warm grey neutral	+\$o		



Description

This component fills the space between a work surface and the interior shelf for a pass-through module. It creates a continuous work area. Attachment hardware is included.

Dimensions



Specification Information

Step 1.

C153

Step 2. Usage

0. for top shelf

1. for interior shelf

Step 3. Width

24	24" wide
30	30" wide
36	36" wide
48	48" wide

Prices for Steps 1-3.	
C1530. 24	\$56
30	\$70
36	\$83
48	\$93
C1531. 24	\$57
30	\$82
36	\$94
48	\$108

Step 4	4. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+ \$o
LU	soft white	+ \$o
WN	warm grey neutral	+\$0



Description


This assembly of 2 access panels attaches to a freestanding module and encloses both sides of the service chase. The panels mount in a flush or recessed position and are removable for service access. Each panel has 2 knockouts for standard duplex electrical outlets. Brackets and mounting hardware are included. Shipped knocked down.

The assembly has 2 30"- or 2 42"-high access panels.

Notes


Order low module (C1112.) or pass-through module (C1114.) separately. For interior shelf mounted at 36" on pass-through module, specify 30"-high assembly.

For interior shelf mounted at 48'' on pass-through module, specify 42''-high assembly.





FRONT



Specification Information

Step 1.

C1412.

Step 2. Height

30 30" high

42 42" high

Step 3. Width

19 19¹/₄" wide

25 25¹/₄" wide

31 31¹/₄" wide

43 43¹/₄" wide

Prices for Steps 1-3.				
	19	25	31	43
C1412. 30	\$519	604	631	651
42	\$581	679	742	818

Step 2	4. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+ \$o
LU	soft white	+\$0
WN	warm grey neutral	+\$0




Description

This assembly of 4 access panels attaches to a freestanding module and encloses both sides of the service chase. The panels mount in a flush or recessed position and are removable for service access. Each panel has 2 knockouts (4 per module side) for standard duplex electrical outlets. Brackets and mounting hardware are included. Shipped knocked down. The 24"-high assembly has 4 12"-high access panels; the 42"-high assembly has 2 12"- and 2 30"-high access panels.

Notes


Order low module (C1112.) or pass-through module (C1114.) separately. For interior shelf mounted 30" high on pass-through module, specify 24"-high assembly.

For interior shelf mounted 48" high on pass-through module, specify 42"-high assembly.





FRONT



Combination Access Panel Assembly continued

Specification Information

Step 1.

C1422.

Step 2. Height

24 24" high

42 42" high

Step 3. Width

19 19¹/₄" wide


25 25¹/₄" wide

31 31¹/₄" wide

43 43¹/₄" wide

Prices for Steps 1-3.				
	19	25	31	43
C1422. 24	\$1093	1146	1200	1243
42	\$1161	1290	1433	1512

Step 2	4. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+ \$o
LU	soft white	+ \$o
WN	warm grey neutral	+ \$o



Description


This assembly attaches to a pass-through module and creates a low cup sink assembly in the lower service chase. It encloses both sides of the service chase and has removable panels for service access. The assembly includes 4 12"-high access panels, 1 shelf support frame, 1 interior shelf, 1 bottom chase cover, and attachment hardware. Each 12"-high panel has 2 knockouts for standard duplex electrical outlets. The assembly has a 15^{1} /4" space for a cup sink; this divides the service chase into 2 12" segments. Shipped knocked down.


Notes

Order pass-through module (C1114.) separately.

To enclose service chase of adjacent module, order center closure panel (C1455.1607) separately.

Customer must order sink and fixtures separately and field cut interior shelf to hold sink.





Low Cup Sink Access Panel Assembly continued


Specification Information

43

Step 1.		
C1432.	42	
Step 2.	Width	
19	19¹/₄" wide	
25	25 ¹ / ₄ " wide	
43	43 ¹ / ₄ " wide	
Prices f	for Steps 1-2.	
C1432.	42 19	\$1739
	25	\$1898

Step 3. Surface Finish		
JD	warm grey neutral Chemsurf®	+\$o
KH	light tone Chemsurf®	+\$o
KR	inner tone light Chemsurf®	+\$o
LE	soft white Chemsurf®	+ \$o

\$2098




Description

This 32"-high, 11"-wide panel encloses the top section of a pass-through module.


Notes

Order pass-through module (C1114.) separately.

Order 2 panel/panel change-of-height finished ends (AO272.) separately.



Speci	fication Information	
Step 1		
C1450.3211		\$124
Step 2	. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o



Description


This assembly encloses the end of a module when a terminal panel is not used

The end closure for a low module has a 48"-high panel and 2 48"-high end caps. The full-height end closure for a pass-through module has an 80"-high panel and 2 80"-high end caps. The half-height end closure for a pass-through module has a 48"-high panel and 2 80"-high end caps.


Notes

Order low module (C1112.) or pass-through module (C1114.) separately.


Dimensions



END CLOSURE FOR LOW MODULE




FULL-HEIGHT END CLOSURE FOR PASS-THROUGH MODULE



HALF-HEIGHT END CLOSURE FOR PASS-THROUGH MODULE

SIDE FRONT

- 1	200 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	
Step 1.		
C145		
C+	T	
Step 2.	туре	
1.4811	end closure for low module	
1.8011	1.8011 full-height end closure for pass-through module	
2.8011	half-height end closure for pass-through module	
Prices f	or Steps 1-2.	
C1451.4	4811	\$254
C1451.8011		\$333
C1452.8011		\$249
Step 3.	Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$o



Description

This panel attaches to the inside of a module run and encloses the end of a service chase. Attachment hardware is included.

Notes


To enclose service chase on either end of cup sink opening or 14"-high openings created by partial chases at panel connections, order 16"-high panel.

To enclose lower interior module chase, order 43"-high panel.

To enclose upper section of pass-through module, order 32"-high panel.

Dimensions





Specification Information

Step 1.


C1455.

C: 11 : 1 :		
Step 2. Height		
1607	16" high	
3207	32" high	

4307 43" high

Prices for Steps 1-2.	
C1455. 1607	\$586
3207	\$593
4307	\$597

Step 3. Surface Finish		
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+ \$o
WN	warm grey neutral	+\$o



Description

This aluminum rail attaches to 1 side of an equal-width module or support panel and provides 1 or 2 cable management channels. It cannot be used with an Action Office® Series 1 or 2 connector.

Notes


Order connection components separately:

- Wire chase rail bridge (C1326. or C1327.) to connect across mid-run support panel or back side of 3-way connector
- 90° inside corner wire chase rail (C1321. or C1322.) to connect around 90° inside corner

Order wire chase rail covers (C1315., C1316., C1317., or C1318.) and rail end cap (C1331. or C1332.) separately.

Receptacles must be customer supplied.

Dimensions



Specification Information

Step 1.

C131

Step 2. Type

- 3. double channel
- 4. single channel

Step 3. Width


- 24" wide
- **30** 30" wide
- **36** 36" wide
- **42** 42" wide A
- **48** 48" wide
- **60** 60" wide
- **72** 72" wide
- **96** 96" wide

Prices for Steps 1-3.	
C1313. 24	\$99
30	\$126
36	\$161
42	\$186
48	\$197
60	\$247
72	\$295
96	\$392
C1314. 24	 \$52
30	\$67
36	\$85
42	\$95
48	\$102
60	\$128
72	\$151
96	\$203

Step 4. Surface Finish

AN clear anodized aluminum

Wire Chase Rail, 90° Inside Corner C1321. C1322.



Product Information


Description

This aluminum rail attaches to 2 modules or support panels at a 90° inside corner. It cannot be used with an Action Office® Series 1 or 2 connector.

Notes

Rail includes 2 6" mitered corner covers. Order wire chase rail covers to cover the remaining length of each rail (C1315., C1316., C1317., or C1318.) separately.

Dimensions



Specification Information

Step 1.

C132

Step 2. Type

1. double channel

2. single channel

Ste	n	3.	W	id	tl	h

24 24" wide

30 30" wide

48 48" wide

Step 4. Depth

For 2	4" wide (24)			
24	24" deep	 	 	
30	30" deep			
48	48" deep			

_	wide (30)
	30" deep


48 48" deep

For 48" wide (48)

48 48" deep

Prices for Steps 1-4.			
	24	30	48
C1321. 24	\$355	381	454
30	_	\$404	478
48	_	_	\$553
C1322. 24	\$238	255	293
30	_	\$267	305
48	_	_	\$343

Step 5. Surface Finish		
AN	clear anodized aluminum	+\$0



Description


This bridge joins wire chase rails across a mid-run support panel or on the back side of a 3-way 90° connector. It allows the continuation of a cable management channel and includes a blank cover. The bridge cannot be used with an Action Office® Series 1 or 2 connector.

Notes

Order 1 of the following products separately:

- 3-way 90° connector (C1231.)
- Mid-run support (C1272. or C1273.)

Dimensions



Specification Information

Step 1.


C132

Step 2. Type

- **6.** double channel
- 7. single channel


Prices for Steps 1-2.	
C1326.	\$95
C1327.	\$95

Step 3	3. Surface Finish	
ΔΝ	clear anodized aluminum	±\$0



Description

This end cap finishes the exposed end of a wire chase rail. Package



Spec	cification Information	
Step :	1.	
C133		
Step :	2. Type	
1.	double channel	
2.	single channel	
Prices	s for Steps 1-2.	
C133	1.	\$86
C133	2.	\$86
Step	3. Surface Finish	
AN	clear anodized aluminum	+ \$o

C1318.

Product Information


Description

This cover encloses a wire chase rail.

Notes

Specify type to cover wire chase rail's power and/or data outlets. Specify covers in 6"- or 12"-wide increments to match rail width.

Dimensions



Specification Information

Step 1.

C131

Step 2. Type

5.06 6"-wide plain cover
5.12 12"-wide plain cover
6. duplex cutout cover
7. GFI cutout cover

8. gang cover


Prices for Steps 1-2.

C1315.06	\$19
C1315.12	\$23
C1316.	\$26
C1317.	\$27
C1318.	\$27

Step 3. Surface Finish

AN clear anodized aluminum +\$o

\$40




Product Information

Description


This 11/2"-wide cover wraps around the outside of a wire chase channel to cover gaps between adjoining channels and covers.

Notes

Order wire chase and covers separately.



Spe	Specification Information			
Step	1.			
C135	58.			
Step	2. Type			
D	double channel			
Т	triple channel			
Price	es for Steps 1-2.			
C135	58. D	\$2	7	



Description

This work surface hangs from a panel, module, or wall strips and is available as laminate or solid surface. The work surface, when supported by work surface supports has a top-loading capacity of up to 200 pounds. The solid surface option, when supported by sink products has a top-loading capacity of 250 pounds. Appropriate sink products and work surface supports must be specified and ordered separately.

Notes

For laminate option (L), order standard duty work surface supports (C2140.S) separately.

For solid surface option (K), order 1 of the following products separately:

- ADA sink support, systems mounted (C2520.)
- Sink module (C2512.)
- Sink unit for storage (C2513.)
- Work surface support, standard duty (C2140.S)

Work surface with solid surface option (K) will not accept a storage case (C5040.) without storage case leg (C5050.).

To support a 30" deep work surface installed on wall strips, a storage case (C5040.) with storage case leg (C5050.) is required. Order separately


Work surface with solid surface option (K) has a matte gloss level and includes a backsplash.

When hanging work surface with laminate option (L) from wall strips, order recommended backsplash (C2910.) separately.

Order optional products separately:

- Drawer bearer package (CO228)
- Gripper package (CO266), for laminate (L) option only
- Keyboard/Mouse Support, Adjustable (Y7710.)
- Lockable drawer bearer (CO396)

Dimensions



Specification Information

Step 1.

C2142.

Step	2.	Depth

24 24" deep30 30" deep

Step 3. Width

24 24" wide30 30" wide

36 36" wide48 48" wide

60 60" wide 72 72" wide

96 96" wide

Step 4. Material

For 24" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), or 60" wide (60)

L laminate


K solid surface

For 72" wide (72) or 96" wide (96)

L laminate

Prices for Steps 1-4.		
	L	K
C2142. 24 24	\$110	1397
30	\$136	1623
36	\$186	2036
48	\$250	2304
60	\$364	2758
72	\$456	_
96	\$587	_
30 24	\$117	1565
30	\$169	1816
36	\$213	2276
48	\$348	2566
60	\$494	3069
72	\$607	_
96	\$727	_

Step	5. Finish	
For la	minate (L)	
HF	inner tone light	+\$0
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$0
For so	olid surface (K)	
1N	canyon	+\$o
CQ	glacier white	+\$o
QB	bisque	+\$o
1E	flint	+\$50
1S	savannah	+\$50
1U	sahara	+\$50
31	mojave	+\$50
41	sandstone	+\$50
IQ	aurora	+\$50
OQ	silt	+\$50
Q1	vanilla	+\$50
7I	raffia	+\$200
BLP	blue pebble	+\$200



Description

This work surface is supported by heavy duty work surface supports, allowing it to hang from a panel, module, or wall strips. It has a top-loading capacity of up to 400 pounds. It can hold drawers and includes drawer bearers. Appropriate heavy duty work surface supports must be specified and ordered separately.

The work surface can be specified as laminate, Chemsurf®, stainless steel, Corian®, or phenolic resin. A 51"-wide surface has a 3" overhang on the left or right side to fill the gap created by a mid-run support, C-series 3-way connector, or interface connector. A 96"-wide surface includes 1 additional support bracket. The 24"-wide surface cannot accept components.

Notes

Order heavy duty work surface supports (C2140.H, C2140.W, C214.R) separately.

For stainless steel option (S), order 1 of the following products separately:

- Heavy duty work surface support, with component adapter (C2140.W)
- Heavy duty work surface support, component adapter ready (C2140.R) Work surface with phenolic resin option (P) or Corian option (K) must use work surface supports (C2140.) heavy duty, with component adaptor (W) or heavy duty, component adaptor ready (R) options. Finish color on phenolic resin option (P) is black.

Work surface with solid surface option (K) has a matte gloss level.

Work surface with phenolic resin option (P) or Corian option (K) includes a backsplash.

Work surface with phenolic resin option (P) is 1" thick and includes spacers for transitioning to adjacent, non-resin surfaces.

To support a 30" deep work surface installed on wall strips, a storage case (C5040.) with storage case leg (C5050.) is required. Order separately.

When hanging a work surface with laminate option (L) or Chemsurf® option (C) from wall strips, order recommended backsplash (C2910.) separately.

Work surface supports the following quantities of drawers, C frames, or storage cabinets:


- 24" wide: 1 drawer or 1 C frame
- 30", 36" wide: 1 drawer, 1 C frame or 1 storage cabinet
- 48", 51", and 60" wide: 2 drawers, 2 C frames, or 2 storage cabinets
- 72" wide: 3 drawers, 3 C frames, or 3 storage cabinets
- 96" wide: 4 drawers, 4 C frames, or 4 storage cabinets

24"-wide work surface will not accept a storage case (C5040.) or storage unit (C5010.).

Work surface with phenolic resin option (P) or stainless steel option (S) will not accept a keyboard/mouse support, adjustable (Y7710.).

Order optional products separately:

- Drawer (CO207, CO208, or CO209)
- Keyboard tray (Y7710.), except for stainless steel (S) and phenolic resin (P) surfaces
- Pencil drawer (Ao48o.)
- Storage case (C5040.) and adapter rail (C1610.), except for phenolic resin (P) surfaces
- Storage unit (C5010.) and component adaptor (C2193.)



Step	1.
C214	1.
Step	2. Depth
24	24" deep
30	30" deep
Step	3. Width
24	24" wide
30	30" wide
36	36" wide
48	48" wide
5L	51" wide, left overhang
5R	51" wide, right overhang
60	60" wide
72	72" wide
96	96" wide
Step	4. Material
For 2	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide
For 2 (60),	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72)
For 2 (60), L	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate
For 2 (60), L C	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf®
For 2 (60), L C S	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf® stainless steel
For 2 (60), L C	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf®
For 2. (60), L C S P	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf® stainless steel A phenolic resin corian
For 2 (60), L C S P K	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf® stainless steel A phenolic resin
For 2 (60), L C S P K	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf® stainless steel A phenolic resin corian 1" wide, left overhang (5L) or 51" wide, right overhang (5R)
For 2 (60), L C S P K	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf® stainless steel A phenolic resin corian 1" wide, left overhang (5L) or 51" wide, right overhang (5R) laminate Chemsurf®
For 2 (60), L C S P K	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf® stainless steel A phenolic resin corian 1" wide, left overhang (5L) or 51" wide, right overhang (5R) laminate
For 2 (60), L C S P K	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf® stainless steel A phenolic resin corian 1" wide, left overhang (5L) or 51" wide, right overhang (5R) laminate Chemsurf® stainless steel A
For 2 (60), L C S P K For 5	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf® stainless steel A phenolic resin corian 1" wide, left overhang (5L) or 51" wide, right overhang (5R) laminate Chemsurf® stainless steel A phenolic resin corian
For 2 (60), L C S P K For 5	4" wide (24), 30" wide (30), 36" wide (36), 48" wide (48), 60" wide or 72" wide (72) laminate Chemsurf® stainless steel A phenolic resin corian 1" wide, left overhang (5L) or 51" wide, right overhang (5R) laminate Chemsurf® stainless steel A phenolic resin


Prices for Steps 1-4.					
	L	C	S	P	K
C2141. 24 24	\$148	200	2021	716	1186
30	\$158	235	2175	846	1327
36	\$206	299	2598	976	1618
48	\$281	383	3103	1235	2182
5L	\$336	443	3380	1365	2477
5R	\$336	443	3380	1365	2477
60	\$391	494	4044	1430	2479
72	\$484	601	4615	1625	3063
96	\$618	753	5169	_	_
30 24	\$190	226	2276	780	1353
30	\$200	306	2568	976	1518
36	\$272	407	3267	1107	1836
48	\$366	518	3844	1430	2441
5L	\$381	534	4158	1625	2780
5R	\$381	534	4158	1625	2780
60	\$504	656	4935	1690	2773
72	\$618	784	5683	1949	3443
96	\$759	959	6067	_	_

Step 5	5. Finish		
For la	minate (L)		
HF	inner tone light	+\$0	
LT	light tone	+ \$o	
LU	soft white	+\$o	
WN	warm grey neutral	+\$0	
For Ch	nemsurf® (C)		
JD	warm grey neutral Chemsurf®	+\$0	
KH	light tone Chemsurf®	+\$0	
KR	inner tone light Chemsurf®	+ \$o	
LE	soft white Chemsurf®	+\$0	

stainless steel A

S

For cor	ian (K)	
CQ	glacier white	+\$0
Q1	vanilla	+\$0
QB	bisque	+\$0
1N	canyon	+\$370
15	savannah	+\$370
1U	sahara	+\$370
31	mojave	+\$370
41	sandstone	+\$370
IQ	aurora	+\$370
0Q	silt	+\$370
1E	flint	+\$740
7I	raffia	+\$740
BLP	blue pebble	+\$740



Description

This set of supports allows a standard- or heavy duty work surface to hang from a panel, module, or wall strips.

Notes

Supports for standard duty surfaces include the following:

Surface Width—Supports

24"-60" wide - 2 supports (1 left and 1 right)

72" wide — 3 supports (1 left, 1 middle, and 1 right)

96" wide - 4 supports (1 left, 2 middle, and 1 right)

Supports for heavy duty surfaces include the following:

Surface Width—Supports


24"-72" wide — 2 supports (1 left and 1 right)

96" wide — 3 supports (1 left, 1 middle, and 1 right)

For standard duty option (S), order standard duty work surface (C2142.) separately.

For heavy duty options (H,W, or R), order appropriate heavy duty work surface (C2141.) separately.

Dimensions



Specification Information

Step 1.

C2140.

Step 2. Type

S standard duty

H heavy duty

W heavy duty, with component adaptor

R heavy duty, component adaptor ready

Step 3. Depth

24 for 24"-deep work surface

30 for 30"-deep work surface

Step 4. Width

24 24" wide

30 30" wide

36 36" wide48" wide

60 60" wide

72 72" wide


96 96" wide

Prices for S	Steps 1-4.						
		24	30	36	48	60	72
C2140. S	24	\$153	153	153	153	153	229
	30	\$153	153	153	153	153	229
Н	24	\$213	213	213	213	213	213
	30	\$228	228	228	228	228	228
W	24	\$385	410	432	454	496	540
	30	\$423	449	469	492	534	577
R	24	\$277	277	277	277	277	304
	30	\$314	314	314	314	314	347

		96
C2140. S	24	\$304
	30	\$304
Н	24	\$319
	30	\$341
W	24	\$747
	30	\$804
R	24	\$455
	30	\$518

Step	5. Finish	
HF	inner tone light	+\$
LT	light tone	+\$
LU	soft white	+\$
WN	warm grey neutral	+\$

\$1159



Product Information

Description


This 24"-deep rectangular work surface hangs on a standard rail. It has predrilled holes for hanging drawers, suspended pedestals, or keyboard trays. The work surface can be installed, relocated, or removed without tools and without disturbing other rail components. The support arm finish matches the surface finish.

Notes

Order optional work surface backsplash (C2910.) separately.

Order optional keyboard/mouse support, adjustable (Y7710.) separately.

Dimensions



Specification Information Step 1. C2410.24

Step 2. Width


72

48 48" wide
72 72" wide

Prices for Steps 1-2.

C2410.24 48 \$936


Step	3. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o




Description

This cover provides a continuous surface for 2 back-to-back squarededge work surfaces hung on an open panel frame. It is 48" wide but can be field cut to fit a narrower opening.

|--|



fication Information	
.48	\$78
. Finish	
black umber	+\$ o
inner tone light	+ \$o
light tone	+ \$o
soft white	+ \$o
warm grey neutral	+\$0
	Finish black umber inner tone light light tone soft white



Description

This assembly attaches to a support panel, a module, or wall strips and provides support for a C frame and components hanging on the C frame. Attachment hardware is included. Shipped knocked down.


The assembly holds the following quantities of C frames:

- 24" and 30" wide: 1 C frame
- 48" and 60" wide: 2 C frames
- 72" wide: 3 C frames

Notes

For 51"-wide work surface, order 48"-wide C frame support.

Dimensions



Specification Information

Step 1.

C2192.

Step	2.	Work	Surf	ace	Depth	ı

for 24"-deep work surface

30 for 30"-deep work surface

Step 3. Width


24	24" wide
30	30" wide
36	36" wide

48 48" wide **60** 60" wide

72 72" wide

Prices for Steps 1-3. 72 24 30 36 48 60 C2192. 24 585 590 600 \$573 578 595 30 \$631 681 688 672 676 694

Step 4.	. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+ \$o
LU	soft white	+ \$o
WN	warm grey neutral	+\$0



Description

This kit attaches to a heavy-duty work surface support and allows C frames to hang below the work surface. Attachment hardware is included.

The kit holds the following quantities of C frames or storage units:

- 24" wide: 1 C frame
- 30" wide: 1 C frame or storage unit
- 48" and 60" wide: 2 C frames or storage units
- 72" wide: 3 C frames or storage units
- 96" wide: 4 C frames or storage units

Notes

The kit attaches to a work surface support (C2140.) with heavy duty, component adapter ready option (R).

For 51"-wide work surface, order 48"-wide kit. 96"-wide kit is shipped in 2 pieces.

 ${\it Order\ storage\ components\ separately:}$

- C frame (CO205)
- Storage unit (C5010.)

Dimensions



Specification Information

Step 1.

C2193.

Step 2. Width					
24	24" wide				
30	30" wide				
48	48" wide				
60	6o" wide				
72	72" wide				
96	96" wide				


Prices for Steps 1-2.	
C2193. 24	\$155
30	\$185
48	\$232
60	\$279
72	\$328
96	\$420

Step	3. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+ \$o
LU	soft white	+ \$o
WN	warm grey neutral	+\$o

WN

warm grey neutral

+\$0




Product Information


Description

This plastic barrier attaches to the back edge of a C-series rectangular or corner work surface or to a process table and creates a backsplash. Package contains 6.

Dimensions



Specification Information Step 1. C2910. Step 2. Width 72" wide Prices for Steps 1-2. C2910. 72 \$348 Step 3. Surface Finish For 72" wide (72) black **+\$**o HF inner tone light +\$o LT light tone **+\$**0 soft white LU **+\$**o



Description

This component attaches to a Casework work surface. It fills the gap created when Casework base cabinets are placed against terminal panels or wall strips. Package contains 6. Attachment hardware is included.

Notes


Order work surface (CG211. or CG212.) or corner work surface (CG218. or CG219.) separately.



Speci	fication Information	
Step 1.		
C2920	96	\$148
Step 2.	Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o

30

\$356



Product Information

Description


This side splash attaches to a Co/Struc phenolic resin work surface. It is used at the end of a run or at a corner. The finish is available in black only.

Notes

Order phenolic work surface (C2141) separately.



Spec	cification Information	
Step	1.	
C293	0.	
Step	2. Depth	
24	24" deep	
30	30" deep	
Step	3. Material	
P	phenolic resin	
Price	s for Steps 1-3.	
		P
C293	0. 24	\$331



Description


This 48"-wide work surface hangs from modules, C-series support panels, or wall strips and provides additional top-loading capacity. It includes 1 drawer bearer package. The work surface cannot attach to Action Office® Series 1 or Series 2 connectors. Attachment hardware is included.


The work surface has a laminate or Chemsurf® top. The support arm finish matches the surface finish.


Notes

When hanging work surface from wall strips, order recommended work surface backsplash (C2910.) separately.

Order optional keyboard/mouse support, adjustable (Y7710.) separately.








Speci	fication Information		
Step 1			
C2215	. A		
Step 2	. Depth		
2448	24" deep A		
3048	30" deep A		
Prices	for Steps 1-2.		
C2215	. 2448		\$767
	3048		\$888
Step 3	. Surface Finish		
Lamina	ate Top		
		2448	3048
HF	inner tone light A	+\$o	0
LT	light tone A	+\$o	0
LU	soft white A	+\$o	0
WN	warm grey neutral A	+\$0	0
Chems	urf® Top		
		2448	3048
JD	warm grey neutral Chemsurf® 🛕	+\$167	194
KH	light tone Chemsurf® A	+\$167	10/
	tight tone chemsun	+\$10/	194
KR	inner tone light Chemsurf® A	+\$167	194

Heavy-Duty Corner Work Surface, C2216. Action Office® Series 1 or 2 90° Corner Connector



Product Information

Description


This 48"-wide work surface hangs from panels or wall strips connected by an Action Office® Series 1 or 2 90° corner connector. It provides additional top-loading capacity and includes 1 drawer bearer package. The support arm finish matches the surface finish. The work surface cannot attach to a C-series 90° corner connector. Attachment hardware is included.


The work surface has a laminate or Chemsurf® top.


Notes

When hanging work surface from wall strips, order recommended work surface backsplash (C2910.) separately.

Order optional keyboard/mouse support, adjustable (Y7710.) separately.








Heavy-Duty Corner Work Surface, Action Office® Series 1 or 2 90° Corner Connector continued

C	C - 1 - 1 - 5 1 - 5 1		
Step 1	fication Information		
C2216			
Step 2	. Depth		
2448	24" deep 🖪		
3048	30" deep A		
Prices	for Steps 1-2.		
C2216	. 2448		\$767
	3048		\$888
Step 3	. Surface Finish		
Lamin	ate Top		
		2448	3048
HF	inner tone light A	+\$o	0
LT	light tone A	+\$o	0
LU	soft white A	+\$o	0
WN	warm grey neutral A	+\$0	0
Chems	urf® Top		
		2448	3048
JD	warm grey neutral Chemsurf® A	+\$167	194
KH	light tone Chemsurf® A	+\$167	194
KR	inner tone light Chemsurf® A	+\$167	194
LE	soft white Chemsurf® A	+\$167	194



Description

This 30"-deep height-adjustable table has a laminate or Chemsurf® top. The electrically driven assembly has a 15" travel range from $26^1/_2$ " to $41^1/_2$ ". The table base is equipped with a motor, height-control box, and $8^1/_2$ ' power cord. It has a load capacity of 400 pounds when weight is evenly distributed. The table is available with leveling glides or 4" casters with brakes. It is UL listed for the U.S. and Canada. Attachment hardware is included. Shipped knocked down.

A riser frame attaches to the table base to provide upright support for hanging components. The table includes 1 stretcher to support a C frame or storage unit. The 48"-wide workstation has 1 set of drawer bearers; 60"- and 72"-wide workstations have 2 sets of drawer bearers.

Notes

60"- and 72"-wide riser frames require a riser frame divider for overhead storage components. See Co/Struc Planning Guide for information.

Order storage components (to hang on stretcher) separately:

- B-style shelf (CO521.)
- C frame (CO205)
- Drawer (CO207, CO208, or CO209)
- Storage unit (C5010.)

Order workstation storage components separately:

Reagent shelf (C4327.)


Order overhead components (to hang on riser frame) separately:

- B-style shelf (C3011.)
- Equipment rail (C1172.)
- Flat panel mount (C4324.)
- Marker board (C4323.)
- Power strip (C4322.)
- Riser frame divider (C4320.)
- Riser frame end cover (C4321.)
- Shelf unit (C5020. or C5021.)

Order Action Office® overhead components separately. See Co/Struc Planning Guide for information.

When specifying a Chemsurf® finish, the corresponding frame finish is as follows:

- For warm grey neutral Chemsurf (JD), warm grey neutral (WN) frame.
- For light tone Chemsurf (KH), light tone (LT) frame.
- For inner tone light Chemsurf (KR), inner tone light (HF) frame.
- For soft white Chemsurf (LE), soft white (LU) frame.




Workstation, Electric Height Adjustment continued

Specification Information

Step 1	l .		
C4310	0.30		
Step 2	2. Width		
48	48" wide		
60	6o" wide		
72	72" wide		
Step 3	3. Surface Material		
L	laminate		
С	Chemsurf®		
Prices	for Steps 1-3.		
		L	
C4310		\$6097	638
	60	\$6178	646
	72	\$6459	674
Step 2	4. Casters/Glides		
GA	glides		+\$
CA	casters		+\$37
Step	5. Top/Frame Finish		
For la	minate (L)		
HF	inner tone light		+\$
LT	light tone		+\$
LU	soft white		+\$
WN	warm grey neutral		+\$
For Ch	nemsurf® (C)		
JD	warm grey neutral Chemsurf®		+\$
KH	light tone Chemsurf®		+\$
KR	inner tone light Chemsurf®		+\$
LE	soft white Chemsurf®		+\$

Workstation, Manual Height Adjustment

C4311.



Product Information

Description

This 30"-deep height-adjustable work station has a laminate or Chemsurf® top. It is pin adjustable in 1" increments and has a 15" travel range from $26^1/_2$ " to $41^1/_2$ ". The table has a load capacity of 400 pounds when weight is evenly distributed. The table is available with leveling glides or 4" casters with brakes. It is UL listed for the U.S. and Canada. Attachment hardware is included. Shipped knocked down.

A riser frame attaches to the table base to provide upright support for hanging components. The table includes 1 stretcher to support a C frame or storage unit. The 48"-wide workstation has 1 set of drawer bearers; 60"- and 72"-wide workstations have 2 sets of drawer bearers.

Notes

60"- and 72"-wide riser frames require a riser frame divider for overhead storage components. See Co/Struc Planning Guide for information.

Order storage components (to hang on stretcher) separately:

- B-style shelf (CO521.)
- C frame (CO205)
- Drawer (CO207, CO208, or CO209)
- Storage unit (C5010.)

Order workstation storage components separately:

• Reagent shelf (C4327.)


Order overhead components (to hang on riser frame) separately:

- Equipment rail (C1172.)
- Flat panel mount (C4324.)
- Marker board (C4323.)
- Power strip (C4322.)
- Riser frame divider (C4320.)
- Riser frame end cover (C4321.)

Order Action Office® overhead components separately. See Co/Struc Planning Guide for information.

When specifying a Chemsurf® finish, the corresponding frame finish is as follows:

- For warm grey neutral Chemsurf (JD), warm grey neutral (WN) frame.
- For light tone Chemsurf (KH), light tone (LT) frame.
- For inner tone light Chemsurf (KR), inner tone light (HF) frame.
- For soft white Chemsurf (LE), soft white (LU) frame.




Workstation, Manual Height Adjustment continued

Spec	ification Information		
Step 1			
C4311	.30		
Cton	. NA/: - 41-		
	2. Width		
48	48" wide		
60	60" wide		
72	72" wide		
Step 3	3. Surface Material		
L	laminate		
С	Chemsurf®		
Prices	for Steps 1-3.		
		L	(
C4311	.30 48	\$4635	4946
	60	\$4742	505
	72	\$4962	5364
Step 2	. Casters/Glides		
GA	glides		+\$0
CA	casters		+\$370
Step	;. Top/Frame Finish		
For la	minate (L)		
HF	inner tone light		+\$0
LT	light tone		+\$0
LU	soft white		+\$0
WN	warm grey neutral		+\$0
For Ch	nemsurf® (C)		
JD	warm grey neutral Chemsurf®		+\$0
KH	light tone Chemsurf®		+\$0
KR	inner tone light Chemsurf®		+\$0
LE	soft white Chemsurf®		+\$0

Work Table, Electric Height Adjustment

C4312.



Product Information

Description

This 2-sided, 42"-deep height-adjustable table has a laminate or Chemsurf® top. The electrically driven assembly has a 15" travel range from $26^{1}/2$ " to $41^{1}/2$ ". The table base is equipped with a motor, heightcontrol box, and $8^{1/2}$ power cord. It has a load capacity of 400 pounds when weight is evenly distributed. The table is available with leveling glides or 4" casters with brakes. It is UL listed for the U.S. and Canada. Attachment hardware is included. Shipped knocked down.

The table includes 2 stretchers (one on each side) to support C frames and/or storage units. 48"- and 60"-wide tables have 4 sets of drawer bearers.


Notes


Order storage components separately:


- B-style shelf (CO521.)
- C frame (CO205)
- Drawer (CO207, CO208, or CO209)
- Storage unit (C5010.)

When specifying a Chemsurf® finish, the corresponding frame finish is as follows:

- For warm grey neutral Chemsurf (JD), warm grey neutral (WN) frame.
- For light tone Chemsurf (KH), light tone (LT) frame.
- For inner tone light Chemsurf (KR), inner tone light (HF) frame.
- For soft white Chemsurf (LE), soft white (LU) frame.








Work Table, Electric Height Adjustment continued

Spec	fication Information		
Step 1			
C4312	.42		
Stop	Width		
48	Width 48" wide		
60	60" wide		
60	60 wide		
Step 3	. Surface Material		
L	laminate		
С	Chemsurf®		
Prices	for Steps 1-3.		
		L	C
C4312	.42 48	\$5615	5910
	60	\$5792	6087
Step 4	. Casters/Glides		
GA	glides		+\$0
CA	casters		+\$370
Step 5	. Top/Frame Finish		
For lai	minate (L)		
HF	inner tone light		+\$0
LT	light tone		+\$0
LU	soft white		+\$0
WN	warm grey neutral		+\$0
For Ch	emsurf® (C)		
JD	warm grey neutral Chemsurf®		+\$0
KH	light tone Chemsurf®		+\$0
KR	inner tone light Chemsurf®		+\$0
LE	soft white Chemsurf®		+\$0

Process Table, Electric Height Adjustment

C4313.



Product Information

Description

This height-adjustable process table has a laminate or Chemsurf® top. The electrically driven assembly has a 15" travel range from $26^{1}/_{2}$ " to $41^{1}/_{2}$ ". The table base is equipped with a motor, height-control box, and $8^{1}/_{2}$ ' power cord. It has a load capacity of 400 pounds when weight is evenly distributed. The table is available with leveling glides or 3" locking casters. It is UL listed for the U.S. and Canada. Attachment hardware is included. Shipped knocked down.

The table includes 1 stretcher to support a C frame or storage unit.

Notes


Order storage components separately:


- B-style shelf (CO521.)
- C frame (CO205)
- Drawer (CO207, CO208, or CO209)
- Reagent shelf (C4326.)
- Storage unit (C5010.)


When specifying a Chemsurf® finish, the corresponding frame finish is

- For warm grey neutral Chemsurf (JD), warm grey neutral (WN) frame.
- For light tone Chemsurf (KH), light tone (LT) frame.
- For inner tone light Chemsurf (KR), inner tone light (HF) frame.
- For soft white Chemsurf (LE), soft white (LU) frame.

Dimensions







Specification Information

Step 1.

C4313.

Step 2.	Depth
24	24" deep
30	30" deep

Step 3. Width		
48	48" wide	
60	6o" wide	

72 72" wide

Step 4. Surface Material


Step 6. Top/Frame Finish

L laminateC Chemsurf®

Prices for Steps 1-4.		
	L	C
C4313. 24 48	\$5385	5679
60	\$5561	5854
72	\$5796	6088
30 48	\$5486	5778
60	\$5601	5894
72	\$5831	6125

Step 5	. Casters/Glides	
GA	glides	+\$0
CA	casters	+\$310

For la	minate (L)	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0
For Cl	hemsurf® (C)	
JD	warm grey neutral Chemsurf®	+\$0
KH	light tone Chemsurf®	+\$0
KR	inner tone light Chemsurf®	+\$0
LE	soft white Chemsurf®	+\$0



Description

This metal shelf attaches to a workstation or process table to provide storage for reagents below the surface.


Notes

Order workstation (C4310. or 4311.) separately.

Order process table (C4313.) separately.

For Workstation

Dimensions



Specification Information

Step 1.

C432

Step 2. Type

for electric process table 6.

7. for workstation

Step 3. Depth

For electric process table (6.)

16 16" deep

For workstation (7.)

18" deep

Step 4. Width


48 48" wide

60 60" wide

72" wide 72

Prices for Steps 1-4.			
	48	60	72
C4326. 16	\$749	764	814
C4327. 18	\$784	800	853

Step 5.	Finish	
HF	inner tone light	+ \$o
LT	light tone	+ \$o
LU	soft white	+ \$o
WN	warm grey neutral	+ \$o



Description

This freestanding table base is available at seated height (29" high) and standing height (36" high), has 2 C-legs, a steel modesty panel, and 1" adjustable glides. Base finishes come in both metallic and standard healthcare finishes. Work table surface and available reagent shelf are ordered separately. The table has a load capacity of 500 pounds when weight is evenly distributed. Shipped knocked down.

Notes

If casters are desired at initial purchase, Options is recommended as a solution to maintain 29" and 36" heights. Casters are also available as a service part if replacing glides with casters is desired in the future.

At lowest glide adjustment, actual height of table with glides and work surface is 283/4" for a 29" seated-height table base and 353/4" for a 36" standing-height table base. Actual height of table with casters and work surface is 31.4" for a 29" seated-height table base and 38.4" for a 36" standing-height table base. Order casters as a service part.

Clearance from front of table to modesty panel is 241/4" for a 30"-deep work surface and 181/4" for a 24"-deep work surface.

Order work table surface (C4700.) separately:


- 24" and 30" depths, and 48", 51" (left and right), and 60" widths are compatible with table bases 24" deep \times 48" wide.
- 60" wide work table surface overhangs on center of a 48"-wide table base.
- 51" wide (left and right) work table surfaces overhang left or right of a 48"-wide table base.
- Order drawer bearers separately, with work table surface.

Order separately:


- Work table reagent shelf (C4900.).
- Cord cleat (G1331.).
- Cable management trough (AO382.).

To accommodate a refrigerator, a 36"-high table base is recommended. Refrigerators 21.34" wide or less can be used with 1 set of adjacent drawer bearers for a 3" or 6" drawer (CO207, CO208) only.

Due to the nature of the paint for metallic finishes (CN, MS), the end finish is more prone to show a light blemish or texture change on the top of the table foot when light hits it directly. This is normal as a result of the metal transition from the weld underneath.



Spec	ification Information		
Step 1			
C460	Α		
Step 2	2. Glides		
G.	glides A		
Step 3	3. Depth		
24	24" deep A		
Step 2	4. Width		
48	48" wide A		
Step 5	5. Height		
29	29" high 🛕		
36	36" high 🛕		
Prices	for Steps 1-5.		
		29	36
C4600	G. 24 48	\$942	994
Step 6	6. Base/Modesty Panel Finish		
HF	inner tone light A		+\$0
LT	light tone A		+\$0
LU	soft white A		+\$0
WN	warm grey neutral A		+\$0
CN	metallic champagne A		+\$20
MS	metallic silver A		+\$20



Description

This freestanding table base is available at seated height (29" high) and standing height (36" high) and has stretchers to support C frames, drawers, and storage units. Table base comes with 1" adjustable glides. Base finishes come in both metallic and standard healthcare finishes. Work table surface and available reagent shelf are ordered separately. The table has a load capacity of 500 pounds when weight is evenly distributed. Shipped knocked down.

Notes

If casters are desired at initial purchase, Options is recommended as a solution to maintain 29" and 36" heights. Casters are also available as a service part if replacing glides with casters is desired in the future.

At lowest glide adjustment, actual height of table with glides and work surface is $28^3/4$ " for a 29" seated-height table base and $35^3/4$ " for a 36" standing-height table base. Actual height of table with casters and work surface is 31.4" for a 29" seated-height table base and 38.4" for a 36" standing-height table base.

Order work table surface (C4700.) separately:

- 24" deep \times 36" wide base is only compatible with a 24" deep \times 36" wide work table surface.
- 24" deep \times 48" wide base is compatible with a 24" deep \times 48", 51" (left and right), 60" wide work table surface.
- 24" deep × 72" wide base is only compatible with a 24" deep × 72" wide work table surface.
- 30" deep \times 36" wide base is only compatible with a 30" deep \times 36" wide work table surface.
- 30" deep \times 48" wide base is compatible with a 30" deep \times 48", 51" (left and right), 60" wide work table surface.
- 30" deep × 485" wide base is only compatible with a 36" deep × 72" wide work table surface.
- 30" deep \times 72" wide base is only compatible with a 30" deep \times 72" wide work table surface.

60" wide work table surface (C4700.) overhangs on center of a 48"-wide table base.

51" wide (left and right) work table surfaces (C4700.) overhang left or right of a 48"-wide table base.

Order drawer bearers separately, with work table surface (C₄700.). Due to the nature of the paint for metallic finishes (CN, MS), the end finish is more prone to show a light blemish or texture change on the top of the table foot when light hits it directly. This is normal as a result of the metal transition from the weld underneath.

Order storage components separately:

- Drawer (CO207, CO208, CO209).
- Storage unit (C5010.).
- Work table reagent shelf (C4900.).
- C frame (CO205).

Maximum C frames or storage units using maximum number of stretchers and a 36"-high table base:

- For 36"-wide base: 2 C frames or 1 storage unit.
- For 48"-wide base: 4 C frames or 2 storage units.
- For 48S"-wide × 30"-deep base: 8 C frames (4 front/4 back) or 4 storage units.
- For 72"-wide base: 6 C frames or 3 storage units.

Order technology support and cord management accessories separately:


- Connect[™] S300 (Y1323.xxxxS).
- Flo® single-screen monitor arm (Y91171.CM).
- Flo plus monitor arm assembly (Y91174.xCM).
- Tablet support (Y91173.).
- Laptop mount (Y91177.).
- LE series keyboard solutions (Y7710.1).
- JAW™ adjustable CPU holder assembly (Y91103.).
- Cord cleat (G1331.).
- Cable management trough (AO382.).
- Cable trough (Y1200.).

Order Y1200.30 for a C461G.3036; order Y1200.42 for a C461G.3048.

Dimensions

Spec	cification Information			
Step	1.			
C461	A			
Step	2. Glides			
G.	glides A			
Step	3. Depth			
24	24" deep 🖪			
30	30" deep A			
Step	4. Width			
For 2	4" deep (24)			
36	36" wide 🗚			
48	48" wide A			
72	72" wide A			
 For 30	o" deep (30)			
36	36" wide 🛕			
48	48" wide 🛕			
72	72" wide 🛕			
485	48" wide, symmetrical support [A]			
Step	5. Height			
29	29" high 🖪			
36	36" high A			
Step	6. Number of Stretchers			
For 2	9" high (29)			
A	upper placement only A			
For 3	6" high (36)			
Α	upper placement only A			
В	upper and lower placement A			
Prices	s for Steps 1-6.		_	
		29A	36A	36B
C461	G. 24 36	\$1064	1127	1312

Step	7. Base/Stretcher Finish	
HF	inner tone light A	+\$0
LT	light tone A	+\$0
LU	soft white A	+\$0
WN	warm grey neutral 🛕	+\$0
CN	metallic champagne 🛕	+\$20
MS	metallic silver A	+\$20



Description

This freestanding 1" pin-adjustable table base has stretchers to support C frames, drawers, and storage units. Table base comes with 1" adjustable glides or 4" locking casters. Base finishes come in both metallic and standard healthcare finishes. Work table surface and available reagent shelf are ordered separately. The table has a load capacity of 500 pounds when weight is evenly distributed. Shipped knocked down.

Notes

Table with glides and work surface adjusts from 27" high to 36" high. Table with casters adjusts from 30" high to 38" high.

Order work table surface (C4700.) separately:

- 24" deep x 36" wide base is only compatible with a 24" deep x 36" wide work table surface.
- 24" deep × 48" wide base is compatible with a 24" deep × 48", 51" (left and right), 60" wide work table surface.
- 24" deep \times 72" wide base is only compatible with a 24" deep \times 72" wide work table surface.
- 30" deep × 36" wide base is only compatible with a 30" deep × 36" wide work table surface.
- 30" deep \times 48" wide base is compatible with a 30" deep \times 48", 51" (left and right), 60" wide work table surface.
- 30" deep \times 48S" wide base is only compatible with a 36" deep \times 72" wide work table surface.
- 30" deep \times 72" wide base is only compatible with a 30" deep \times 72" wide work table surface.

 60° wide work table surface (C4700.) overhangs on center of a 48° -wide table base.

51" wide (left and right) work table surfaces (C4700.) overhang left or right of a 48"-wide table base.

Order drawer bearers separately, with work table surface (C4700.).

- Order storage components separately: Drawer (CO207, CO208, CO209).
- Storage unit (C5010.).
- Work table reagent shelf (C4900.).
- C frame (CO205).

Maximum C frames or storage units using maximum number of stretchers and a 36"-high table base:

- For 36"-wide base: 2 C frames or 1 storage unit.
- For 48"-wide base: 4 C frames or 2 storage units.
- For 485"-wide × 30"-deep base: 8 C frames (4 front/4 back) or 4 storage units.
- For 72"-wide base: 6 C frames or 3 storage units.


For tables adjusted below 32" high, C frames hang from top stretcher only, and storage units cannot mount.

Due to the nature of the paint for metallic finishes (CN, MS), the end finish is more prone to show a light blemish or texture change on the top of the table foot when light hits it directly. This is normal as a result of the metal transition from the weld underneath.

Order technology support and cord management accessories separately:

- Connect™ S300 (Y1323.xxxxS).
- Flo® single-screen monitor arm (Y91171.CM).
- Flo plus monitor arm assembly (Y91174.xCM).
- Tablet support (Y91173.).
- Laptop mount (Y91177.).
- LE series keyboard solutions (Y7710.1).
- JAW™ adjustable CPU holder assembly (Y91103.).
- Cord cleat (G1331.).
- Cable management trough (AO382.).
- Cable trough (Y1200.).

Order Y1200.30 for a C462x.3036; order Y1200.42 for a C462x.3048.




Specification Information

A
. Casters/Glides
glides A
casters A
. Depth
24" deep A
30" deep A
. Width
" deep (24)
36" wide A
48" wide A
72" wide A
/2 wide A
" deep (30)
36" wide A
48" wide A
72" wide A
48" wide, symmetrical support
. Number of Stretchers
upper placement only A
upper and lower placement A

Ston 6	. Base Finish	
		φ.
HF	inner tone light A	+\$ 0
LT	light tone A	+\$o
LU	soft white A	+\$o
WN	warm grey neutral 🛕	+\$o
CN	metallic champagne A	+\$20
MS	metallic silver A	+\$20

Prices for Steps	1-5.	
	Α	В
C462G. 24 36	\$1695	1886
48	\$1707	1899
72	\$1803	2079
30 36	\$1712	1901
48	\$1723	1914
72	\$1818	2094
485	\$1970	2214
	Α	В
C462C. 24 36	\$1762	1953
48	\$1773	1966
72	\$1870	2146
30 36	\$1777	1967
48	\$1790	1981
72	\$1885	2161
485	\$2036	2280



Description

This work table surface attaches to work table bases (C460G., C461G., C462G., C462C.). Predrilled holes accept drawer bearers ordered with the work surface. Surface material is available in laminate, Chemsurf®, Corian®, or phenolic resin.

Notes

Drawer bearers must be field installed.


Only used on work table bases C460G., C461G., C462G., C462C., ordered separately. See table base information for work surface compatibilities and top load weight capacities.

Work table surface supports the following maximum pairs of drawer bearers (each pair has 1 left and 1 right bearer):

- For 36"-wide table: 1 pair.
- For 48"-wide table: 2 pair.
- For 485"-wide table: 4 pair (2 pair on either side).
- For 60"-wide table: 2 pair.
- For 72"-wide table: 3 pair.

Phenolic resin is 1" thick. When ordering adjustable work table base (C462C.) with phenolic resin top (C4700.xxxxP) that needs to be level with tables (C462C.) using laminate, Chemsurf, or Corian tops (C4700.L, C4700.K, C4700.C), order phenolic shim service part.

See Co/Struc® System planning guide for capacity information.



Step 1. **C4700.** A

Specification Information

C-17 0 C					
Step 2	2. Depth				
24	24" deep 🖪				
30	30" deep 🛕				
36	36" deep 🖪				
Step 3	3. Width				
For 24	i" deep (24)				
36	36" wide 🛕				
48	48" wide 🛕				
51L	51" wide, left overhang A				
51R	51" wide, right overhang A				
60	60" wide 🛕				
72	72" wide 🗚				
For 30	o" deep (30)				
36	36" wide 🛕				
48	48" wide 🛕				
51L	51" wide, left overhang 🛕				
51R	51" wide, right overhang A				
60	60" wide 🛕				
72	72" wide 🗚				
 For 36	5" deep (36)				
72	72" wide 🗚				
Step 2	4. Surface Material				
L	laminate 🖪				
C	Chemsurf® A				
K	corian 🖪				
P	phenolic resin A				
Prices	for Steps 1-4.				
		L	С	K	P
C4700). 24 36	\$214	281	726	905
	48	\$233	306	909	1203
	51L	\$295	386	1092	1458
	51R	\$295	386	1092	1458
	60	\$303	395	1093	1460
	72	\$327	485	1278	1760


	30 36	\$240	307	857	1116
	48	\$269	336	1067	1480
	51L	\$331	420	1286	1798
	51R	\$331	420	1286	1798
	60	\$327	483	1279	1791
	72	\$370	517	1494	2159
	36 72	\$491	802	1745	3299
Step 5	. Top/Edge Finish				
For lar	minate (L)				
HF	inner tone light 🛕				+\$o
LT	light tone A				+\$o
LU	soft white A				+\$o
WN	warm grey neutral A				+\$ 0
For lar	ninate (L)				
28	canyon 🖪				+\$ 0
29	misted A				+\$0
For Ch	emsurf® (C)				
JD	warm grey neutral Chemsurf® [Α			+\$ 0
KH	light tone Chemsurf® A				+\$ 0
KR	inner tone light Chemsurf® A				+\$0
LE	soft white Chemsurf® A				+\$0
For co	rian (K)				
CQ	glacier white [A]				+\$0
Q1	vanilla 🖪				+\$0
QB	bisque A				+\$0
1N	canyon A				+\$400
15	savannah 🗚				+\$400
1U	sahara 🗚				+\$400
31	mojave A				+\$400
41	sandstone A				+\$400
IQ	aurora A				+\$400
0Q	silt A				+\$400
1E	flint A				+\$525
7I	raffia 🖪				+\$525
BLP	blue pebble A				+\$525

For phenolic resin (P)

BK black A

+\$0

Step	6. Pairs of Drawer Bearers	
For 36	5" wide (36)	
NB	without drawer bearers A	+\$0
1D	1 pair drawer bearers A	+\$15
For 48	8" wide (48), 51" wide, left overhang (51L), 51" wide, right	
overh	ang (51R), or 60" wide (60)	
NB	without drawer bearers A	+\$o
1D	1 pair drawer bearers A	+\$15
2D	2 pair drawer bearers A	+\$30
For 22	q" deep (24) or 30" deep (30) with 72" wide (72)	
NB	without drawer bearers A	+\$o
1D	1 pair drawer bearers A	+\$15
2D	2 pair drawer bearers A	+\$30
3D	3 pair drawer bearers A	+\$45
For 36	5" deep (36) with 72" wide (72)	
NB	without drawer bearers A	+\$0
1D	1 pair drawer bearers 🛕	+\$15
2D	2 pair drawer bearers A	+\$30
3D	3 pair drawer bearers A	+\$45
4D	4 pair drawer bearers A	+\$60



Description

This high-pressure laminate shelf attaches to work table bases to provide storage below the surface.

Notes


Compatible with basic work table base (C460G.), standard work table base (C461G.xxxxxxxxA), and adjustable work table base with upper stretcher placement only (C462x.xxxxxxA).

Compatible only with 36 high table bases, due to installation clearance needed.

16" deep shelf (C4900.16) is compatible with 24"- and 30"-deep table bases. 18" deep shelf (C4900.18) is compatible with only 30"-deep table bases.

Order warm grey neutral (WN) reagent shelf finish to complement a metallic champagne (CN) table base. Order cool grey neutral (CL) reagent shelf finish to complement a metallic silver (MS) table base. Reagent shelf (C4900.) rests on table base legs and is not height adjustable.

Dimensions



Specification Information

Step 1.

C4900. A

Step	2.	Depth
------	----	-------

16" deep A

18 18" deep A

Step 3. Width

36 36" wide A

48 48" wide A

72 72" wide A


Step 4. Surface Material

L laminate A

Prices for Steps 1-4.	
	L
C4900. 16 36	\$376
48	\$425
72	\$615
18 36	\$395
48	\$443
72	\$635

Step	5. Top/Edge Finish	
HF	inner tone light A	+\$o
LT	light tone A	+\$o
LU	soft white A	+\$o
WN	warm grey neutral A	+\$o
CL	cool grey neutral A	+\$o

\$346



Product Information


Description

This piece divides a workstation riser frame into smaller sections and provides upright riser support for hanging components. The dividers bolt to the riser frame at predrilled attachment locations. Finish is black.

Notes

Order workstation (C4310. or C4311.) separately.


Dimensions



Specification Information

Step 1.

C4320.



Description


This metal cover finishes the exposed hardware on the end of a riser frame

0	^	+	_	N
	ρ			


Specify 2 end covers per riser frame.

Specify end cover finish to match the corresponding finish of workstation (C4310. or C4311.) to which cover will be attached.

Order workstation (C4310. or C4311.) separately.



Specif	ication Information	
Step 1.		
C4321.		\$90
Step 2.	Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$0



Description

This power strip mounts horizontally on hanger clips and provides single-circuit, 15-amp, 100-volt outlets on a riser frame. It has a lighted on/off switch with circuit breaker and an 8' cord. It is UL listed. Mounting hardware is included.

Power strips have the following number of outlets:

Width-Outlets

24"-4 outlets

30"-5 outlets

36"−6 outlets

48''-8 outlets

60"-10 outlets

72"—12 outlets

Notes


To mount power strip inside riser frame (vertically), specify 36" width.


Specify power strip finish to match the corresponding finish of workstation (C4310. or C4311.) to which it will be attached.

Order workstation (C4310. or C4311.) separately.

For lab modules (C1112., C1114.) and support panels (C1150., C1151.), power strip can be specified as part of the product.

Dimensions





Specification Information


Step 1.

C4322.

Step	2. Width			
24	24" wide			
30	30" wide			
36	36" wide			
48	48" wide			
60	6o" wide			
72	72" wide			

Prices for Steps 1-2.	
C4322. 24	\$387
30	\$416
36	\$446
48	\$460
60	\$519
72	\$549

Step 3.	Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+ \$o




Description

This 2-sided white board mounts between 2 riser frame dividers or a riser frame and a riser frame divider. Both sides are magnetic with white erasable writing surfaces. Includes a black-texture marker tray on one side only. Mounting hardware is included.

Notes

Order workstation (C4310. or C4311.) separately.

Dimensions



Specification Information


Step 1.

C4323.

Step	2. Width			
24	24" wide			
30	30" wide			
36	36" wide			

Prices for Steps 1-2.	
C4323. 24	\$335
30	\$371
36	\$395

\$503



Product Information


Description


This arm mounts to a workstation riser frame and holds a flat panel display monitor. The arm meets VESA standards (Video Electronic Standards Association - Mounting Interface Standards) with a 100 mm or 75 mm mounting pattern. It provides 360° monitor rotation and 60° horizontal and vertical tilt range. Finish is black. Mounting hardware is included. Finish is black.

Notes

Order workstation (C4310. or C4311.) separately.

Dimensions






Specification Information

Step 1.

C4324.



Description

This component hangs from a module, a support panel, or wall strips and holds a 12"-high sink. It includes a front panel and side panels to conceal the bottom of the sink and left and right work surface supports. Attachment hardware is included. Shipped knocked down.

The module accepts the following cutout dimensions:

Depth:


- 24": 181/4"-deep maximum cutout
- 30": $24^{1/4}$ "-deep maximum cutout Width:
- 24": 19⁵/₈"-wide maximum cutout
- 30": 25⁵/₈"-wide maximum cutout
- 48": 435/8"-wide maximum cutout
- 60": 55⁵/₈"-wide maximum cutout

Notes

Sink must be customer supplied and field installed.

Order standard duty work surface (C2142.) with solid surface option (K) separately.

Dimensions



Specification Information

Step 1.

C2512.

Step 2. Depth

24" deep

30 30" deep

Step 3. Width

24" wide

30 30" wide


36 36" wide

48 48" wide

60 60" wide

Prices for Steps 1-3.					
	24	30	36	48	60
C2512. 24	\$510	595	625	654	712
30	\$607	712	750	783	878

C+	Confere Finish	
Step 2	4. Surface Finish	
91	white	+\$o
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o
нм	natural maple	+\$o
76	light brown walnut	+\$o
нх	aged cherry	+\$o



Description

This 24"-deep, 48"-wide structure hangs from 2 24"-wide panels or freestanding modules or from 3 wall strips. It supports drawer storage on the left or right side and a sink unit on the other side. The maximum cutout for the sink is 18^{1} /4" deep × 19^{5} /8" wide. The surface includes a front panel and side panels to conceal the bottom of the sink and left and right work surface supports. Attachment hardware is included. Shipped knocked down.

Notes


Order storage components for open space opposite sink unit separately:

- Drawer bearers (CO228 or CO396)
- Drawer caddy with casters (CO692. or CO693.)
- Gripper package (CO266) and C frame (CO205) for Co/Struc® drawers
- L cart (CO212) and single counter top (CO555) for C frames and drawers when surface is mounted 38" above floor
- Storage case (C5040.2422) and 24"-wide Co/Struc adapter rail (C1610.24)


Sink must be customer supplied and field installed.

Storage unit (C5010.) cannot mount under sink unit.

Order standard duty work surface (C2142.) with solid surface option (K) separately.



	cification Information	
Step	1.	
C251	3.24	\$611
Step	2. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0
91	white	+\$75



Description


This sink support hangs from a module or wall strips. It supports a $5^{1}/2$ "-deep sink and a solid surface or alternative material top. The sink support has a 4"-high laminate front trim and a removable front panel to allow access to plumbing. Attachment hardware is included. Shipped knocked down.

Notes

Order standard duty work surface (C2142.) with solid surface option (K) separately.

Sink must be customer supplied and field installed. A Just SL-ADA-1921-A-GR sink, with drain located in rear center, fits in a 24" \times 24" sink support.

Dimensions



Specification Information

Step 1.

60

C2520.

Step 2. Depth					
24	24" deep				
30	30" deep				

Step 3.	Width
24	24" wide
30	30" wide
48	48" wide

60" wide

Prices for Steps 1-3.				
	24	30	48	60
C2520. 24	\$967	1088	1208	1329
30	\$1027	1147	1329	1450

Sten	4. Surface Finish	
91	white	+\$120
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+ \$o
нм	natural maple	+\$120
76	light brown walnut	+\$120
нх	aged cherry	+\$120
OG	honey maple	+\$180

91

нм

НΧ

OG

white

natural maple

aged cherry


honey maple

+\$60

+\$60

+\$60

+\$90



Product Information

Description

This cabinet hangs below a heavy-duty work surface or process table. It has 2 side panels, a back panel, and a top and bottom panel.


Notes

Order heavy-duty work surface (C2141.), process table (C4313.), standard work table base (C461G.), adjustable lab work table base (C462G., C462C.), work table (C4312.), or workstation (C4310. or C4311.) separately.


Order shelves (C5015.1622) separately; shelves adjust in $1^1\!/_a{}^n$ increments.

Order A-size drawers (CO207) or B-size drawers (CO208) and lockable drawer bearer (C5048.) separately.

Order storage unit door (C5016. or C5017.) separately.



Spe	cification Information	
Step		
C501	0.1622	\$616
Step	2. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+\$ o
LU	soft white	+\$o
WN	warm grey neutral	+\$o
76	light brown walnut	+\$60




Description


This shelf mounts inside a storage unit and adjusts in 1¹/₄" increments.

Votes

Order storage unit (C5010.) separately.



Speci	Specification Information			
Step 1.				
C5015	.1622	\$177		
Sten 2	. Surface Finish			
HF	inner tone light	+\$0		
LT	light tone	+\$0		
LU	soft white	+ \$o		
WN	warm grey neutral	+ \$o		
76	light brown walnut	+\$15		
91	white	+\$15		
нм	natural maple	+\$15		
нх	aged cherry	+\$15		
OG	honey maple	+\$25		



Description

This door attaches to the front of a storage unit. Attachment hardware is included

Notes


Co/Struc® drawers cannot be specified behind doors.

When specifying 19"- or 22"-high lockable door that has drawers above it, specify interior filler (C5009.) to prevent access to cabinet below.

For general information on locks, see Keyed-Alike Information in

Appendices.

Dimensions



Specification Information

Step 1.

C501

Step	2.	Hinge	Location

6. left-hand hinges

7. right-hand hinges

Step 3. Height

1922 19" high

2222 22" high

2522 25" high

Step 4. Lock

N no lock


L keyed differently

K keyed alike

Prices for Steps 1-4.			
	N	L	K
C5016. 1922	\$229	284	284
2222	\$262	318	318
2522	\$306	361	361
C5017. 1922	\$229	284	284
2222	\$262	318	318
2522	\$306	361	361

Step 5	. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o
76	light brown walnut	+\$30
91	white	+\$30
нм	natural maple	+\$30
нх	aged cherry	+\$30
OG	honey maple	+\$45

Step 6	Pull Option/Finish	
NN	none	+\$0
BK	black	+\$ o
NH	brushed nickel	+\$0





Description


This filler mounts inside a storage unit and provides closure to the cabinet below when Co/Struc drawers are installed. It mounts at either 19" or 22" high.

Notes	
Order storage unit (C5010.) separately.	

Filler is used with A-size drawers (CO207) or B-size drawers (CO208).








Description

These drawer bearers attach Co/Struc® A-, B-, or C-size drawers in a storage unit. They have integral stops to prevent drawers from falling when pulled out. Mounting hardware is included. Finish is black. \\ Package contains 1 pair.

Order A-size drawers (CO207), B-size drawers (CO208), or C-size drawers (CO209) separately.

Can be used with locking or non-locking drawers.

Dimensions




Specification Information

Step 1.

C5048.

\$75



Description

This case mounts under a heavy-duty rectangular or corner work surface. It has 2 side panels, a back panel, and a top and bottom panel. Interior components can be height adjusted in $1^1/4^n$ increments.

Notes

Order heavy-duty work surface (C2141.) or heavy-duty corner work surface (C2215. or C2216.) separately.


To support case mounted under rectangular work surface, order adapter rail (C1610.). To support case mounted under corner work surface, order corner storage case rail (C1611.) separately.

Order optional products separately:

- Storage case bottle drawer (C5045.)
- Storage case door (C5052. or C5053.)
- Storage case drawer (C5043. or C5044.)
- Storage case pullout shelf (C5042.)
- Storage case shelf (C5041.)

For top slot of case, specify 16"-deep drawer or shelf to clear support beam under heavy-duty work surface.

Dimensions



Specification Information

Step 1.


C5040.

Step 2.	Depth
2422	23" deep

3022 29" deep

Prices for Steps 1-2.	
C5040. 2422	\$697
3022	\$767

Step 3.	Surface Finish	
HF	inner tone light	+\$ o
LT	light tone	+ \$o
LU	soft white	+ \$o
WN	warm grey neutral	+ \$o
76	light brown walnut	+\$70
91	white	+\$70
нм	natural maple	+\$70
НХ	aged cherry	+\$70
OG	honey maple	+\$105




Description

This shelf mounts inside a storage case and adjusts in $1^{1}/_{4}$ " increments.

Notes

To clear support beam on heavy-duty work surface, shelf must be positioned at least $2^1\!/2^n$ below the surface.

Dimensions



Specification Information


Step 1.

C5041.

Step 2	. Depth			
2422	23" deep			
3022	29" deep			
Duisse for Change of				

Prices for Steps 1-2.	
C5041. 2422	\$200
3022	\$211

Step 3. Surface Finish					
HF	inner tone light	+\$ 0			
LT	light tone	+\$ o			
LU	soft white	+ \$o			
WN	warm grey neutral	+ \$o			
76	light brown walnut	+\$20			
91	white	+\$20			
нм	natural maple	+\$20			
нх	aged cherry	+\$20			
OG	honey maple	+\$30			




Description


This pullout shelf is $21^1/_4$ " deep and mounts inside a storage case. It adjusts in $1^1/_4$ " increments.

Notes

To clear support beam on heavy-duty work surface, shelf must be positioned at least $2^1\!/2^n$ below the surface.



Spec	cification Information			
Step	1.			
C504	C5042.2422			
Step	2. Surface Finish			
HF	inner tone light	+\$o		
LT	light tone	+\$o		
LU	soft white	+\$o		
WN	warm grey neutral	+\$o		
76	light brown walnut	+\$30		
91	white	+\$30		
нм	natural maple	+\$30		
нх	aged cherry	+\$30		
OG	honey maple	+\$45		



Description


This drawer mounts inside a storage case and is available with or without a lock. The 22"-deep drawer cannot be used in the top position of a storage case.

Notes

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions





Specification Information

Step 1.

C504

Step 2. Depth

3. 16" deep

4. 22" deep

Step 3. Height

For 16" deep (3.)

0422 small, 5" high

0822 medium, 7¹/₂" high

For 22" deep (4.)

0422 small, 5" high

0822 medium, 7¹/₂" high

1222 large, 10" high

Step 4. Lock

N no lock


L keyed differently

K keyed alike

Prices for Steps 1-4.			
	N	L	K
C5043. 0422	\$259	428	428
0822	\$324	489	489
C5044. 0422	\$284	374	374
0822	\$340	504	504
1222	\$359	522	522

Step 5.	. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o
76	light brown walnut	+\$40
91	white	+\$40
нм	natural maple	+\$40
нх	aged cherry	+\$40
OG	honey maple	+\$60


Step	6. Pull Option/Finish	
NN	none	+\$
BK	black	+\$
NH	brushed nickel	+\$



Description

This pullout drawer mounts in a storage case and has a partial-height front for easy access to contents. It includes 3 dividers.

Dimensions



Specification Information

Step 1.


C5045.

Step 2. Height						
	5" high					
1222	7 ¹ / ₂ " high					

Prices for Steps 1-2.	
C5045. 0822	\$609
1222	\$649

Step 3. Surface Finish					
HF	inner tone light	+\$o			
LT	light tone	+\$o			
LU	soft white	+\$o			
WN	warm grey neutral	+\$o			
76	light brown walnut	+\$60			
91	white	+\$60			
нм	natural maple	+\$60			
нх	aged cherry	+\$60			
OG	honey maple	+\$90			


Step 4	Pull Option/Finish	
NN	none	+\$o
BK	black	+\$ 0
NH	brushed nickel	+ \$o



Description

This divider fits inside a storage case bottle drawer. Finish is soft white. Package contains 3.

Dimensions



Specification Information

Step 1.

C5049.

Step 2. Height

0824 5¹/₄" high

1224 9¹/₈" high


Prices for Steps 1-2.

C5049. 0824

1224

\$94

\$77



Description

This door attaches to the front of a storage case. Attachment hardware is included.


Notes

Co/Struc® storage case drawers and storage case bottle drawers cannot be specified behind doors.

When specifying 17"- or 20"-high lockable door that has drawers above it, specify storage case shelf (C5041.) to prevent access to cabinet below.

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

C505

Step	2.	H	ing	ge	Lo	C	a	tic	on
			-						

2. left-hand hinges

3. right-hand hinges

Step 3. Height

1722 17" high

2022 20" high

2522 25" high

Step 4. Lock

N no lock


L keyed differently

K keyed alike

Prices for Steps 1-4.			
	N	L	K
C5052. 1722	\$214	268	268
2022	\$231	286	286
2522	\$292	347	347
C5053. 1722	\$214	268	268
2022	\$231	286	286
2522	\$292	347	347

Step 5	. Surface Finish	
91	white	+\$30
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+ \$o
WN	warm grey neutral	+ \$o
нм	natural maple	+\$30
76	light brown walnut	+\$30
нх	aged cherry	+\$30
OG	honey maple	+\$45

Step	6. Pull Option/Finish	
NN	none	+\$0
BK	black	+ \$o
NH	brushed nickel	+\$o




Description


This leg supports a storage case mounted under a heavy-duty work surface. It provides support when the work surface holds heavy loads. The legs adjust in 1/2" increments to allow a work surface height of 34" to 38". Leg finish is medium tone.

Notes

Order storage case (C5040.) separately.

Dimensions






Specification Information

Step 1.

C5050. \$170





Description

This support rail attaches to wall strips or to 24"-wide panels or modules joined by an Action Office® Series 1 or 2 corner connector or a C-series corner connector. It supports a storage case hung from a 48"-wide, 24"-or 30"-deep heavy-duty corner work surface.


Notes

Order storage case (C5040.) separately.





Spec	cification Information	
Step :	1.	
C1611.48		\$230
Step	2. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+\$0
LU	soft white	+\$o
WN	warm grey neutral	+ \$o



Description


This unit hangs on a module, support panel, or wall strips above a work surface. It has sliding glass doors. Shipped knocked down.

Notes

Order optional task light separately:

• Utility task light (G6136.)

Dimensions



Specification Information


Step 1.

C5014.16

Step 2.	. Width
30	30" wide
48	48" wide

Prices for Steps 1-2.	
C5014.16 30	\$1161
48	\$1358

Step 3	3. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o
76	light brown walnut	+\$110
91	white	+\$110
нм	natural maple	+\$110
нх	aged cherry	+\$110
OG	honey maple	+\$160



Description

This 3- or 5-shelf unit hangs from a module, panel, or wall strips. It has $9^{\prime\prime}$ - or $16^{\prime\prime}$ -deep shelves that adjust in $2^{1}/_{2}^{\prime\prime}$ increments. Shelves are flat or slanted at either 11° or 22° angles for gravity-feed dispensing. Shipped knocked down.

Notes

For additional shelf, order add-on shelf (C5025.) separately.

To cover shelf unit with dust cover, order add-on shelf (C5025.) and place in top position of unit.

To enclose back of unit, order back panel (C5026.) separately.

To hang unit on architectural wall and enclose back of unit, order shelf unit wall attachment panel (C5036.) separately.

To enclose front of unit, order shelf unit door (C5027., C5028., C5032., C5037., or C5038.) separately.


To ease cleaning and prevent the stacking of materials on top of shelf unit, order sloped shelf unit top (CO397.) separately.

Order optional task light for 16"-deep unit separately:


- Utility task light (G6136.)
- Energy-efficient task light (G6120.)

To attach task light to slanted shelf, order light bracket (C5051.) separately.





3 Shelves



5 Shelves

Specification Information

Step 1.

C502

Step 2. Shelves/Height

0. 3 shelves and 24" high

1. 5 shelves and 34" high

Step 3. Depth

09 9" deep

16" deep

Step 4. Width

24 24" wide


30 30" wide

36 36" wide48 48" wide

D: 6 6:

Prices for Steps 1-4.				
	24	30	36	48
C5020. 09	\$438	462	472	498
16	\$549	577	610	683
C5021. 09	\$624	644	672	716
16	\$748	800	846	938

Step	5. Side Panel/Shelf Finish	
HF	inner tone light	+\$o
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$o
91	white	+\$60



Description

This unit hangs from 24"-wide panels connected by an Action Office® 2-way 90° connector or wall strips used in a corner position. Shelves adjust in $2^1/2^n$ increments and are flat or slanted at either 11° or 22° angles for gravity-feed dispensing. Shipped knocked down.

Notes


For additional shelf in center section, order 24"-wide add-on shelf (C5025.1624). For additional shelf in side section, order 14"-wide add-on shelf (C5025.1614).


To cover back of center unit, order back panel (C5026.) separately.

Order optional 24"-wide task light for center section separately:


- Utility task light (G6136.24)
- Energy-efficient task light (G6120.24)

To attach task light to slanted center shelf, order light bracket (C5051.) separately.





3 Shelves



5 Shelves

Specification Information

Step 1.

C502

Step 2. Shelves/Height

2.1648 3 shelves and 24" high


3.1648 5 shelves and 34" high

Prices for Steps 1-2.

C5022.1648 \$1742

C5023.1648 \$2340

Step 3	3. Side Panel/Shelf Finish	
91	white	+\$160
HF	inner tone light	+\$ o
LT	light tone	+\$ o
LU	soft white	+\$ o
WN	warm grey neutral	+\$o




Description


This 72"-high shelf unit hangs from 24"- or 48"-wide panels connected by 1 or 2 Action Office® 2-way 90° connectors or hangs from wall strips used in 1 or 2 90° corners. Shelves adjust in $2^1/_2$ " increments and are flat or slanted at either 11° or 22° angles for gravity-feed dispensing. Shipped knocked down.

Notes


For additional shelf, order add-on shelf (C5025.1631) separately.

To cover shelf unit with dust cover, order add-on shelf (C5025.1631) and place in top position of unit.






Specification Information			
Step	1.		
C5024.1648		\$3185	
Step	2. Side Panel/Shelf Finish		
HF	inner tone light	+ \$o	
LT	light tone	+\$0 +\$0	
LU	soft white	+\$0	
WN	warm grey neutral	+\$0	
91	white	+\$300	



Description

This shelf fits inside a shelf unit, diagonal corner shelf unit, or 90° corner shelf unit. It is used as an additional shelf or dust cover.

Dimensions



Specification Information

Step 1.

C5025.

~ .		_		. 1
Ste	กว	- 11	Δr	٦th

09 9" deep

16 16" deep

Step 3. Width

For 9" deep (09)

24 22³/₄" wide

30 28³/₄" wide

36 34³/₄" wide

48 46³/₄" wide

For 16" deep (16)

14 13¹/₂" wide

24 22³/₄" wide

30 28³/₄" wide


31 30¹/₂" wide

36 34³/₄" wide

48 46³/₄" wide

Prices for Steps 1-3.						
	14	24	30	31	36	48
C5025. 09	_	\$88	94	_	98	113
16	\$01	112	115	110	121	130

Step	4. Shelf Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0





Description

This 3"-wide clip fastens to the edge of a metal shelf and holds content labels. Stick-on labels can be applied to the front of the clip or slide-in labels can be inserted from the back. Package contains 36 clear clips.










Description

This panel attaches to a 3- or 5-shelf unit or diagonal corner shelf unit and closes the back. The shelf unit's top and bottom shelves must be flat.

Notes

Order shelf unit (C5020. or C5021.) or diagonal corner shelf unit (C5022. or C5023.) separately.

Dimensions



Specification Information

Step 1.

C5026.

24 24" high

34 34" high

Step 3. Width

24 22³/₄" wide


30 28³/₄" wide

36 34³/₄" wide

48 46³/₄" wide

Prices for Steps 1-3.				
	24	30	36	48
C5026. 24	\$165	179	191	228
34	\$199	218	235	272

Step 2	4. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$ o
WN	warm grey neutral	+\$0



Description

These brackets fasten a task light to the bottom of an overhead cabinet, a flat or slanted shelf in a 3- or 5-shelf unit, or a center shelf in a diagonal corner shelf unit. Package contains 1 pair. Finish is black.

Notes

For 16"-deep shelf, order energy-efficient task light (G6120.) separately.

Dimensions



Step 1.


C5051.


Step 2. Usage

16 for 16"-deep shelf

Prices for Steps 1-2.

C5051. 16 \$11





Description


This 24"-wide door encloses a shelf unit. Attachment hardware is included.

Notes

For 24"-high door, order 3-shelf unit (C5020.) separately. For 34"-high door, order 5-shelf unit (C5021.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

C502

Step 2. Hinge Location

7. left-hand hinges

8. right-hand hinges

Step 3. Height

2424 24" high

3424 34" high

Step 4. Lock

N no lock


L keyed differently

K keyed alike

Prices for Steps 1-4.			
	N	L	K
C5027. 2424	\$340	440	440
3424	\$449	548	548
C5028. 2424	\$340	440	440
3424	\$449	548	548

Step	5. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$ o
LU	soft white	+\$ o
WN	warm grey neutral	+\$ o
76	light brown walnut	+\$40
91	white	+\$40
нм	natural maple	+\$40
нх	aged cherry	+\$40
OG	honey maple	+\$60

Step 6.	Pull Option/Finish	
NN	none	+\$o
BK	black	+\$o
NH	brushed nickel	+\$o



Description

These doors enclose a shelf unit. The doors are available with or without a lock. Attachment hardware is included.


Notes

For 24"-high doors, order 3-shelf unit (C5020.) separately.

For 34"-high doors, order 5-shelf unit (C5021.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

C5032.

Step 2	. Height			
24	24" high			
34	34" high			
Step 3	. Width			
30	30" wide			

Step 4. Lock

36

48

no lock

L keyed differently

36" wide


48" wide

K keyed alike

Prices for Steps 1-4.			
	N	L	K
C5032. 24 30	\$358	458	458
36	\$419	519	519
48	\$514	614	614
34 30	\$464	564	564
36	\$546	646	646
48	\$716	816	816

Step 5	. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0
76	light brown walnut	+\$50
91	white	+\$50
нм	natural maple	+\$50
нх	aged cherry	+\$50
OG	honey maple	+\$75

Step 6.	Pull Option/Finish	
NN	none	+\$0
BK	black	+\$0
NH	brushed nickel	+ \$0




Description


This top attaches to a shelf unit. It is sloped to ease cleaning and prevent the stacking of materials on top of the cabinet. Attachment hardware is included.

Notes

Order shelf unit (C5020. or C5021.) separately.

Dimensions





Specification Information

Step 1.

CO397.

Step 2. Depth

09 10¹/₂" deep

16 17" deep

Step 3. Width

24 24" wide


30 30" wide

36 36" wide

48 48" wide

Prices for Steps 1-3.				
	24	30	36	48
CO397. 09	\$258	287	311	341
16	\$287	307	341	406

Step.	4. Surface Finish	
HF	inner tone light	+\$0
LT	light tone	+ \$o
LU	soft white	+ \$o
WN	warm grey neutral	+\$o




Description

This panel encloses the back of a 3- or 5-shelf unit and allows the shelf unit to hang on an architectural wall.

Notes

Order 3-shelf unit (C5020.) or 5-shelf unit (C5021.) separately.

Dimensions



Specification Information

Step 1.


C5036.

Step 2	. Height			
24	24" high			
34	34" high			
Step 3	. Width			
2/1	24" wide			

24 24" wide 30 30" wide 36 36" wide 48 48" wide

Prices for Steps 1-3.				
	24	30	36	48
C5036. 24	\$184	205	247	288
34	\$199	213	253	301

Step 4.	. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+\$ o
LU	soft white	+\$ o
WN	warm grey neutral	+\$ o
76	light brown walnut	+\$25
91	white	+\$25
HM	natural maple	+\$25
нх	aged cherry	+\$25
OG	honey maple	+\$40



Description

This $15^{1/2}$ "-high shelf with a lip hangs from a standard rail or adapter rail. It can hold up to 2 add-on shelves. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.

Notes

For enclosed storage, order B-style flipper door (A3310., A3312., or A3313.) separately.


For fully enclosed storage, order A-/B-style flipper door back panel (A3390.) and B-style flipper door (A3310., A3312., or A3313.) separately. Order optional add-on shelves of equal depth and width separately:

- Add-on shelf (AO521.1524 or AO521.1548)
- Add-on shelf with lip (C3019.)

Order optional task light separately:

- Utility task light (G6136.)
- Energy-efficient task light (G6120.)

Dimensions



Specification Information

Step 1.

CO521.

Step 2. Depth

13 13" deep

16 16" deep

Step 3. Width

24" wide

30 30" wide


36" wide

42 42" wide

48 48" wide

Prices for Steps 1-3.					
	24	30	36	42	48
C0521. 13	\$410	421	440	450	474
16	\$568	581	610	625	678

Step 4	. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+ \$o



Description

This $15^{1}/2$ "-high shelf with a lip hangs from a panel, a freestanding module, or wall strips. It can hold up to 2 add-on shelves. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.

Notes

For enclosed storage, order B-style flipper door (A3310., A3312., or A3313.) separately.

For fully enclosed storage of 24"- to 48"-wide shelf, order A-/B-style flipper door back panel (A3390.) and B-style flipper door (A3310., A3312., or A3313.) separately.


Order optional add-on shelves of equal depth and width separately:


- Add-on shelf (AO521.1524 or AO521.1548)
- Add-on shelf with lip (C3019.)


For 60"-wide shelf, order 60"-wide task light.

Order optional task light separately:

- Utility task light (G6136.)
- Energy-efficient task light (G6120.)







FRONT

Specification Information

Step 1.

C3011.

Step 2.	Depth
13	13" deep
16	16" deer

Step 3. Width

24 24" wide 30 30" wide


36 36" wide

42" wide 42 48 48" wide

60" wide

Prices for Steps 1-3.						
	24	30	36	42	48	60
C3011. 13	\$233	245	258	269	284	314
16	\$332	345	364	376	414	498

Step 4	. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+\$0
LU	soft white	+ \$o
WN	warm grey neutral	+ \$o



Description

This 16"-deep shelf hangs from a panel, a freestanding module, or wall strips and has a sloped metal top. A filler piece is included for enclosed storage applications. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.

Note:


For enclosed storage, order B-style flipper door (A3310., A3312., or A3313.) separately.


To enclose back of shelf, order A-/B-style flipper door back panel (A3390.) separately.


Order optional task light separately:

- Utility task light (G6136.)
- Energy-efficient task light (G6120.)

Dimensions







Specification Information


Step 1.

C3091.16

Step 2	2. Width				
24	24" wide				
30	30" wide				
42	42" wide				
48	48" wide				

Prices for S	teps 1-2.	
C3091.16	24	\$458
	30	\$480
	42	\$559
	48	\$586

Step 3.	Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o



Description


This shelf fits into a B-style shelf or B-style rail-hanging shelf to add a second shelf. Attachment hardware is included.


Notes

Order support components separately:

- B-style (C3011.) shelf
- B-style (CO521.) rail-hanging shelf

Dimensions





Specification Information

Step 1.

C3019.

Step 2.	Depth
13	13" deep


16 16" deep

Step 3. Width

24	24"	wide


Prices for Steps 1-3.						
	24	30	36	42	48	60
C3019. 13	\$134	141	157	168	181	202
16	\$196	210	228	240	278	349

Step	4. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$o





Description

This 3"-wide clip fastens to the edge of a metal shelf and holds content labels. Stick-on labels can be applied to the front of the clip or slide-in labels can be inserted from the back. Package contains 36 clear clips.









Description

This enclosed storage unit hangs on rails for stationary storage, or it is used with the TR $_3$ cart for mobile storage. The locker has 20 drawer slots in $_3$ " increments and holds interchangeable components.

The locker has the following lock options: no lock, keyed lock, or keyless lock. The keyless lock mounts on the right-hand side and front of the locker. It includes a back-up keyed lock and has 89 potential user codes and 11 potential supervisor codes. The keyless lock is shipped with 6 AA batteries.

Notes


Order rail components separately:


• Standard rail assembly (CO₂82)

Order transport/storage components separately:

- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO₅6₅)
- Drawer (CO207, CO208, or CO209)
- Pullout catheter rack (CO469.)
- TR3 cart (CO342)


For general information on locks, see Keyed-Alike Information in Appendices.





15	Specif	ication Information		13	key number 13	+\$0
16	Step 1.			14	key number 14	+\$0
Step 2_ Lock	CO561			15	key number 15	+ \$o
FF				16	key number 16	+ \$o
FL keyed lock, keyed differently	Step 2.	Lock		17	key number 17	+ \$o
KA keyed lock, keyed alike EL keyless lock, keyed differently EK keyless lock, keyed differently EK keyless lock, keyed differently EX keynumber 22 22 keynumber 23 33 keynumber 23 450 CO561 FF \$16,44 \$51791 26 keynumber 25 450 KA \$1791 27 keynumber 25 450 KA \$1791 27 keynumber 26 450 KA \$1791 27 keynumber 27 450 KA \$1791 27 keynumber 28 450 EK \$2647 29 keynumber 28 450 EK \$2647 29 keynumber 30 8eynumber 30 8eynumber 30 8eynumber 31 8eynumber 31 8eynumber 31 8eynumber 31 8eynumber 31 8eynumber 31 8eynumber 32 450 Step 3. Locker Finish 31 keynumber 31 8eynumber 33 8eynumber 36 8eynumber 36 8eynumber 36 8eynumber 37 450 Step 4. Door Finish 36 keynumber 38 450 A keynumber 39 450 A keynumber 49 450 A keynu	FF	no lock		18	key number 18	+ \$o
EL keyless lock, keyed differently 21 key number 21 450 EK keyless lock, keyed alike 22 key number 32 450 Prices for Steps 1:2. 24 key number 23 450 COS61 FF \$1644 25 key number 26 key number 36 450 KA \$1791 26 key number 27 450 KA \$1791 27 key number 28 450 EK \$2647 28 key number 27 450 EK \$2647 28 key number 27 450 EK \$2647 28 key number 28 450 EK \$2647 29 key number 29 450 Step 3. Locker Finish 31 key number 31 450 HF inner tone light 450 32 key number 33 450 LT light tone 450 33 key number 34 450 WN war grey neutral 450 34 key number 34 450 </th <th>FL</th> <th>keyed lock, keyed differently</th> <th></th> <th>19</th> <th>key number 19</th> <th>+\$0</th>	FL	keyed lock, keyed differently		19	key number 19	+\$0
EK keyless lock, keyed allike 22 key number 22 +\$0 Prices for Steps 1:-2. 24 key number 23 +\$0 CO561 FF \$1644 25 key number 25 +\$0 FL \$1791 26 key number 26 +\$0 KA \$1791 27 key number 28 +\$0 EL \$2647 28 key number 28 +\$0 EK \$2647 29 key number 28 +\$0 FER \$2647 29 key number 29 +\$0 Step 3. Locker Finish 31 key number 30 +\$0 Step 3. Locker Finish 31 key number 30 +\$0 Step 3. Locker Finish 31 key number 32 +\$0 LU soft white +\$0 33 key number 32 +\$0 LU soft white +\$0 33 key number 32 +\$0 LU soft white +\$0 33 key number 32 +\$0 LU soft white +\$0	KA	keyed lock, keyed alike		20	key number 20	+\$0
Sep number 23	EL	keyless lock, keyed differently		21	key number 21	+\$0
Prices for Steps 1-2. 24 key number 24 \$\$0.	EK	keyless lock, keyed alike		22	key number 22	+\$0
\$\frac{1}{4} \begin{align*} \begin				23	key number 23	+\$0
FL	Prices f	or Steps 1-2.		24	key number 24	+\$0
KA \$1791 27 key number 27 +\$0 EL \$2647 28 key number 28 450 EK \$2647 29 key number 30 450 Step 3. Locker Finish 31 key number 31 450 HF inner tone light 450 32 key number 32 450 LU soft white 450 34 key number 32 450 LU soft white 450 34 key number 32 450 LU soft white 450 34 key number 34 450 WN warm grey neutral 450 34 key number 35 450 Step 4. Door Finish 37 key number 36 450 450 Step 4. Door Finish 37 key number 36 450 Bb beachglass 450 38 key number 38 450 Bb beachglass 450 38 key number 38 450 BC pL 450 450 450 450 450	CO561	FF	\$1644	25	key number 25	+ \$o
FL \$2647 28 key number 28 +\$0 EK \$2647 29 key number 29 +\$0 \$20 key number 30 \$450 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50 \$50		FL	\$1791	26	key number 26	+\$0
EK \$2647 29 key number 29 +\$o 30 key number 30 +\$o HF inner tone light +\$o 31 key number 31 +\$o HF inner tone light +\$o 32 key number 32 +\$o LU soft white +\$o 34 key number 33 +\$o LU soft white +\$o 34 key number 33 +\$o WN warm grey neutral +\$o 35 key number 36 +\$o Step 4. Door Finish 37 key number 36 +\$o 35 Abe beachglass +\$o 38 key number 37 +\$o 45 AB beachglass +\$o 38 key number 37 +\$o 45 AB beachglass +\$o 37 key number 38 +\$o 45 AB beachglass +\$o 39 key number 38 +\$o 45 BB beachglass +\$o 40 key number 38 +\$o <td></td> <td>KA</td> <td>\$1791</td> <td>27</td> <td>key number 27</td> <td>+\$0</td>		KA	\$1791	27	key number 27	+\$0
Step 3, Locker Finish		EL	\$2647	28	key number 28	+\$0
Step 3. Locker Finish		EK	\$2647	29	key number 29	+\$0
HF inner tone light +\$o 32 key number 32 +\$o 4\$o LT light tone +\$o 33 key number 34 +\$o 35 key number 34 +\$o LU soft white +\$o 35 key number 36 +\$o 4\$o Step 4. Door Finish 37 key number 36 +\$o 38 key number 36 +\$o 3B beachglass +\$o 38 key number 38 +\$o AC 450 38 key number 38 +\$o BF inner tone light +\$o 39 key number 38 +\$o BF inner tone light +\$o 39 key number 39 +\$o 4\$o LU soft white +\$o 40 key number 40 +\$o 45 41 key number 43 +\$o LU soft white +\$o 43 key number 43 +\$o 45 key number 44 +\$o LU soft white +\$o 44 key n				30	key number 30	+\$0
LT light tone	Step 3.	Locker Finish		31	key number 31	+\$0
LU soft white +\$o 34 key number 34 +\$o 4\$o WN warm grey neutral +\$o 35 key number 35 +\$o 4\$o Step 4. Door Finish 37 key number 37 +\$o 38 key number 37 +\$o 3B beachglass +\$o 38 key number 38 +\$o DA deep blue +\$o 39 key number 39 +\$o HF inner tone light +\$o 40 key number 40 +\$o LU soft white +\$o 41 key number 40 +\$o LU soft white +\$o 42 key number 43 +\$o LU soft white +\$o 43 key number 43 +\$o Kear +\$o 43 key number 43 +\$o Key Number 45o 44 key number 43 +\$o Key Pumber 45 Key number 45 key number 45 +\$o Key Pumber 45 key number 46 +\$o <th< td=""><td>HF</td><td>inner tone light</td><td>+\$o</td><td>32</td><td>key number 32</td><td>+\$0</td></th<>	HF	inner tone light	+\$o	32	key number 32	+\$0
WN warm grey neutral +\$o 35 key number 35 +\$o +\$o Step 4. Door Finish 37 key number 37 +\$o 3B beachglass +\$o 38 key number 39 +\$o 3B beachglass +\$o 39 key number 39 +\$o 4B blue +\$o 40 key number 40 +\$o 4B blue +\$o 41 key number 40 +\$o 4B blue +\$o 41 key number 40 +\$o 4B blue +\$o 41 key number 42 +\$o 4B blue +\$o 43 key number 43 +\$o 4B blue +\$o 43 key number 44 +\$o 4B key number 44 +\$o 45 key number 44 +\$o 4B key number 45 +\$o 45 key number 46 +\$o 5Cep 5. Key Number 46 key number 46 +\$o 45	LT	light tone	+\$o	33	key number 33	+\$0
Step 4. Door Finish 37 key number 36 4\$0 38 key number 37 4\$0 38 key number 38 4\$0 38 key number 39 4\$0 39 key number 39 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0 4\$0	LU	soft white	+\$o	34	key number 34	+\$0
Step 4. Door Finish 37 key number 37 +\$o 3B beachglass +\$o 38 key number 38 +\$o DA deep blue +\$o 39 key number 39 +\$o HF inner tone light +\$o 40 key number 40 +\$o LT light tone +\$o 41 key number 41 +\$o LU off white +\$o 42 key number 42 +\$o TR clear +\$o 43 key number 43 +\$o WN warm grey neutral +\$o 44 key number 43 +\$o WN warm grey neutral +\$o 44 key number 44 +\$o VSO No 44 key number 45 +\$o Step 5. Key Number 46 key number 46 +\$o For keyed lock, keyed alike (KA) or keyless lock, keyed alike (EK) 47 key number 47 +\$o 10 key number 01 +\$o 48 key number 49 +\$o 10 <td>WN</td> <td>warm grey neutral</td> <td>+\$o</td> <td>35</td> <td>key number 35</td> <td>+\$0</td>	WN	warm grey neutral	+\$o	35	key number 35	+\$0
3B beachglass +\$o 38 key number 38 +\$o DA deep blue +\$o 39 key number 39 +\$o 4\$o key number 40 +\$o +\$o 4\$o key number 41 +\$o 4\$o key number 42 +\$o 4\$o key number 42 +\$o 4\$o key number 42 +\$o 4\$o key number 43 +\$o 4\$o key number 43 +\$o 4\$o key number 44 +\$o 4\$o key number 44 +\$o \$o \$o M*o				36	key number 36	+\$0
DA deep blue +\$o 39 key number 39 +\$o 4\$o HF inner tone light +\$o 40 key number 40 +\$o 4\$o LT light tone +\$o 41 key number 41 +\$o 4\$o LU soft white +\$o 42 key number 42 +\$o 4\$o TR clear +\$o 43 key number 43 +\$o 4\$o WN warm grey neutral +\$o 44 key number 43 +\$o 4\$o Step 5. Key Number 45 key number 45 key number 45 +\$o 4\$o key number 45 +\$o Step 5. Key Number 46 key number 46 *\$o 4\$o key number 45 +\$o 4\$o key number 47 *\$o 4\$o 4\$o key number 47 *\$o 4\$o	Step 4.	Door Finish		37	key number 37	+ \$o
HF inner tone light +\$o 40 key number 4o +\$o LT light tone +\$o 41 key number 41 +\$o LU soft white +\$o 42 key number 42 +\$o +\$o 43 key number 43 +\$o +\$o 40 key number 44 +\$o 40 key number 44 +\$o 40 key number 45 +\$o 40 key number 45 +\$o 40 key number 46 +\$o 40 key number 46 +\$o 40 key number 47 +\$o 40 key number 47 +\$o 40 key number 48 +\$o 40 key number 49 +\$o 40 key number 49 +\$o 40 key number 50 45 45 45 45 45 45 45 45 45 45 45 45 45	3B	beachglass	+ \$0	38	key number 38	+ \$o
LT light tone +\$o 41 key number 41 +\$o 42 key number 42 +\$o 45o 7c 4clear +\$o 43 key number 43 +\$o 45o 44 key number 44 +\$o 45o 44 key number 44 +\$o 45o key number 45 +\$o 46 key number 45 +\$o 46 key number 46 +\$o 45o 46 key number 47 +\$o 45o 46 key number 47 +\$o 45o 45o <t< th=""><th>DA</th><th>deep blue</th><th>+\$0</th><th>39</th><th>key number 39</th><th>+\$0</th></t<>	DA	deep blue	+\$ 0	39	key number 39	+\$0
LU soft white +\$o 42 key number 42 +\$o 45 TR clear +\$o 43 key number 43 +\$o 4\$o WN warm grey neutral +\$o 44 key number 44 +\$o Very Number 45 key number 45 +\$o Step 5. Key Number 46 key number 45 +\$o For keyed alike (KA) or keyless lock, keyed alike (EK) 47 key number 46 +\$o For keyed alike (KA) or keyless lock, keyed alike (EK) 47 key number 46 +\$o For keyed alike (KA) or keyless lock, keyed alike (EK) 47 key number 47 +\$o To key number 01 +\$o 48 key number 47 +\$o 49 key number 48 +\$o To key number 02 +\$o 49 key number 50 key number 50 +\$o 45o	HF	inner tone light	+\$o	40	key number 40	+\$0
TR clear +\$o 43 key number 43 +\$o 45 WN warm grey neutral +\$o 44 key number 44 +\$o Step 5. Key Number 45 key number 45 +\$o Step 5. Key Number 46 key number 46 +\$o For keyed lock, keyed alike (KA) or keyless lock, keyed alike (EK) 47 key number 47 +\$o 01 key number 01 +\$o 48 key number 48 +\$o 02 key number 02 +\$o 49 key number 49 +\$o 03 key number 03 +\$o 50 key number 50 +\$o 04 key number 04 +\$o 51 key number 51 +\$o 05 key number 05 +\$o 52 key number 52 +\$o 06 key number 07 +\$o 54 key number 53 +\$o 07 key number 08 +\$o 55 key number 55 key number 56 +\$o 08 key number 10 +\$o 56	LT	light tone	+\$o	41	key number 41	+\$0
WN warm grey neutral +\$o 44 key number 44 +\$o 45 key number 45 +\$o 45 key number 45 +\$o 46 key number 46 +\$o 48 key number 47 +\$o 48 key number 47 +\$o 48 key number 48 +\$o 20 key number 02 4\$o 49 key number 49 +\$o 4\$o 4\$	LU	soft white	+\$ 0	42	key number 42	+ \$o
Step 5. Key Number 45 key number 45 +\$o For keyed lock, keyed alike (KA) or keyless lock, keyed alike (EK) 47 key number 47 +\$o 01 key number o1 +\$o 48 key number 48 +\$o 02 key number o2 +\$o 49 key number 49 +\$o 03 key number o3 +\$o 50 key number 50 +\$o 04 key number o4 +\$o 51 key number 51 +\$o 05 key number o5 +\$o 52 key number 52 +\$o 06 key number o6 +\$o 53 key number 53 +\$o 07 key number o7 +\$o 54 key number 54 +\$o 08 key number o8 +\$o 55 key number 55 +\$o 09 key number o9 +\$o 56 key number 57 +\$o 10 key number o 57 key number 57 +\$o 11 key number o 58 key number 58 +\$o	TR	clear	+\$ 0	43	key number 43	+ \$o
Step 5. Key Number 46 key number 46 +\$0 For keyed lock, keyed alike (KA) or keyless lock, keyed alike (EK) 47 key number 47 +\$0 01 key number o1 +\$0 48 key number 48 +\$0 02 key number o2 +\$0 49 key number 49 +\$0 03 key number o3 +\$0 50 key number 50 +\$0 04 key number o4 +\$0 51 key number 51 +\$0 05 key number o5 +\$0 52 key number 52 +\$0 06 key number o6 +\$0 53 key number 53 +\$0 07 key number o7 +\$0 54 key number 54 +\$0 08 key number o8 +\$0 55 key number 55 +\$0 09 key number o9 +\$0 56 key number 56 +\$0 10 key number o +\$0 57 key number 57 +\$0 11 key number o 58 key number 58 +\$0	WN	warm grey neutral	+\$o	44	key number 44	+\$0
For keyed lock, keyed alike (KA) or keyless lock, keyed alike (EK) 47 key number 47 +\$0 01 key number o1 +\$0 48 key number 48 +\$0 02 key number o2 +\$0 49 key number 49 +\$0 03 key number o3 +\$0 50 key number 50 +\$0 04 key number o4 +\$0 51 key number 51 +\$0 05 key number o5 +\$0 52 key number 52 +\$0 06 key number o6 +\$0 53 key number 53 +\$0 07 key number o7 +\$0 54 key number 54 +\$0 08 key number o8 +\$0 55 key number 55 +\$0 09 key number o9 +\$0 56 key number 56 +\$0 10 key number o +\$0 57 key number 57 +\$0 11 key number o 58 key number 58 +\$0				45	key number 45	+\$0
01 key number o1 +\$0 48 key number 48 +\$0 02 key number o2 +\$0 49 key number 49 +\$0 03 key number o3 +\$0 50 key number 50 +\$0 04 key number o4 +\$0 51 key number 51 +\$0 05 key number o5 +\$0 52 key number 52 +\$0 06 key number o6 +\$0 53 key number 53 +\$0 07 key number o7 +\$0 54 key number 54 +\$0 08 key number o8 +\$0 55 key number 55 +\$0 09 key number o9 +\$0 56 key number 56 +\$0 10 key number o +\$0 57 key number 57 +\$0 11 key number o 58 key number 58 +\$0	Step 5.	Key Number		46	key number 46	+\$0
02 key number o2 +\$o 49 key number 49 +\$o 03 key number o3 +\$o 50 key number 50 +\$o 04 key number o4 +\$o 51 key number 51 +\$o 05 key number o5 +\$o 52 key number 52 +\$o 06 key number o6 +\$o 53 key number 53 +\$o 07 key number o7 +\$o 54 key number 54 +\$o 08 key number o8 +\$o 55 key number 55 +\$o 09 key number o9 +\$o 56 key number 56 +\$o 10 key number o +\$o 57 key number 57 +\$o 11 key number o 58 key number 58 +\$o	For key	red lock, keyed alike (KA) or keyless lock, keyed alike (EK)		47	key number 47	+\$0
03 key number 03 +\$0 50 key number 50 +\$0 04 key number 04 +\$0 51 key number 51 +\$0 05 key number 05 +\$0 52 key number 52 +\$0 06 key number 06 +\$0 53 key number 53 +\$0 07 key number 07 +\$0 54 key number 54 +\$0 08 key number 08 +\$0 55 key number 55 +\$0 09 key number 09 +\$0 56 key number 56 +\$0 10 key number 10 +\$0 57 key number 57 +\$0 11 key number 11 +\$0 58 key number 58 +\$0	01	key number 01	+\$ o	48	key number 48	+\$0
04 key number 04 +\$0 51 key number 51 +\$0 05 key number 05 +\$0 52 key number 52 +\$0 06 key number 06 +\$0 53 key number 53 +\$0 07 key number 07 +\$0 54 key number 54 +\$0 08 key number 08 +\$0 55 key number 55 +\$0 09 key number 09 +\$0 56 key number 56 +\$0 10 key number 10 +\$0 57 key number 57 +\$0 11 key number 11 +\$0 58 key number 58 +\$0	02	key number 02	+\$o	49	key number 49	+\$0
05 key number 05 +\$0 52 key number 52 +\$0 06 key number 06 +\$0 53 key number 53 +\$0 07 key number 07 +\$0 54 key number 54 +\$0 08 key number 08 +\$0 55 key number 55 +\$0 09 key number 09 +\$0 56 key number 56 +\$0 10 key number 10 +\$0 57 key number 57 +\$0 11 key number 11 +\$0 58 key number 58 +\$0	03	key number 03	+\$ o	50	key number 50	+\$0
06 key number o6 +\$0 53 key number 53 +\$0 07 key number o7 +\$0 54 key number 54 +\$0 08 key number o8 +\$0 55 key number 55 +\$0 09 key number o9 +\$0 56 key number 56 +\$0 10 key number 10 +\$0 57 key number 57 +\$0 11 key number 11 +\$0 58 key number 58 +\$0	04	key number 04	+\$o	51	key number 51	+\$0
07 key number o7 +\$0 54 key number 54 +\$0 08 key number o8 +\$0 55 key number 55 +\$0 09 key number o9 +\$0 56 key number 56 +\$0 10 key number 10 +\$0 57 key number 57 +\$0 11 key number 11 +\$0 58 key number 58 +\$0	05	key number 05	+\$o	52	key number 52	+\$0
08 key number 08 +\$0 55 key number 55 +\$0 09 key number 09 +\$0 56 key number 56 +\$0 10 key number 10 +\$0 57 key number 57 +\$0 11 key number 11 +\$0 58 key number 58 +\$0	06	key number o6	+\$o	53	key number 53	+\$0
09 key number 09 +\$0 56 key number 56 +\$0 10 key number 10 +\$0 57 key number 57 +\$0 11 key number 11 +\$0 58 key number 58 +\$0	07	key number 07	+\$o	54	key number 54	+\$0
10 key number 10 +\$0 57 key number 57 +\$0 11 key number 11 +\$0 58 key number 58 +\$0	08	key number o8	+\$o	55	key number 55	+\$0
11 key number 11 +\$0 58 key number 58 +\$0	09	key number o9	+\$o	56	key number 56	+\$0
	10	key number 10	+\$o	57	key number 57	+\$0
12 key number 12 +\$0 59 key number 59 +\$0	11	key number 11	+\$o	58	key number 58	+\$ o
	12	key number 12	+\$0	59	key number 59	+\$0

60	key number 60	+\$o
61	key number 61	+\$o
62	key number 62	+ \$0
63	key number 63	+\$o
64	key number 64	+\$o
65	key number 65	+\$o
66	key number 66	+\$o
67	key number 67	+\$o
68	key number 68	+\$o
69	key number 69	+\$o
70	key number 70	+\$o
71	key number 71	+\$o
72	key number 72	+\$o
73	key number 73	+\$ 0
74	key number 74	+\$ 0
75	key number 75	+\$ 0
76	key number 76	+\$0
77	key number 77	+\$0
78	key number 78	+\$0
79	key number 79	+\$ 0
80	key number 80	+\$0
81	key number 81	+\$0
82	key number 82	+\$0
83	key number 83	+\$0
84	key number 84	+\$0
85	key number 85	+\$0
86	key number 86	+\$0
87	key number 87	+\$0
88	key number 88	+\$0
89	key number 89	+\$0
90	key number 90	+\$0
91	key number 91	+\$0
92	key number 92	+\$0
93	key number 93	+\$0
94	key number 94	+\$0
95	key number 95	+\$0
96	key number 96	+\$0
97	key number 97	+\$0
98	key number 98	+\$0
99	key number 99	+\$0



Description

This enclosed storage unit hangs on rails for stationary storage, or it is used with the TR3 cart for mobile storage. The locker has 20 drawer slots in 3'' increments and holds interchangeable components. Door is hinged at the right.

Available with no lock, keyless, or a keyed lock.

Notes

Door frame color matches locker color.

Order rail components separately:


• Standard rail assembly (CO282)


Order transport/storage components separately:

- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO₅6₅)
- Drawer (CO207, CO208, or CO209)
- Pullout catheter rack (CO469.)
- TR3 cart (CO342)

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions





Specification Information

Step 1.

CO570

Step 2	. Lock
FF	no lock
FL	keyed lock, keyed differently
KA	keyed lock, keyed alike
KP	keyless lock


or Steps 1-2.	
FF	\$2122
FL	\$2288
KA	\$2288
KP	\$2645
	FF FL KA

Step	3. Locker Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0

Step 4	. Door Insert Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
TR	clear	+\$o
WN	warm grey neutral	+\$o

Step	5. Key Number	
For k	eyed lock, keyed alike (KA)	
01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number o6	+\$0
07	key number 07	+\$0
80	key number o8	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0
15	key number 15	+\$0
16	key number 16	+\$ 0

17	key number 17	+\$o	64	key number 64	+\$o
18	key number 18	+\$o	65	key number 65	+\$o
19	key number 19	+\$o	66	key number 66	+\$o
20	key number 20	+\$o	67	key number 67	+\$o
21	key number 21	+\$o	68	key number 68	+\$o
22	key number 22	+\$o	69	key number 69	+\$o
23	key number 23	+\$o	70	key number 70	+\$o
24	key number 24	+\$o	71	key number 71	+\$o
25	key number 25	+\$o	72	key number 72	+\$o
26	key number 26	+\$o	73	key number 73	+\$0
27	key number 27	+\$o	74	key number 74	+\$0
28	key number 28	+\$o	75	key number 75	+\$0
29	key number 29	+\$o	76	key number 76	+\$0
30	key number 30	+\$o	77	key number 77	+\$0
31	key number 31	+\$o	78	key number 78	+\$0
32	key number 32	+\$o	79	key number 79	+\$0
33	key number 33	+\$o	80	key number 80	+\$0
34	key number 34	+\$o	81	key number 81	+\$0
35	key number 35	+\$ 0	82	key number 82	+\$ 0
36	key number 36	+\$ o	83	key number 83	+\$0
37	key number 37	+\$o	84	key number 84	+\$0
38	key number 38	+\$o	85	key number 85	+\$0
39	key number 39	+\$o	86	key number 86	+\$0
40	key number 40	+\$o	87	key number 87	+\$0
41	key number 41	+\$o	88	key number 88	+\$0
42	key number 42	+\$o	89	key number 89	+\$0
43	key number 43	+\$o	90	key number 90	+\$0
44	key number 44	+\$o	91	key number 91	+\$o
45	key number 45	+\$o	92	key number 92	+\$o
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0 +\$0	95	key number 95	+\$0 +\$0
49 50	key number 49	+\$0 +\$0	96 07	key number 96 key number 97	+\$0 +\$0
50 51	key number 50	+\$0 +\$0	97 98	key number 98	+\$0
52	key number 51 key number 52	+\$0 +\$0	99	key number 99	+\$0
53	key number 53	+\$0	"	key number 99	+40
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			
62	key number 62	+\$0			
63	key number 63	+\$0			
	,	. 40			



Description

This storage unit hangs on rails for stationary storage, or it is used with the TR $_3$ cart for mobile storage. The locker has 20 drawer slots in $_3$ " increments and holds interchangeable components.

Notes


Order rail components separately:


• Standard rail assembly (CO₂82)

Order transport/storage components separately:


- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO₅6₅)
- Drawer (CO207, CO208, or CO209)
- Pullout catheter rack (CO469.)
- TR3 cart (CO342)

Order optional locker tambour door (CO219) separately.





Specification Information					
Step 1	Step 1.				
CO562FF		\$1450			
Step 2	. Surface Finish				
HF	inner tone light	+ \$o			
LT	light tone	+ \$o			
LU	soft white	+\$0			
WN	warm grey neutral	+\$0			



Description

This enclosed, mobile storage unit holds interchangeable components to store and transport materials. The locker has 20 drawer slots in 3" increments and a counterweighted base with 5"-diameter casters (2 swivel and 2 swivel with brake). The tambour door slides vertically for access to locker materials. The locker does not include a handle; it is pushed or pulled by grasping the sides. It cannot be linked for multiple towing and cannot be transported in a truck.


The locker has the following lock options: no lock, keyed lock, or keyless lock. The keyless lock mounts on the right-hand side and front of the locker. It includes a back-up keyed lock and has 89 potential user codes and 11 potential supervisor codes. The keyless lock is shipped with 6 AA batteries.

Notes

Order internal storage components separately:

- A-, B-, or C-size drawers (CO207, CO208, or CO209)
- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO₅6₅)
- Pullout catheter rack (CO469.)

For general information on locks, see Keyed-Alike Information in Appendices.




Specif	fication Information		13	13 key number 13
Step 1.			14	, ,
CO564			15	
			16	, -
Step 2.	. Lock		17	
FF	no lock		18	· ·
FL	keyed lock, keyed differently		19	
KA	keyed lock, keyed alike		20	, ,
EL	keyless lock, keyed differently		21	•
EK	keyless lock, keyed alike		22	,
			23	,
	for Steps 1-2.		24	,
CO564		\$2295	25	, -
	FL	\$2407	26	key number 26
	KA	\$2407	27	27 key number 27
	EL	\$3300	28	28 key number 28
	EK	\$3300	29	29 key number 29
			30	30 key number 30
Step 3.	. Locker Finish		31	
HF	inner tone light	+\$o	32	, ,
LT	light tone	+\$0	33	, ,
LU	soft white	+\$0	34	, , , , , , , , , , , , , , , , , , , ,
WN		+\$0	35	
AA IA	warm grey neutral	+30		,
~.			36	, -
	. Door Finish		37	, .,
3B	beachglass	+\$0	38	, ,
DA	deep blue	+\$o	39	39 key number 39
HF	inner tone light	+\$o	40	40 key number 40
LT	light tone	+\$ 0	41	41 key number 41
LU	soft white	+\$o	42	42 key number 42
TR	clear	+\$0	43	key number 43
WN	warm grey neutral	+\$o	44	
	3 1, 1 11 11		45	,
Sten E	. Key Number		46	,
	yed lock, keyed alike (KA) or keyless lock, keyed alike (EK)		47	,
01	key number o1		48	, .,
	•	+\$o		,
02	key number o2	+\$o	49	, ,,
03	key number 03	+\$0	50	
04	key number 04	+\$o	51	
05	key number 05	+\$o	52	52 key number 52
06	key number o6	+\$ o	53	key number 53
07	key number 07	+\$o	54	54 key number 54
08	key number o8	+\$0	55	55 key number 55
09	key number 09	+\$0	56	
10	key number 10	+\$0	57	
11	key number 11	+\$0	58	
12	key number 12	+\$0	59	
14	KCy HUHIDEI 12	τ⊅∪	37	Ney Hulling 59

60	key number 60	+ \$o
61	key number 61	+\$0
62	key number 62	+ \$o
63	key number 63	+\$0
64	key number 64	+\$0
65	key number 65	+\$0
66	key number 66	+\$0
67	key number 67	+\$0
68	key number 68	+\$0
69	key number 69	+\$0
70	key number 70	+\$0
71	key number 71	+\$0
72	key number 72	+\$0
73	key number 73	+\$0
74	key number 74	+\$0
75	key number 75	+\$ 0
76	key number 76	+\$ 0
77	key number 77	+\$ 0
78	key number 78	+\$0
79	key number 79	+\$0
80	key number 80	+\$0
81	key number 81	+\$0
82	key number 82	+\$0
83	key number 83	+\$0
84	key number 84	+\$0
85	key number 85	+\$0
86	key number 86	+\$ 0
87	key number 87	+\$ 0
88	key number 88	+\$ 0
89	key number 89	+\$ 0
90	key number 90	+\$ 0
91	key number 91	+\$ 0
92	key number 92	+\$0
93	key number 93	+\$0
94	key number 94	+\$0
95	key number 95	+\$0
96	key number 96	+\$0
97	key number 97	+\$0
98	key number 98	+\$0
99	key number 99	+\$0

Locker on Wheel Base, with Hinged Door

CO574



Product Information

Description

This enclosed, mobile storage unit holds interchangeable components to store and transport materials. The locker has 20 drawer slots in 3" increments and a counterweighted base with 5"-diameter casters (2 swivel and 2 swivel with brake). The door is hinged at the right and swings open for access to locker materials. The locker does not include a handle; it is pushed or pulled by grasping the sides. It cannot be linked for multiple towing and cannot be transported in a truck.

Available with no lock, keyless, or a keyed lock.

Notes


Door frame color matches locker color.

Order internal storage components separately:

- A-, B-, or C-size drawers (CO207, CO208, or CO209)
- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO₅6₅)
- Pullout catheter rack (CO469.)

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CO574


Step	2. Lock	
FF	no lock	
FL	keyed lock, keyed differently	
KA	keyed lock, keyed alike	
KP	keyless lock	
Price	s for Steps 1-2.	
CO57	74 FF	\$2906
	FL	\$3070
	KA	\$3070
	KP	\$3429

Step 3	Step 3. Locker Finish			
HF	inner tone light	+\$o		
LT	light tone	+\$o		
LU	soft white	+\$o		
WN	warm grey neutral	+\$o		

Step 4	. Door Finish	
HF	inner tone light	+\$o
LT	light tone	+ \$o
LU	soft white	+\$0
TR	clear	+\$o
WN	warm grey neutral	+\$o

Step 5	. Key Number	
For ke	yed lock, keyed alike (KA)	
01	key number 01	+\$o
02	key number 02	+\$o
03	key number 03	+\$o
04	key number 04	+\$o
05	key number 05	+\$o
06	key number o6	+\$o
07	key number 07	+\$o
08	key number o8	+\$o
09	key number 09	+\$o
10	key number 10	+\$o
11	key number 11	+\$o
12	key number 12	+\$o
13	key number 13	+\$o
14	key number 14	+\$o
15	key number 15	+\$0
16	key number 16	+\$0

17	key number 17	+ \$0	64	key number 64	+ \$0
18	key number 18	+\$ 0	65	key number 65	+\$o
19	key number 19	+\$ 0	66	key number 66	+\$o
20	key number 20	+ \$0	67	key number 67	+\$o
21	key number 21	+\$o	68	key number 68	+\$o
22	key number 22	+\$0	69	key number 69	+\$o
23	key number 23	+ \$0	70	key number 70	+\$o
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36		+\$0 +\$0		•	+\$0
	key number 36	+\$0 +\$0	83	key number 83 key number 84	+\$0 +\$0
37 38	key number 37	+\$0 +\$0	84 85	key number 85	+\$0 +\$0
	key number 38		85	-	
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$o	89	key number 89	+\$o
43	key number 43	+\$o	90	key number 90	+\$o
44	key number 44	+\$o	91	key number 91	+\$o
45	key number 45	+\$o	92	key number 92	+\$o
46	key number 46	+\$o	93	key number 93	+\$o
47	key number 47	+\$o	94	key number 94	+\$o
48	key number 48	+\$o	95	key number 95	+\$o
49	key number 49	+\$o	96	key number 96	+\$0
50	key number 50	+\$ 0	97	key number 97	+\$o
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+ \$0			
57	key number 57	+ \$o			
58	key number 58	+ \$o			
59	key number 59	+\$o			
60	key number 60	+\$ 0			
61	key number 61	+\$o			
62	key number 62	+\$o			
63	key number 63	+\$o			



Description

This mobile storage unit holds interchangeable components to store and transport materials. The locker has 20 drawer slots in 3" increments and a counterweighted base with 5"-diameter casters (2 swivel and 2 swivel with brake). The locker does not include a handle; it is pushed or pulled by grasping the sides. It cannot be linked for multiple towing and cannot be transported in a truck.


Notes

Order internal storage components separately:


- A-, B-, or C-size drawers (CO207, CO208, or CO209)
- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO565)
- Pullout catheter rack (CO469.)

Order chest tube rack (CO466) separately.

Order optional locker tambour door (CO219) and locker lock kit (CO277) separately.




Spec	ification Information	
Step 1	1.	
CO56	3FF	\$2269
Step 2	2. Surface Finish	
HF	inner tone light	+\$ o
LT	light tone	+\$ o
LU	soft white	+\$ o
WN	warm grey neutral	+\$0




Description

This vertical door provides closure on a locker. It opens by sliding into the locker bottom.

Dimensions



Specification Information Step 1.		
CO21	9FF	\$221
Step 2	2. Surface Finish	
3B	beachglass	+\$o
DA	deep blue	+\$o
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
TR	clear	+\$o
WN	warm grey neutral	+\$o




Description

This door provides closure on a locker. It is hinged at the right.

Attachment hardware included.

Available with no lock, keyless, or a keyed lock.

Dimensions



Specification Information

Step 1.

CO290

14

15

16

key number 14

key number 15

key number 16


Step 2	2. Lock	
FF	no lock	
FL	keyed lock, keyed differently	
KA	keyed lock, keyed alike	
KP	keyless lock	
Prices	s for Steps 1-2.	
CO29	0 FF	\$658
	FL	\$827
	KA	\$827
	КР	\$1180
Step 3	3. Door Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
WN	warm grey neutral	+\$0
Step	4. Insert Finish	
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TR	clear	+\$0
WN	warm grey neutral	+\$0
Step	5. Key Number	
For ke	eyed lock, keyed alike (KA)	
01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number o6	+\$0
07	key number 07	+\$0
80	key number o8	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
		and the second s

+\$o

+\$o

+\$o

17	key number 17	+\$o	64	key number 64	+\$o
18	key number 18	+\$o	65	key number 65	+\$o
19	key number 19	+\$o	66	key number 66	+\$o
20	key number 20	+ \$0	67	key number 67	+\$0
21	key number 21	+\$o	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
		+\$o			+\$0
23	key number 23		70	key number 70	
24	key number 24	+\$o	71	key number 71	+\$0
25	key number 25	+\$o	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+ \$0	74	key number 74	+\$ 0
28	key number 28	+\$o	75	key number 75	+\$o
29	key number 29	+\$o	76	key number 76	+\$0
30	key number 30	+\$o	77	key number 77	+\$o
31	key number 31	+\$o	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35		+\$0	82	key number 82	+\$0
	key number 35				
36	key number 36	+\$o	83	key number 83	+\$0
37	key number 37	+\$o	84	key number 84	+\$0
38	key number 38	+\$o	85	key number 85	+\$0
39	key number 39	+ \$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$o	88	key number 88	+\$ 0
42	key number 42	+\$o	89	key number 89	+\$0
43	key number 43	+\$o	90	key number 90	+\$ 0
44	key number 44	+\$o	91	key number 91	+\$o
45	key number 45	+\$o	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$o	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
		+\$o			+\$0
50	key number 50		97	key number 97	
51	key number 51	+\$o	98	key number 98	+\$0
52	key number 52	+\$o	99	key number 99	+\$0
53	key number 53	+\$o			
54	key number 54	+ \$0			
55	key number 55	+ \$0			
56	key number 56	+\$o			
57	key number 57	+ \$0			
58	key number 58	+\$o			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			
62	key number 62	+\$0			
63	key number 63	+\$0 +\$0			
65	key number 03	-φU			



Description


This kit includes a lock and converts a locker with tambour door to a lockable storage unit.

Notes

For general information on locks, see Keyed-Alike Information in Appendices.

Kit must be field installed.

Dimensions



Spe	Specification Information			
Step	1.			
CO27	77			
Step	2. Lock			
FF	keyed differently			
KA	keyed alike			
Price	es for Steps 1-2.			
	77 FF	\$154		
CO2,	KA	\$154		
	NA.	V±34		
Step	3. Surface Finish			
HF	inner tone light	+\$0		
LT	light tone	+\$ o		
LU	soft white	+\$ o		
WN	warm grey neutral	+\$0		
Step	4. Key Number			
For k	keyed alike (KA)			
01	key number 01	+\$0		
02	key number 02	+\$0		
03	key number 03	+\$0		
04	key number 04	+\$0		
05	key number 05	+\$0		
06	key number o6	+\$ o		
07	key number 07	+\$ o		
80	key number o8	+\$ o		
09	key number 09	+\$ o		
10	key number 10	+\$ o		
11	key number 11	+\$ o		
12	key number 12	+\$ o		
13	key number 13	+\$ o		
14	key number 14	+\$ o		
15	key number 15	+\$0		
16	key number 16	+\$ o		
17	key number 17	+\$0		
18	key number 18	+\$0		
19	key number 19	+\$0		
20	key number 20	+\$0		
21	key number 21	+\$0		
22	key number 22	+\$ o		

23

24

25

26

27

key number 23

key number 24

key number 25

key number 26

key number 27

+\$o


+\$o

+\$o

+\$o

+\$o

28	key number 28	+\$ 0	75	key number 75	+\$ 0
29	key number 29	+\$ o	76	key number 76	+\$o
30	key number 30	+\$0	77	key number 77	+\$o
31	key number 31	+\$o	78	key number 78	+\$o
32	key number 32	+\$o	79	key number 79	+ \$ 0
33	-	+ \$o	80	key number 80	+\$o
34		+\$0	81	key number 81	+\$o
35		+\$0	82	key number 82	+\$0
36		+\$0	83	key number 83	+\$0
37		+\$0	84	key number 84	+\$0
38		+\$0	85	key number 85	+\$0
39		+\$0	86	key number 86	+\$0
40		+\$0	87		+\$0 +\$0
	•			key number 87	
41		+\$o	88	key number 88	+\$o
42		+\$o	89	key number 89	+\$o
43		+\$o	90	key number 90	+\$o
44		+\$0	91	key number 91	+\$0
45	,	+ \$o	92	key number 92	+\$o
46	46 key number 46	+ \$0	93	key number 93	+\$o
47	key number 47	+\$ 0	94	key number 94	+\$ 0
48	key number 48	+\$ 0	95	key number 95	+\$ 0
49	key number 49	+\$o	96	key number 96	+\$o
50	key number 50	+\$ o	97	key number 97	+\$ o
51	key number 51	+\$ o	98	key number 98	+\$o
52	key number 52	+\$0	99	key number 99	+ \$o
53	key number 53	+\$o			
54		+\$o			
55		+ \$o			
56		+ \$o			
57		+ \$o			
58		+ \$o			
59		+\$0			
60		+\$0			
61	•	+\$0			
62	•	+\$0			
63	•	+\$0			
64		+\$0			
	•	+\$0			
65					
66	•	+\$o			
67		+\$o			
68	•	+\$o			
69		+\$o			
70		+\$o			
71	•	+\$0			
72	•	+ \$o			
73		+ \$0			
74	74 key number 74	+\$o			



Description


This keyless lock converts a locker with tambour door to a lockable storage unit. It mounts on the right-hand side and front of the locker. The keyless lock includes a back-up keyed lock and has 89 potential user codes and 11 potential supervisor codes. The keyless lock is shipped with 6 AA batteries.


Notes

For general information on locks, see Keyed-Alike Information in Appendices.

Kit must be field installed.

Dimensions





Specification Information

Step 1.

CO279E


Step	2. Lock
L	keyed differently
K	keyed alike
Price	os for Stans 1-2

\$1049
\$1049

Step 3	. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$o
LU	soft white	+\$o
WN	warm grey neutral	+\$o

Step.	4. Key Number	
For ke	eyed alike (K)	
01	key number 01	+\$o
02	key number 02	+\$o
03	key number 03	+\$o
04	key number 04	+\$o
05	key number 05	+\$o
06	key number o6	+\$o
07	key number 07	+\$o
80	key number o8	+\$o
09	key number 09	+\$o
10	key number 10	+\$o
11	key number 11	+\$o
12	key number 12	+\$o
13	key number 13	+\$o
14	key number 14	+\$0
15	key number 15	+ \$0
16	key number 16	+\$0
17	key number 17	+\$o
18	key number 18	+\$0
19	key number 19	+ \$0
20	key number 20	+ \$0
21	key number 21	+ \$0
22	key number 22	+ \$0
23	key number 23	+ \$0
24	key number 24	+\$0
25	key number 25	+\$0
26	key number 26	+\$0
27	key number 27	+\$0

28	key number 28	+\$ 0	75	key number 75	+\$ 0
29	key number 29	+\$ o	76	key number 76	+\$o
30	key number 30	+\$0	77	key number 77	+\$o
31	key number 31	+\$o	78	key number 78	+\$o
32	key number 32	+\$o	79	key number 79	+ \$ 0
33	-	+ \$o	80	key number 80	+\$o
34		+\$0	81	key number 81	+\$o
35		+\$0	82	key number 82	+\$0
36		+\$0	83	key number 83	+\$0
37		+\$0	84	key number 84	+\$0
38		+\$0	85	key number 85	+\$0
39		+\$0	86	key number 86	+\$0
40		+\$0	87		+\$0 +\$0
	•			key number 87	
41		+\$o	88	key number 88	+\$o
42		+\$o	89	key number 89	+\$o
43		+\$o	90	key number 90	+\$o
44		+\$0	91	key number 91	+\$0
45	,	+ \$0	92	key number 92	+\$o
46	46 key number 46	+ \$0	93	key number 93	+\$o
47	key number 47	+\$ 0	94	key number 94	+\$ 0
48	key number 48	+\$ 0	95	key number 95	+\$ 0
49	key number 49	+\$o	96	key number 96	+\$o
50	key number 50	+\$ o	97	key number 97	+\$ o
51	key number 51	+\$ o	98	key number 98	+\$o
52	key number 52	+\$0	99	key number 99	+ \$o
53	key number 53	+\$o			
54		+\$o			
55		+ \$o			
56		+ \$o			
57		+ \$o			
58		+ \$o			
59		+\$0			
60		+\$0			
61	•	+\$0			
62	•	+\$0			
63	•	+\$0			
64		+\$0			
	•	+\$0			
65					
66	•	+\$o			
67		+\$o			
68	•	+\$o			
69		+\$o			
70		+\$o			
71	•	+\$0			
72	•	+ \$o			
73		+ \$0			
74	74 key number 74	+\$o			



Description

This 92"-high, enclosed cabinet attaches to an architectural wall and must stand on the floor. It has a tambour door that slides vertically for access to cabinet materials. The cabinet is available with or without a lock and holds components that support endoscopes and catheters. Attachment hardware is included.


Notes

Order interior components separately:

- Cabinet floor (CV314.)
- Shelf support kit (CV317.)
- Distal tube bracket (CV311.)
- Endoscope bracket (CV310.)

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



24" (610)

Specification Information

keyed alike

Step 1.

CV300.9214

Step 2. Lock

For no lock (*), skip this step.

L keyed differently

Prices for Steps 1-2.

CV300.9214 \$2620

CV300.9214 L \$2781

K \$2781

Step 3. Surface/Door Finish

LU soft white +\$0


C+		
	4. Key Number	
	eyed alike (K)	
01	key number 01	+\$0
02	key number 02	+\$ 0
03	key number 03	+\$ 0
04	key number 04	+\$ 0
05	key number 05	+\$0
06	key number o6	+\$ 0
07	key number 07	+\$0
80	key number o8	+\$0
09	key number 09	+\$ 0
10	key number 10	+\$ 0
11	key number 11	+\$ 0
12	key number 12	+\$ 0
13	key number 13	+\$ 0
14	key number 14	+\$ 0
15	key number 15	+\$ 0
16	key number 16	+\$ 0
17	key number 17	+\$ 0
18	key number 18	+ \$o
19	key number 19	+ \$o
20	key number 20	+ \$o
21	key number 21	+ \$o
22	key number 22	+\$ 0
23	key number 23	+ \$o
24	key number 24	+ \$0
25	key number 25	+\$o
26	key number 26	+\$o
27	key number 27	+\$o
28	key number 28	+\$o

Endoscopy/Catheter Cabinet

continued

Co/Struc® Storage

29	key number 29	+\$o	76	key number 76	+\$o
30	key number 30	+\$o	77	key number 77	+\$0
31	key number 31	+\$o	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$o	80	key number 80	+\$0
34	key number 34	+ \$o	81	key number 81	+\$0
35	key number 35	+\$ 0	82	key number 82	+\$0
36	key number 36	+\$o	83	key number 83	+\$0
37	key number 37	+\$o	84	key number 84	+\$o
38	key number 38	+\$o	85	key number 85	+\$0
39	key number 39	+\$o	86	key number 86	+\$0
40	key number 40	+\$o	87	key number 87	+\$0
		+\$0	88	key number 88	+\$0
41	key number 41				
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$ 0	90	key number 90	+\$0
44	key number 44	+\$o	91	key number 91	+\$0
45	key number 45	+\$ 0	92	key number 92	+\$o
46	key number 46	+\$o	93	key number 93	+\$0
47	key number 47	+\$o	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$ 0	98	key number 98	+\$0
52	key number 52	+\$ 0	99	key number 99	+\$0
53	key number 53	+\$ o			
54	key number 54	+\$ 0			
55	1	+\$o			
	key number 55	140			
56	key number 55 key number 56				
56 57	key number 56	+ \$o			
57	key number 56 key number 57	+\$0 +\$0			
57 58	key number 56 key number 57 key number 58	+\$0 +\$0 +\$0			
57 58 59	key number 56 key number 57 key number 58 key number 59	+\$0 +\$0 +\$0 +\$0			
57 58 59 60	key number 56 key number 57 key number 58 key number 59 key number 60	+\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66 key number 66 key number 66	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66 67	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66 key number 67 key number 67	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66 67 68 69	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66 key number 66 key number 67 key number 68 key number 68	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66 67 68 69 70	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66 key number 67 key number 68 key number 69 key number 70	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66 67 68 69	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66 key number 66 key number 67 key number 68 key number 68	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66 67 68 69 70	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66 key number 67 key number 68 key number 69 key number 70	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66 67 68 69 70 71	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66 key number 67 key number 68 key number 69 key number 70 key number 71 key number 72	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66 67 68 69 70 71 72 73	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66 key number 67 key number 68 key number 69 key number 70 key number 71 key number 72 key number 73	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			
57 58 59 60 61 62 63 64 65 66 67 68 69 70 71	key number 56 key number 57 key number 58 key number 59 key number 60 key number 61 key number 62 key number 63 key number 64 key number 65 key number 66 key number 67 key number 68 key number 69 key number 70 key number 71 key number 72	+\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0 +\$0			



Description


This bracket mounts at any height in a cabinet. It has 5 slots and adjustable cradles for vertical storage of endoscopes. Each slot is $3^3/4^{\prime\prime}$ deep. Mounting hardware is included.


Notes


Order following products separately:


- Cabinet (CV300.)
- Distal tube bracket (CV311.)

Recommended installation height is $5^{1}/_{4}"$ below top of cabinet back panel.









Description


This bracket mounts at any height in a cabinet. It has 5 slots for vertical storage of an endoscope distal tube. Each slot is $^7/_8$ " deep. Mounting hardware is included.

Notes


Order following products separately:

- Cabinet (CV300.)
- Endoscope bracket (CV310.)

Recommended installation height is $44^{1}\!/_{\!4}"$ below top of cabinet back panel.



Spec	ification Information	
Step:	1.	
CV31	1.	\$208
Step:	2. Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$ 0
LU	soft white	+\$ 0
WN	warm grey neutral	+\$0



Description


This kit allows the use of shelves inside the endoscopy/catheter cabinet. It includes a pair of bearers and spacers.

Notes

Order shelves separately:

- CS tray/shelf (CO204)
- CS wire shelf (CO253)

Dimensions



Specification Information


Step 1.

CV317.

\$195

WN

warm grey neutral



Product Information

Description


This component provides storage in an endoscopy/catheter cabinet. It has 2 usable surfaces: 1 surface has raised edges to prevent items from rolling off; the reverse surface is flat for applications where raised edges would interfere with material storage.

Notes

Order shelf support kit (CV317.) separately.

Tray/shelf surface with raised edges can be configured with subcontainers; order subcontainers separately.


Dimensions





Spec	cification Information	
Step	1.	
CO204FF		
Step	2. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+\$0
LU	soft white	+\$o

+\$0





Description


This metal floor encloses the bottom of an endoscopy/catheter cabinet. Attachment hardware is included.

Notes

Order cabinet (CV300.) separately.







Description

This frame hangs on a standard rail or adapter rail, or it is supported by an L cart, process table, counter top, or heavy-duty storage work surface. It has 4 slots at 3" intervals to enclose and support storage components and accessories. The C frame cannot mount under a 24"-wide work surface and it is not recommended for overhead storage.


Notes


Order support components separately:


- Adapter rail (C1610.)
- Counter top (CO555, or CO556)
- L cart (CO212)
- Process table (C4115., C4215., or C4313.)
- Lab work table (C461G., C462G., or C462C.)
- Standard rail (CO282)
- Work table (C4312.)
- Workstation (C4310. or C4311.)

Order storage components separately:


- C tray/shelf (CO203)
- Drawer (CO207, CO208, or CO209)







\$190



Product Information


Description


This surface mounts on a C frame to provide a work surface or storage surface. It includes a gripper package for attaching the top to a C frame. The single or extended counter top fits on 1 C frame. The extended counter top extends 4" past each C frame side.


Notas

Order C frame (CO205) separately.

Dimensions







Specification Information


Step 1.

CO556FF

co

Step 2. Configuration/Width	
555FF single	
556FF extended	
Prices for Steps 1-2.	
CO555FF	\$156

Step 3. Surface Finish			
HF	inner tone light	+\$o	
LT	light tone	+\$o	
LU	soft white	+\$o	
WN	warm grey neutral	+\$o	




Description

This lid provides a top for a C frame. It has raised side edges to keep items from sliding off.

Notes


Use lid when storing L cart below 38"-high process table or cantilevered work surface; existing drawer bearers must be removed.





FRONT

Specif	ication Information	
Step 1.		
CO201FF		\$149
Step 2.	Surface Finish	
HF	inner tone light	+\$o
LT	light tone	+\$0
LU	soft white	+ \$o
WN	warm grey neutral	+\$0



Description


This component provides storage within a C frame, locker, or interior cabinet adapter. It is designed with 2 usable surfaces: 1 surface has raised edges to prevent items from rolling off; the reverse surface is flat for applications where raised edges would interfere with material storage.


Notes

Order support component separately:


- C frame (CO205)
- Interior cabinet adapter (CG591.)
- Locker or locker on wheel base (CO561, CO562, CO563, or CO564)

 Tray/shelf surface with raised edges can be configured with subcontainers or subdivider vane. Order subdivider vane (CO242) or subcontainers separately.







Specification Information			
Step 1.			
CO203FF		\$100	
Step 2. Surface Finish			
HF	inner tone light	+\$o	
LT	light tone	+\$0	
LU	soft white	+\$o	
WN	warm grey neutral	+\$0	





Description

This component organizes IV bags, forms, tape rolls, and other small items within a drawer or C tray/shelf.









Description


This wire shelf fits into the slots of a C frame, locker, or interior cabinet adapter and stores items requiring air circulation. It cannot be used with a subdivider vane or subdivider.

Notes

Order support component separately:

- C frame (CO205)
- Interior cabinet adapter (CG591.)
- Locker or locker on wheel base (CO561, CO562, CO563, or CO564)

Dimensions






FRONT



Step 2. Surface Finish

LU soft white +\$0





Description

This hardware package mounts a C frame under a work surface. 1 gripper package is required to mount each C frame to a work surface. Attachment screws are included.

Notes


C frame cannot mount with grippers under heavy-duty work surface.





FRONT

Spe	cification Information	
Step	1.	
CO266FF		\$60
Step	2. Surface Finish	
HF	inner tone light	+ \$o
LT	light tone	+\$ o
LU	soft white	+\$0
WN	warm grey neutral	+\$0



Description

This mobile caddy fits below a hanging or freestanding work surface. It has full-width pulls, casters, 2 3"-high drawers, 1 6"-high drawer, and 1 9"-high drawer. The drawers can be removed and reconfigured within the caddy or used in other products that accept 20"-wide plastic drawers.


The 9"-high drawer is available with a lock.

Notes

For caddy without drawers, order drawers separately.

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CO69

Step 2. Configuration				
2.20	no drawers			
3.20	drawers			

Step 3. Lock

For no lock (*), skip this step.

For no drawers (2.20), skip this step.

For drawers (3.20)

L 9"-high drawer keyed differently

K 9"-high drawer keyed alike


Prices for Steps 1-3.	\$1362
600)2.20	¥1)02
C0693.20	\$1813
CO693.20 L	\$1914
K	\$1914

Step 4. Surface Finish			
HF	inner tone light	+\$o	
LT	light tone	+\$o	
LU	soft white	+\$o	
WN	warm grey neutral	+\$o	

Step 5	5. Drawer Finish	
For dr	awers (3.20)	
1B	khaki	+\$0
2B	juniper	+\$0
3B	beachglass	+\$0
4B	plum	+\$0
5B	java	+\$0
5E	aquamarine	+\$0
6F	midnight	+\$0
DA	deep blue	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
VR	vivid red	+\$0
WN	warm grey neutral	+\$0
YT	yellow	+\$0
ZP	terra cotta	+\$0

tep	6. Key Number		46	key number 46	+
or 9'	"-high drawer keyed alike (K)		47	key number 47	+
1	key number 01	+\$0	48	key number 48	+
2	key number 02	+ \$0	49	key number 49	+
3	key number 03	+\$0	50	key number 50	+
4	key number 04	+ \$0	51	key number 51	+
15	key number 05	+ \$0	52	key number 52	+
6	key number o6	+ \$o	53	key number 53	+
7	key number 07	+ \$0	54	key number 54	+
8	key number o8	+ \$o	55	key number 55	+
9	key number 09	+ \$o	56	key number 56	+
0	key number 10	+ \$o	57	key number 57	+
1	key number 11	+ \$0	58	key number 58	4
2	key number 12	+ \$o	59	key number 59	+
3	key number 13	+ \$o	60	key number 60	+
4	key number 14	+ \$o	61	key number 61	4
5	key number 15	+ \$o	62	key number 62	=
6	key number 16	+ \$o	63	key number 63	-
7	key number 17	+ \$o	64	key number 64	
3	key number 18	+ \$0	65	key number 65	
•	key number 19	+\$0	66	key number 66	
)	key number 20	+\$0	67	key number 67	
l	key number 21	+\$0	68	key number 68	
	key number 22	+\$0	69	key number 69	
3	key number 23	+\$0	70	key number 70	
ŀ	key number 24	+\$0	71	key number 71	
,	key number 25	+\$o	72	key number 72	
•	key number 26	+\$o	73	key number 73	
,	key number 27	+\$o	74	key number 74	
3	key number 28	+\$o	75	key number 75	
)	key number 29	+\$o	76	key number 76	
)	key number 30	+\$o	77	key number 77	
l	key number 31	+\$o	78	key number 78	
2	key number 32	+\$o	79	key number 79	
3	key number 33	+\$o	80	key number 80	
ŀ	key number 34	+\$o	81	key number 81	
5	key number 35	+\$o	82	key number 82	
5	key number 36	+\$o	83	key number 83	
7	key number 37	+\$o	84	key number 84	
3	key number 38	+\$o	85	key number 85	
)	key number 39	+\$0	86	key number 86	
)	key number 40	+\$0	87	key number 87	
	key number 41	+ \$o	88	key number 88	
2	key number 42	+ \$o	89	key number 89	
3	key number 43	+ \$o	90	key number 90	
4	key number 44	+ \$o	91	key number 91	
5	key number 45	+\$o	92	key number 92	

93	key number 93	+\$o
94	key number 94	+\$o
95	key number 95	+\$o
96	key number 96	+\$o
97	key number 97	+\$o
98	key number 98	+\$o
99	kev number 99	+\$0



Description

This drawer fits into a locker or C frame and has stops to prevent accidental removal. It can also be mounted under a work surface using drawer bearers, ordered separately. The drawer is available with no lock, a keyed lock, or keyless lock.


Notes

To mount drawer under table or work surface, order drawer bearer package (CO228) or lockable drawer bearer (CO396) separately.

Order optional A-size subcontainer (CO231, CO232, CO233, CO234, CO235, or CO369) or B-size subcontainer (CO236, CO238, or CO240) separately.


Drawer with keyless lock option (KP) cannot be used in locker with tambour door.

Keyed lock is located on right side; keyless lock is located on left side. For general information on locks, see Keyed-Alike Information in Appendices.



Caracia	instinu Information		C+	= Var Nirmhar	
	ication Information			5. Key Number	
Step 1.				eyed alike (KA)	Φ.
CO20			01	key number 01	+\$0
Cı	C: /II: 14		02	key number 02	+\$0
	Size/Height		03	key number 03	+\$0
7	a size		04	key number 04	+\$0
8	b size		05	key number 05	+\$0
9	c size		06	key number o6	+\$0
C4	I = -l.		07	key number 07	+\$0
Step 3.			80	key number o8	+\$0
FF 	no lock		09	key number 09	+\$0
FL	keyed differently		10	key number 10	+\$0
KA	keyed alike		11	key number 11	+\$0
KP	keyless lock		12 13	key number 12	+\$0
Dricoc f	or Stone 4.2		13 14	key number 13	+\$o +\$o
C O207	or Steps 1-3.	\$92	15	key number 14 key number 15	+\$0
CU2U/	FL	\$92 \$192	16	key number 16	+\$0
	KA	\$192 \$192	17	key number 17	+\$0
	KP	\$192 \$617	18	key number 18	+\$0
	Kr	\$017	19	key number 19	+\$0
CO208	FF	\$118	20	key number 20	+\$0
0200	FL	\$219	21	key number 21	+\$0
	KA	\$219	22	key number 22	+\$0
	KP	\$643	23	key number 23	+\$0
		4049	24	key number 24	+\$0
CO209	FF	\$149	25	key number 25	+\$0
	FL	\$250	26	key number 26	+\$0
	KA	\$250	27	key number 27	+\$0
	KP	\$674	28	key number 28	+\$0
		1.71	29	key number 29	+\$0
Step 4.	Surface Finish		30	key number 30	+\$0
1B	khaki	+\$o	31	key number 31	+\$0
2B	juniper	+\$0	32	key number 32	+\$0
3B	beachglass	+\$0	33	key number 33	+\$0
4B	plum	+\$0	34	key number 34	+\$0
5B	java	+\$0	35	key number 35	+\$0
5E	aquamarine	+\$0	36	key number 36	+\$0
6F	midnight	+\$0	37	key number 37	+\$0
DA	deep blue	+\$o	38	key number 38	+\$0
HF	inner tone light	+\$o	39	key number 39	+\$0
LT	light tone	+\$o	40	key number 40	+\$0
LU	soft white	+\$0	41	key number 41	+\$0
VR	vivid red	+\$o	42	key number 42	+\$0
	warm grey neutral	+\$o	43	key number 43	+\$0
ΥT	yellow	+\$o	44	key number 44	+\$0
ZP	terra cotta	+\$o	45	key number 45	+\$0

46	key number 46	+ \$o	93	key number 93	+\$ 0
47	key number 47	+ \$o	94	key number 94	+\$ 0
48	key number 48	+ \$0	95	key number 95	+\$0
49	key number 49	+ \$0	96	key number 96	+\$0
50	key number 50	+ \$0	97	key number 97	+\$0
51	key number 51	+ \$0	98	key number 98	+\$0
52	key number 52	+ \$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+ \$0			
56	key number 56	+\$ 0			
57	key number 57	+\$ 0			
58	key number 58	+\$o			
59	key number 59	+\$o			
60	key number 60	+\$o			
61	key number 61	+\$o			
62	key number 62	+\$o			
63	key number 63	+\$o			
64	key number 64	+\$o			
65	key number 65	+\$o			
66	key number 66	+\$o			
67	key number 67	+\$o			
68	key number 68	+\$o			
69 70	key number 69	+\$0 +\$0			
70 71	key number 70	+\$0 +\$0			
71 72	key number 71 key number 72	+\$0 +\$0			
72 73	key number 73	+\$0 +\$0			
74	key number 74	+\$0 +\$0			
75	key number 75	+\$0			
76	key number 76	+\$0			
77	key number 77	+\$0			
78	key number 78	+\$0			
79	key number 79	+\$0			
80	key number 80	+\$o			
81	key number 81	+\$o			
82	key number 82	+\$0			
83	key number 83	+\$o			
84	key number 84	+\$o			
85	key number 85	+\$o			
86	key number 86	+\$o			
87	key number 87	+\$o			
88	key number 88	+\$o			
89	key number 89	+ \$o			
90	key number 90	+ \$o			
91	key number 91	+\$o			
92	key number 92	+\$o			



Description


This drawer has a hinged steel lid that can be locked. The lid can open when the drawer is pulled out to a stop position; the lid cannot be removed.


The drawer is 3", 6", or 9" high.


Notes


For B-size drawer with lockable lid, order A-size subcontainer only. For C-size drawer with lockable lid, order A- or B-size subcontainers only. For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions









Specification Information

Step 1.

CO37

Step 2.	Size/Height	
0	a size	
2	b size	
4	c size	
Step 3.	Lid Lock	
FF	keyed differently	
KA	keyed alike	
Prices f	or Steps 1-3.	
CO370	FF	\$317
	1/ A	ф

Trices for Steps 1-3.	
CO370 FF	\$317
КА	\$317
C0372 FF	\$344
КА	\$344
C0374 FF	\$376
КА	\$376

Step 4	Step 4. Drawer Finish			
1B	khaki	+\$o		
2B	juniper	+\$o		
3B	beachglass	+\$o		
4B	plum	+\$o		
5B	java	+\$o		
5E	aquamarine	+\$o		
6F	midnight	+\$o		
DA	deep blue	+\$o		
HF	inner tone light	+\$o		
LT	light tone	+\$o		
LU	soft white	+\$o		
VR	vivid red	+\$o		
WN	warm grey neutral	+\$o		
YT	yellow	+\$o		
ZP	terra cotta	+\$ o		

Step 5. Lid Finish

soft white


LU

+\$0

03	key number 03	+\$0	50	key number 50	+\$0
04	key number 04	+\$0	51	key number 51	+\$0
05	key number 05	+\$0	52	key number 52	+\$0
06	key number o6	+\$0	53	key number 53	+\$0
07	key number 07	+\$0	54	key number 54	+\$0
80	key number o8	+\$0	55	key number 55	+\$0
09	key number 09	+\$0	56	key number 56	+\$0
10	key number 10	+\$0	57	key number 57	+\$0
11	key number 11	+\$0	58	key number 58	+\$0
12	key number 12	+\$0	59	key number 59	+\$0
13	key number 13	+\$0	60	key number 60	+\$0
14	key number 14	+\$0	61	key number 61	+\$0
15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$ 0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$ 0	96	key number 96	+\$0

Drawer with Lockable Lid continued

97	key number 97	+\$0
98	key number 98	+\$0
99	key number 99	+\$0



Description

This drawer fits into a locker and has stops to prevent accidental removal. It has a non-removable lid that can open when the drawer is pulled out. The drawer is available with a keyed lock or keyless lock.


Notes

Locked drawer with lid cannot be used with a C frame.

Drawer with keyless lock option (KP) cannot be used in locker with tambour door.

Keyed lock is located on right side; keyless lock is located on left side. For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CO5

Step	2. Size/Height
49	a size

50 b size

51 c size

Step 3. Drawer Lock

FL keyed differentlyKA keyed alike

KP keyless lock

or Steps 1-3.	
FL	\$320
KA	\$320
KP	\$749
FL	\$346
KA	\$346
KP	\$774
FL	\$378
KA	\$378
KP	\$806
	FL KA KP FL KA KP

Step 2	4. Drawer Finish	
1B	khaki	+\$0
2B	juniper	+\$0
3B	beachglass	+\$0
4B	plum	+\$0
5B	java	+\$0
5E	aquamarine	+\$0
6F	midnight	+\$0
DA	deep blue	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
VR	vivid red	+\$0
WN	warm grey neutral	+\$0
ΥT	yellow	+\$0
ΖP	terra cotta	+\$0


Step	5. Drawer Key Number	
For k	keyed alike (KA)	
01	key number 01	+\$0

02	key number 02	+ \$0	49	key number 49	+\$o
03	key number 03	+ \$o	50	key number 50	+\$o
04	key number 04	+ \$0	51	key number 51	+\$o
05	key number 05	+ \$0	52	key number 52	+\$o
06	key number o6	+\$o	53	key number 53	+\$o
07	key number 07	+\$o	54	key number 54	+\$o
08	key number o8	+\$0	55	key number 55	+\$o
09	key number 09	+\$0	56	key number 56	+\$o
10	key number 10	+\$0	57	key number 57	+\$o
11	key number 11	+\$0	58	key number 58	+\$o
12	key number 12	+\$0	59	key number 59	+ \$0
13	key number 13	+\$0	60	key number 60	+\$o
14	key number 14	+\$0	61	key number 61	+\$o
15	key number 15	+\$0	62	key number 62	+\$o
16	key number 16	+\$0	63	key number 63	+\$o
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$o	70 77	key number 77	+\$0
31	key number 31	+\$0 +\$0	77 78	key number 78	+\$0 +\$0
32	key number 32	+\$0 +\$0	78 79	key number 79	+\$0 +\$0
33	•	+\$0 +\$0	80		+\$0 +\$0
34	key number 33	+\$0 +\$0	81	key number 80 key number 81	+\$0 +\$0
35	key number 34	+\$0 +\$0	82	•	+\$0 +\$0
36	key number 35			key number 82	
	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$o	85	key number 85	+\$o
39	key number 39	+\$o	86	key number 86	+\$o
40	key number 40	+\$o	87	key number 87	+\$o
41	key number 41	+\$o	88	key number 88	+\$o
42	key number 42	+\$o	89	key number 89	+\$o
43	key number 43	+\$o	90	key number 90	+\$o
44	key number 44	+\$o	91	key number 91	+\$o
45	key number 45	+\$0	92	key number 92	+\$o
46	key number 46	+\$o	93	key number 93	+\$o
47	key number 47	+\$o	94	key number 94	+\$o
48	key number 48	+\$ 0	95	key number 95	+\$o

96	key number 96	+\$
97	key number 97	+\$
98	key number 98	+\$
99	key number 99	+\$0
Step	6. Lid Finish	
LU	soft white	+\$(

Locked Drawer with Lockable Lid

CO371 CO373 CO375



Product Information

Description

This drawer fits into a locker and has stops to prevent accidental removal. It has a lockable, non-removable lid that can open when the drawer is pulled out. The drawer is available with a keyed lock or keyless lock.


Notes

Locked drawer with lockable lid cannot be used with a C frame.


Drawer with keyless lock option (KP) cannot be used in locker with tambour door.


Keyed lock is located on right side; keyless lock is located on left side. For general information on locks, see Keyed-Alike Information in Appendices.


Dimensions



TOP








Speci	fication Information		02	key number
Step 1.			03	key number o
037			04	key number oz
			05	key number o
	. Size/Height		06	key number of
1	a size		07	key number o
3	b size		08	key number o
5	c size		09	key number og
۵.			10	key number 10
	. Drawer/Lid Lock		11	key number 11
FF	keyed differently		12	key number 12
KA	keyed alike		13	key number 13
KP	keyless lock		14	key number 14
	· · ·		15	key number 15
	for Steps 1-3.	.	16	key number 16
CO371		\$410	17	key number 17
	KA	\$410	18	key number 18
	КР	\$840	19	key number 19
			20	key number 20
CO373		\$437	21	key number 21
	KA	\$437	22	key number 22
	KP	\$866	23	key number 23
			24	key number 22
CO375		\$468	25	key number 25
	KA	\$468	26	key number 26
	KP	\$897	27	key number 27
			28	key number 28
	. Drawer Finish		29	key number 29
1B	khaki	+ \$o	30	key number 30
2B	juniper	+ \$o	31	key number 31
3B	beachglass	+ \$o	32	key number 32
4B	plum	+ \$o	33	key number 33
5B	java	+\$0	34	key number 32
5E	aquamarine	+\$ 0	35	key number 35
6F	midnight	+\$0	36	key number 36
DA	deep blue	+\$ 0	37	key number 37
HF	inner tone light	+\$ 0	38	key number 38
LT	light tone	+\$0	39	key number 39
LU	soft white	+ \$o	40	key number 40
VR	vivid red	+ \$o	41	key number 41
WN	warm grey neutral	+ \$o	42	key number 42
ΥT	yellow	+\$0	43	key number 43
ZP	terra cotta	+\$ o	44	key number 42
			45	key number 45
Step 5.	. Drawer Key Number		46	key number 46
	ved alike (KA)		47	key number 47
			48	

49	key number 49	+\$o	96	key number 96	+\$ 0
50	key number 50	+\$o	97	key number 97	+\$o
51	key number 51	+\$o	98	key number 98	+\$o
52	key number 52	+\$o	99	key number 99	+\$o
53	key number 53	+\$o		,	
54	key number 54	+\$0	Sten	6. Lid Finish	
55	key number 55	+\$0	LU	soft white	+\$0
56			LU	Soft writte	+\$0
	key number 56	+\$o	6.	III III III	
57	key number 57	+\$o		7. Lid Key Number	
58	key number 58	+\$o		eyed alike (KA) or keyless lock (KP)	
59	key number 59	+\$o	01	key number 01	+\$o
60	key number 60	+\$o	02	key number 02	+\$o
61	key number 61	+\$o	03	key number 03	+\$o
62	key number 62	+\$o	04	key number 04	+\$o
63	key number 63	+\$o	05	key number 05	+\$o
64	key number 64	+\$o	06	key number o6	+\$o
65	key number 65	+\$0	07	key number 07	+\$0
66	key number 66	+\$0	08	key number o8	+\$0
67	key number 67	+\$o	09	key number 09	+\$0
68	key number 68	+\$0	10	key number 10	+\$0
69	key number 69	+\$o	11	key number 11	+\$o
70	key number 70	+\$o	12	key number 12	+\$o
71	key number 71	+\$o	13	key number 13	+\$o
72	key number 72	+\$o	14	key number 14	+\$o
73	key number 73	+\$o	15	key number 15	+\$o
74	key number 74	+\$o	16	key number 16	+\$o
75	key number 75	+\$o	17	key number 17	+\$o
76	key number 76	+\$o	18	key number 18	+\$o
77	key number 77	+\$0	19	key number 19	+\$0
78	key number 78	+\$0	20	key number 20	+\$0
79	key number 79	+\$o	21	key number 21	+\$0
80	key number 80	+\$o	22	key number 22	+\$o
81	key number 81	+\$o	23	key number 23	+ \$0
82	key number 82	+\$o	24	key number 24	+\$o
83	key number 83	+\$o	25	key number 25	+\$o
84	key number 84	+\$o	26	key number 26	+\$ o
85	key number 85	+\$o	27	key number 27	+\$o
86	key number 86	+\$o	28	key number 28	+\$o
87	key number 87	+\$o	29	key number 29	+\$o
88	key number 88	+\$0	30	key number 30	+\$0
89	key number 89	+\$0	31	key number 31	+\$0
90	key number 90	+\$0	32	key number 32	+\$0
91	key number 91	+\$o	33	key number 33	+\$0
92	key number 92	+\$0	34	key number 34	+\$0
93	key number 93	+\$ 0	35	key number 35	+\$0
94	key number 94	+\$o	36	key number 36	+\$o
95	key number 95	+\$o	37	key number 37	+\$o

38	key number 38	+\$0	85	key number 85	+\$o
39	key number 39	+\$o	86	key number 86	+\$o
40		+\$o		key number 87	+\$o
41		+\$0		key number 88	+\$o
42		+\$0		key number 89	+\$o
43		+\$ 0		key number 90	+\$o
44	key number 44	+\$0	91	key number 91	+\$o
45	key number 45	+\$0	92	key number 92	+\$o
46	key number 46	+\$o	93	key number 93	+\$o
47		+\$o		key number 94	+\$o
48		+\$0		key number 95	+\$0
49		+\$0		key number 96	+\$0
				•	
50		+\$0		key number 97	+\$o
51		+\$ 0		key number 98	+\$o
52	key number 52	+\$0	99	key number 99	+\$o
53	key number 53	+\$ o			
54	key number 54	+\$o			
55		+\$o			
56		+\$0			
57		+\$0			
58		+\$0			
59		+\$ 0			
60	key number 60	+\$0			
61	. key number 61	+\$0			
62	key number 62	+\$o			
63	key number 63	+\$o			
64		+\$o			
65		+\$0			
66		+\$0			
67		+\$0			
68		+\$0			
69		+\$ 0			
70	key number 70	+\$ 0			
71	key number 71	+\$ 0			
72	key number 72	+ \$o			
73		+\$o			
74		+\$0			
75		+\$0			
76		+\$0			
77		+\$0			
78		+\$0			
79		+\$ 0			
80	key number 80	+ \$0			
81	key number 81	+ \$o			
82		+\$o			
83		+\$0			
84		+\$0			
04	r key ilullibel 04	+\$0			




Description

This clear lid fits on a Co/Struc drawer and provides top closure.

Notes

Order drawer (CO207, CO208, or CO209) separately.

Dimensions






Specification Information

Step 1.


CO508.20


\$137




Description


This package is used to mount a nonlocking drawer under a work surface. The drawer bearer has integral stops to prevent drawers from falling when pulled out. Mounting hardware is included.







+\$o





Product Information

Description

This component is used to mount a drawer with lock under a work surface. The lockable drawer bearer has integral stops to prevent drawers from falling when pulled out. Mounting hardware is included.


Dimensions







black




Description

This assembly holds drawers in place on an L cart and indicates that the drawers are unaltered. It must be used with 2 C frames.

Notes


Order breakaway tabs (CT709.) separately.

Spec	Specification Information			
Step 1.				
CO473FF				
<u> </u>	0. 6. 57.1.			
Step :	2. Surface Finish			
LU	soft white	+ \$o		




Description

This clip snaps onto a drawer, tray/shelf, or flipper to hold content labels. The clip accepts pressure-sensitive labels from the front or slide-in labels from the back. Finish is clear. Package contains 12 clips; labels not included.



Specification Information					
Step 1.					
co					
Step 2. Width					
•					
27427 6" wide					
462 19" wide					
Prices for Steps 1-2.					
CO27427	\$27				
CO462	\$27 \$50				





Description

This assembly fits inside a locker or C frame or fits below a work surface with a drawer bearer. It has a wire cassette, cassette bins, and subdividers. The assembly fills 2 slots in a locker or C frame.

Notes


Order optional cassette cover (CO₅66) separately.







Specification Information	
Step 1.	
CO565FF	\$474
Step 2. Surface Finish	
•	. # -
LU soft white	+5




Description

This cover fits on a cassette assembly and provides top closure.


Notes

Order cassette assembly (CO_56_5) separately.

Dimensions




TOP





Step 2.	Surface Finish		
LU	soft white	+5	\$o





Description

This bin fits inside a cassette assembly. Package contains 10 bins, 20 dividers, and 10 label clips.

Notes

Order cassette assembly (CO₅6₅) separately.






Spe	cification Information	
Step	1.	
CO567FF		\$168
Step	2. Surface Finish	
LU	soft white	+\$o

Dispensing Rail

CO365 CO366 CO367 CO368



Product Information

Description

This rail attaches to a panel, a freestanding module, or wall strips to support subcontainers at a 15° angle for gravity-feed dispensing. Package contains 2.

The A-size rail holds A-size subcontainers; the B-size rail holds B-size subcontainers.

Notes

Dispensing rails hold the following number of subcontainers:

A-Size Rail Width—Number of Subcontainers

24"-5 A-1, A-2, A-3 or A-4

24"-2 A-5

24"-10 A-8

48"-10 A-1, A-2, A-3, or A-4

48"-5 A-5

48"-20 A-8

 $\hbox{B-Size Rail Width-Number of Subcontainers}$


24"-5 B-1 or B-3

24"-2 B-5

48"-10 B-1 or B-3

48"-5 B-5

Order optional A-size subcontainer (CO231, CO232, CO233, CO234, CO235, or CO369) or B-size subcontainer (CO236, CO238, or CO240) separately.



Specification Information Step 1. C036 Step 2. Size/Width 5FF a size and 24" wide 6FF a size and 48" wide 7FF b size and 24" wide 8FF b size and 48" wide Prices for Steps 1-2. CO365FF \$336 CO366FF \$442 CO367FF \$373 CO368FF \$447 Step 3. Surface Finish soft white **+\$**o

Subcontainer C0231
C0232
C0233
C0234
C0235
C0236
C0238
C0240
C0369

Product Information

Description

This storage container fits into a drawer or onto a dispensing rail to organize and store small items. Package contains 12.

The A-size subcontainer is used with an A-size drawer (without lid) or A-size dispensing rail. It is $2^3/_{16}$ " high; depth and width dimensions are: Type—Depth—Width

$$A-1-4^9/_{16}"-4^9/_{16}"$$

$$A-2-9^3/_{16}"-4^9/_{16}"$$

A-4-
$$18^{7}/_{16}$$
"- $4^{9}/_{16}$ "

$$A-5-9^3/_{16}"-9^3/_{16}"$$

The B-size subcontainer is used with a B-size drawer (without lid), B-size dispensing rail, or C-size drawer. It is $5^3/_{16}$ " high; depth and width dimensions are:

Type—Depth—Width

B-1-
$$4^9/_{16}$$
"- $4^9/_{16}$ "


$$B-5-9^3/_{16}"-9^3/_{16}"$$


Notes

Order following products separately:

- Drawer (CO207, CO208, or CO209)
- Dispensing rail (CO365, CO366, CO367, or CO368)

Dimensions






Specification Information Step 1. СО Step 2. Size **231FF** a-1 **232FF** a-2 **233FF** a-3 **234FF** a-4 **235FF** a-5 **369FF** a-8 **236FF** b-1 **238FF** b-3 **240FF** b-5 Prices for Steps 1-2. CO231FF \$73 CO232FF \$73 CO233FF \$80 CO234FF \$80 CO235FF \$80 CO369FF \$75 CO236FF \$80 CO238FF \$110 CO240FF \$110 Step 3. Surface Finish LU soft white +\$o

A-Size Drawer Organizer Kit

C0540
C0541
C0542
C0543
C0544



Product Information

Description

This kit fits inside an A-size drawer (without lid) and has 6 combinations of subcontainers and subcontainer dividers.

A-size subcontainer is $2^{3}/_{16}$ " high. Width and depth dimensions are:

Type-Width-Depth

$$A-1-4^9/_{16}"-4^9/_{16}"$$

$$A-2-4^9/_{16}"-9^3/_{16}"$$


$$A-3-4^9/_{16}"-13^{13}/_{16}"$$

$$A-4-4^9/_{16}"-18^7/_{16}"$$

$$A-5-9^3/_{16}"-9^3/_{16}"$$

$$A-8-2^{1/4}$$
" $-13^{13/16}$ "

Dimensions



Specification Information

Step 1.


CO545

CO54

Step 2	Step 2. Components			
OFF	4 a-3 subcontainers, 6 a-3 dividers			
1FF	1 a-3 subcontainer, 6 a-8 subcontainers, 2 a-3 dividers, 7 a-8 dividers			
2FF	2 a-1 subcontainers, 2 a-2 subcontainers, 1 a-3 subcontainer, 2 a-8 subcontainers, 4 a-3 dividers, 3 a-8 dividers			
3FF	3 a-3 subcontainers, 2 a-8 subcontainers, 4 a-3 dividers, 3 a-8 dividers			
4FF	4 a-2 subcontainers, 1 a-4 subcontainer, 8 a-3 dividers			
5FF	2 a-1 subcontainers, 2 a-3 subcontainers, 1 a-5 subcontainer, 4 a-3 dividers			

Prices for Steps 1-2.	
CO540FF	\$67
CO541FF	\$71
CO542FF	\$71
CO543FF	\$67
CO544FF	\$67
CO545FF	\$67

Step 3.	Surface Finish	
LU	soft white	+\$0



Description

This kit fits inside a B-size drawer (without lid) or C-size drawer and has 2 combinations of subcontainers and subcontainer dividers.

The B-size subcontainer is $5^3/_{16}"$ high. Width and depth dimensions are:


Type-Width-Depth


B-1-4⁹/₁₆"-4⁹/₁₆"

 $B-3-4^9/_{16}"-13^{13}/_{16}"$

 $B-5-9^3/_{16}"-9^3/_{16}"$

Dimensions





Specification Information

Step 1.

CO54

Step 2. Components

6FF 4 b-3 subcontainers, 4 b-3 dividers

7FF 2 b-1 subcontainers, 2 b-3 subcontainers, 1 b-5 subcontainer, 4

b-3 dividers

Prices for Steps 1-2.


CO546FF \$82 **CO547FF** \$82

Step 3. Surface Finish

U soft white +\$o

Subcontainer Divider

CO378 CO379 CO380




Product Information

Description

This divider fits into subcontainers to separate and organize stored items. A divider can be adjusted by sliding it to the desired position. The divider is not available for A-5 and B-5 subcontainers. Package contains

The divider for an A-8 subcontainer is $1^5/8$ " high and $2^1/8$ " wide; the divider for an A-1, A-2, A-3, or A-4 subcontainer is $1^5/8$ " high and $4^1/4$ " wide; and the divider for a B-1 or B-3 subcontainer is 4" high and $4^1/4$ " wide.

Dimensions



Specification Information


Step 1.

CO3

Step 2.	Usage
78FF	for a-8 subcontainer
79FF	for a-1, a-2, a-3, or a-4 subcontainers
80FF	for b-1 or b-3 subcontainers

Prices for Steps 1-2.	
CO378FF	\$68
CO379FF	\$68
CO380FF	\$70

Step	3. Surface Finish	
LU	soft white	+\$0



Description

This snap-on handle has a surface for identification labels in a picking station. It also can convert a small or large subcontainer into a drawer when used on the flat surface of a C tray. The handle cannot be used with a subcontainer in an A-size drawer. Package contains 12.


Votes

Handle accepts following label:

- $2^{1}/_{8}$ "-wide handle: $^{3}/_{4}$ " × 1" gummed label
- $4^{1/2}$ "-wide handle: 3/4" × 3" drop-in label

For $4^{1}/2^{"}$ -wide handle, order label (CO459) separately.

Dimensions



Specification Information

Step 1.

CO37

Step 2. Width

6FF 21/8" wide

7FF 4¹/₂" wide

Prices for Steps 1-2.


CO376FF

CO377FF \$44

\$35

Step 3. Surface Finish


LU soft white +\$o



Description

This label slides into a $4^1/2^n$ -wide subcontainer handle. Package contains 100 cards with 10 perforated labels per card (1000 labels total).

Specification Information	
Step 1.	
CO459	\$29




Description

This flag snaps onto subcontainers that are hung from a dispensing rail. It provides a visible surface for identification labels. Package contains 12.

The flag for an A-8 subcontainer is $1^1/4''$ wide; the flag for an A-3 or B-3 subcontainer is $3^1/2''$ wide.

Dimensions



Specification Information

Step 1.

CO38

Step 2. Usage


1FF for a-8 subcontainer2FF for a-3 or b-3 subcontainer

Prices for Steps 1-2.

CO381FF \$100 CO382FF \$111

Step 3. Surface Finish


LU soft white +\$0




Description

This 20"-wide, powder-coated wire rack fits into the slots of a locker to organize catheters. It is available with 3 or 4 independent, pullout slides. Each pullout slide has 4 prongs, 4 side label clips, and 1 front label holder. Each prong holds up to 10 catheters. Finish is soft white.

Dimensions





Specification Information


Step 1.

CO469.20

Step	2. Slides			
03	3 slides			
04	4 slides			

Prices for Steps 1-2.


C0469.20 03 \$732
04 \$815



Description

This height-adjustable laptop cart delivers point-of-care technology and holds a 17" laptop. The surface provides security by locking a laptop into place during use. The cart holds a maximum of 30 pounds and has 4"-high casters; 2 locking and 2 nonlocking. It has an 18" height-adjustment range and is adjusted with a paddle on the front of the surface. The cart includes a swivel mouse tray. Shipped fully assembled. The surface is Cloud Surf laminate; the base finish is polished aluminum.

Dimensions



Specification Information

Step 1.

CO600.

Step 2. Power Option

N no power stripP power strip

Step 3. Bar Code Scanner Holder

N no bar code scanner holderB bar code scanner holder

Step 4. Utility Basket

N no utility basketU utility basket


Step 5. Drawer

N no drawer

Drings for Change . =

D drawer

Prices for Steps 1-5.				
	NN	ND	UN	UD
CO600. N N	\$2662	3043	2864	3245
В	\$2764	3144	2966	3345
P N	\$3000	3379	3202	3581
В	\$3101	3480	3303	3682




Description

This height-adjustable cart delivers point-of-care technology. It holds a flat panel monitor arm and includes a 10"-wide column-mounted CPU holder with straps. The arm meets VESA standards (Video Electronic Standards Association - Mounting Interface Standards). The cart holds a maximum of 30 pounds and has 4"-high casters; 2 locking and 2 nonlocking. It has an 18" height-adjustment range and is adjusted with a paddle on the front of the surface. The cart includes a swivel mouse tray and 2 plates: 75mm and 100mm. Shipped fully assembled.

The surface is Cloud Surf laminate; the base finish is polished aluminum; and the CPU holder is black.

Dimensions



Specification Information

Step 1.

CO601.

Step 2. Power Option

N no power strip

P power strip

Step 3. Bar Code Scanner Holder

no bar code scanner holder

B bar code scanner holder

Step 4. Utility Basket

N no utility basket

U utility basket


Step 5. Drawer

N no drawer

D drawer

Prices for S	Steps 1-5.				
		NN	ND	UN	UD
CO601. N	N	\$2662	3043	2864	3245
	В	\$2764	3144	2966	3345
P	N	\$3000	3379	3202	3581
	В	\$3101	3480	3303	3682





Description

This hospital-grade power cord has a 1.5' coiled cord that extends to 7.5'. Finish is black.

Dimensions


	C*		
Speci	fication	a Intori	memon
Opcc.	cacio.		illudi Oli

Step 1.

CO602.

\$174







Description


This perforated basket attaches to the column of a laptop or flat-panel cart and holds a utility basket. It can be used to hold charts, paper products, and other supplies. Finish is black.


Dimensions












Description

This holder attaches to the surface of a laptop or flat-panel cart. It manages cables and wires and includes reusable cable ties. Finish is black.

Notes

Laptop-cart cable manager has 1 solid holder. Flat-panel-cart cable manager has 2 pieces that attach to each side of the flat-screen arm holder and go down through the work surface.

Dimensions



Specification Information

Step 1.

CO608.


Step 2. Usage

01 for laptop cart

02 for flat panel cart

Prices for Steps 1-2.


C0608. 01 \$162 **02** \$162





Description


This cart is used to maneuver and transport lockers. It has a foot-activated lever that lifts a locker from a rail to a cart or from a cart to a rail. The cart also has a towing handle and a train connector for linking multiple carts. A maximum of 3 carts can be linked when pulled manually; a maximum of 7 carts can be linked when pulled with a tug. Shipped knocked down.

Dimensions









Description

This mobile cart has a platform, 4 rail profiles for hanging frames, and 4 casters (2 swivel and 2 swivel with brake) that are $4^{3}/4''$ hard-rubber zinc.


Notes


To hang drawer or shelf storage components, order C frame (CO205) separately.

Order optional components separately:

• Breakaway lock assembly (CO₄₇₃)

Dimensions





FRONT SIDE

Specification Information

Step 1.

LU

CO212

Step 2. Application

standard

Prices for Steps 1-2.

CO212 FF \$722


Step 3. Surface Finish HF inner tone light LT light tone

warm grey neutral

+\$0 soft white **+**\$o

+\$o

+\$0



Co/Struc® System	
2-Way 90° Corner Connector	page(s) 32
2-Way Module Corner Connector	48
3-Way 90° Corner Connector	33
3-Way Module Corner Connector	49
Access Panel Assembly	58
Adapter Rail	7
ADA Sink Support, Systems Mounted	118
Add-On Shelf, Shelf Unit	140
Add-On Shelf with Lip	152
Adjustable Work Table Base	105
Anchor Support Kit	51
A-Size Drawer Organizer Kit	215
Backsplash Filler	86
Basic Work Table Base	100
Bottle Drawer Divider	130
Breakaway Lock Assembly	206
B-Size Drawer Organizer Kit	216
B-Style Rail-Hanging Shelf	148
B-Style Shelf	149
B-Style Shelf with Sloped Top	151
Cabinet Floor	179
Cable Manager	226
Cassette Assembly	208
Cassette Bin	210
Cassette Cover	209
Ceiling Entry Support Panel Center Closure Panel	27
C Frame	67
C Frame Support	180 83
Clear Lid	203
C Lid	182
Coiled Power Cord	224
Combination Access Panel Assembly	60
Component Adapter Kit	84
Counter Top	181
CS Tray/Shelf	178
C Tray/Shelf	183
C Wire Shelf	185
Diplomat Panel	45
Dispensing Rail	211
Distal Tube Bracket	176
Drawer	190
Drawer Bearer Package	204
Drawer Caddy with Casters	187
Drawer with Lockable Lid	193
Draw Rod, Seismic	13

End Cap for Standard or Seismic Rail	10
End Closure Assembly	65
End Closure Panel	64
Endoscope Bracket	175
Endoscopy/Catheter Cabinet	173
Equipment Rail, Systems Mounted	9
Equipment Rail, Wall Mounted	8
Filler Extrusion	57
Flat Panel Cart	223
Flat Panel Mount	115
Floor Anchor Adapter Kit	17
Floor Anchor Bracket, Action Office® Series 1	15
Floor Anchor Bracket, Action Office® Series 2	16
Floor Anchor Bracket, Co/Struc® Lab Module	14
Floor Entry Support Panel	25
Frame Stiffener, Ethospace® Top Cap	19
Gripper Package	186
Heavy-Duty Corner Work Surface	88
Heavy-Duty Corner Work Surface, Action Office® Series 1 or 2	90° Corner
Connector	90
Heavy-Duty Work Surface	76
nterface Connector	34
nterior Filler	122
nterior Shelf	53
nterior Shelf and Support Kit	55
abel Clip	207
aptop Cart	222
- Cart	228
ight Bracket	143
ockable Drawer Bearer	205, 123
ocked Drawer with Lid	196
ocked Drawer with Lockable Lid	199
ocker	159
ocker Hinged Door	167
ocker Lock Kit	169
ocker Lock Kit, Keyless	171
ocker on Wheel Base	165
ocker on Wheel Base, with Hinged Door	163
ocker on Wheel Base, with Tambour Door	160
ocker Tambour Door	166
ocker, with Hinged Door	157
ocker, with Tambour Door	15/
ow Cup Sink Access Panel Assembly	62
ow Module	39
ow Module Terminal Panel	42
Marker Board	114
Mid-Run Sunnort Double Sided	

Mid-Run Support, Single Sided	36
Module Support Leg	50
Module Wall Start	47
Open Support Panel	21
Panel Connector	35
Panel Stiffener, Action Office® Series 1 or Series 2 Top Cap	18
Panel Wall Start	31
Pass-Through Cover	82
Pass-Through Module	40
Pass-Through Module Terminal Panel	43
Power Strip	113
Process Table, Electric Height Adjustment	98
Pullout Catheter Rack	221
Rail-Hanging Work Surface	81
Rail Hardware	6
Reagent Shelf	99
Riser Frame Divider	111
Riser Frame End Cover	112
Shelf Backsplash	56
Shelf End Cap	54
Shelf Label Clip	153, 141
Shelf Support Kit	177
Shelf Unit	135
Shelf Unit, 90° Corner	139
Shelf Unit Back Panel	142
Shelf Unit, Diagonal Corner	142 137
Shelf Unit, Diagonal Corner Shelf Unit Door	
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors	137
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel	137 144 145 147
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module	137 144 145
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage	137 144 145 147 116
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top	137 144 145 147 116 117
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface	137 144 145 147 116 117 146
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail	137 144 145 147 116 117 146 74
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base	137 144 145 147 116 117 146 74 5
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case	137 144 145 147 116 117 146 74 5
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Bottle Drawer	137 144 145 147 116 117 146 74 5 102
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Bottle Drawer Storage Case Door	137 144 145 147 116 117 146 74 5 102 124 129
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Bottle Drawer Storage Case Door Storage Case Drawer	137 144 145 147 116 117 146 74 5 102 124 129 131
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Storage Case Door Storage Case Drawer Storage Case Leg	137 144 145 147 116 117 146 74 5 102 124 129 131 127
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Bottle Drawer Storage Case Door Storage Case Leg Storage Case Leg Storage Case Pullout Shelf	137 144 145 147 116 117 146 74 5 102 124 129 131 127 132
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Bottle Drawer Storage Case Door Storage Case Door Storage Case Leg Storage Case Pullout Shelf Storage Case Rail, Corner	137 144 145 147 116 117 146 74 5 102 124 129 131 127 132 126
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Bottle Drawer Storage Case Door Storage Case Drawer Storage Case Leg Storage Case Rail, Corner Storage Case Shelf	137 144 145 147 116 117 146 74 5 102 124 129 131 127 132 126 133
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Storage Case Bottle Drawer Storage Case Door Storage Case Drawer Storage Case Leg Storage Case Pullout Shelf Storage Case Shelf Storage Unit	137 144 145 147 116 117 146 74 5 102 124 129 131 127 132 126 133 125
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Storage Case Bottle Drawer Storage Case Door Storage Case Drawer Storage Case Leg Storage Case Pullout Shelf Storage Case Shelf Storage Case Shelf Storage Unit Storage Unit	137 144 145 147 116 117 146 74 5 102 124 129 131 127 132 126 133 125 119
Shelf Unit, Diagonal Corner Shelf Unit Door Shelf Unit Double Doors Shelf Unit Wall Attachment Panel Sink Module Sink Unit for Storage Sloped Shelf Unit Top Standard-Duty Work Surface Standard Rail Standard Work Table Base Storage Case Storage Case Storage Case Bottle Drawer Storage Case Door Storage Case Drawer Storage Case Leg Storage Case Pullout Shelf Storage Case Shelf Storage Unit	137 144 145 147 116 117 146 74 5 102 124 129 131 127 132 126 133 125

Subcontainer	213
Subcontainer Divider	217
Subcontainer Handle	218
Subcontainer Handle Label	219
Subcontainer Label Flag	220
Subdivider Vane	184
Support Panel	23
Support Panel Top Cap	38
Ferminal Panel	29
Гор Shelf	52
TR3 Cart	227
Jtility Basket	225
/ertical Chase Cover	46
Wall Anchor Bracket	12
Nall Rail, Seismic	11
Nall Strip Cover	20
Nire Chase Overlay	73
Nire Chase Rail	68
Nire Chase Rail, 90° Inside Corner	69
Nire Chase Rail Bridge	70
Nire Chase Rail Cover	72
Nire Chase Rail End Cap	71
Norkstation, Electric Height Adjustment	92
Norkstation, Manual Height Adjustment	94
Nork Surface Backsplash	85
Nork Surface Sidesplash, Resin	87
Nork Surface Support	79
Nork Table, Electric Height Adjustment	96
Nork Table Reagent Shelf	110
Nork Table Surface, Rectangular	107

C460G.	Basic Work Table Base	page(s) 100
C461G.	Standard Work Table Base	102
C462C.	Adjustable Work Table Base	105
C462G.		
C1112.	Low Module	39
C1114.	Pass-Through Module	40
C1122.	Low Module Terminal Panel	42
C1122.	Pass-Through Module Terminal Panel	43
C1124.		
C1146.	Diplomat Panel	45
C1150.	Open Support Panel	21
C1151.	Support Panel	23
C1155.	Floor Entry Support Panel	25
C1156.	Ceiling Entry Support Panel	27
C1161.	Terminal Panel	29
C1170.	Equipment Rail, Wall Mounted	8
C1172.	Equipment Rail, Systems Mounted	9
C1211.	Panel Wall Start	31
C1212.	Module Wall Start	47
C1220.	2-Way Module Corner Connector	48
C1221.	2-Way 90° Corner Connector	32
C1230.	3-Way Module Corner Connector	49
C1231.	3-Way 90° Corner Connector	33
C1261.	Interface Connector	34
C1272.	Mid-Run Support, Single Sided	36
C1273.	Mid-Run Support, Double Sided	37
C1274.	Support Panel Top Cap	38
C1275.	Module Support Leg	50
C1276.	Anchor Support Kit	51
C1281.	Panel Connector	35
C1313.	Wire Chase Rail	68
C1314.		
C1315.	Wire Chase Rail Cover	72
C1316.		
C1317.		
C1318.		
C1321.	Wire Chase Rail, 90° Inside Corner	69
C1322.		
C1326.	Wire Chase Rail Bridge	70
C1327.		
C1331.	Wire Chase Rail End Cap	71
C1332.		
C1358.	Wire Chase Overlay	73
C1412.	Access Panel Assembly	58
C1422.	Combination Access Panel Assembly	60
C1432.	Low Cup Sink Access Panel Assembly	62
C1450.	End Closure Panel	64
• · · · · · · · · · · · · · · · · · · ·		

C1451.	End Closure Assembly	65
C1452.		
C1455.	Center Closure Panel	67
C1460.	Vertical Chase Cover	46
C1510.	Shelf End Cap	54
C1510.	Top Shelf	52
C1511.	Interior Shelf	53
C1511.	Shelf End Cap	54
C1515.	Interior Shelf and Support Kit	55
C1520.	Shelf Backsplash	56
C1530.	Filler Extrusion	57
C1531.		
C1610.	Adapter Rail	7
C1611.	Storage Case Rail, Corner	133
C2140.	Work Surface Support	79
C2141.	Heavy-Duty Work Surface	76
C2142.	Standard-Duty Work Surface	74
C2192.	C Frame Support	83
C2193.	Component Adapter Kit	84
C2215.	Heavy-Duty Corner Work Surface	88
C2216.	Heavy-Duty Corner Work Surface, Action Office® Se	eries 1 or 2
	90° Corner Connector	90
C2399.	Pass-Through Cover	82
C2410.	Rail-Hanging Work Surface	81
C2512.	Sink Module	116
C2513.	Sink Unit for Storage	117
C2520.	ADA Sink Support, Systems Mounted	118
C2910.	Work Surface Backsplash	85
C2920.	Backsplash Filler	86
C2930.	Work Surface Sidesplash, Resin	87
C3011.	B-Style Shelf	149
C3019.	Add-On Shelf with Lip	152
C3091.	B-Style Shelf with Sloped Top	151
C4310.	Workstation, Electric Height Adjustment	92
C4311.	Workstation, Manual Height Adjustment	94
C4312.	Work Table, Electric Height Adjustment	96
C4313.	Process Table, Electric Height Adjustment	98
• • • • • • • • • • • • • • • • • •	Riser Frame Divider	111
• • • • • • • • • • • • • • • • • • • •	Riser Frame End Cover	112
• • • • • • • • • • • • • • • • • • • •	Power Strip	113
• • • • • • • • • • • • • • • • • • • •	Marker Board	114
•	Flat Panel Mount	115
• • • • • • • • • • • • • • • • • • • •	Reagent Shelf	99
C4327.		
	Work Table Surface, Rectangular	
• • • • • • • • • • • • • • • • • • • •	Work Table Reagent Shelf	107
• • • • • • • • • • • • • • • • • • • •	Interior Filler	122

C5010.	Storage Unit	119
C5014.	Storage Unit with Sliding Doors	134
C5015.	Storage Unit Shelf	120
C5016.	Storage Unit Door	121
C5017.		
C5020.	Shelf Unit	135
C5021.		-55
C5022.	Shelf Unit, Diagonal Corner	137
C5023.	, 0	,
C5024.	Shelf Unit, 90° Corner	139
C5025.	Add-On Shelf, Shelf Unit	140
C5026.	Shelf Unit Back Panel	142
C5027.	Shelf Unit Door	144
C5027.	Shell dille Boot	-44
C5032.	Shelf Unit Double Doors	145
C5036.	Shelf Unit Wall Attachment Panel	147
C5040.	Storage Case	
C5040.	Storage Case Shelf	124
•		125
C5042.	Storage Case Pullout Shelf	126
C5043.	Storage Case Drawer	127
C5044.	Ct C D-Hl- D	
C5045.	Storage Case Bottle Drawer	129
C5048.	Lockable Drawer Bearer	123
C5049.	Bottle Drawer Divider	130
C5050.	Storage Case Leg	132
C5051.	Light Bracket	143
C5052.	Storage Case Door	131
C5053.		
CO201	C Lid	182
CO203	C Tray/Shelf	183
CO204	CS Tray/Shelf	178
CO205	C Frame	180
CO207	Drawer	190
CO208		-
CO209		
CO212	L Cart	228
	Draw Rod, Seismic	13
•••••		166
	Drawer Bearer Package	• • • • • • • • • • • • • • • • • • • •
•	Subcontainer	204
CO231	Subcontaillei	213
CO232		
CO233		
CO234		
CO235		
CO236		
CO238		
CO2/10		

CO242	Subdivider Vane	184
CO252	C Wire Shelf	185
CO266	Gripper Package	186
CO274	Label Clip	207
CO277	Locker Lock Kit	169
CO279	Locker Lock Kit, Keyless	171
CO282	Standard Rail	5
CO284	Rail Hardware	6
CO285		
CO290	Locker Hinged Door	167
CO338	End Cap for Standard or Seismic Rail	10
CO342	TR3 Cart	227
CO365	Dispensing Rail	211
CO366		
CO367		
CO368		
CO369	Subcontainer	213
CO370	Drawer with Lockable Lid	193
CO371	Locked Drawer with Lockable Lid	199
CO372	Drawer with Lockable Lid	193
CO373	Locked Drawer with Lockable Lid	199
CO374	Drawer with Lockable Lid	193
CO375	Locked Drawer with Lockable Lid	199
CO376	Subcontainer Handle	218
CO377		
CO378	Subcontainer Divider	217
CO379		
CO380		
CO381	Subcontainer Label Flag	220
CO382		
CO396	Lockable Drawer Bearer	205
CO397.	Sloped Shelf Unit Top	146
CO398	Shelf Label Clip	153, 141
CO459	Subcontainer Handle Label	219
CO462	Label Clip	207
CO469.	Pullout Catheter Rack	221
CO473	Breakaway Lock Assembly	206
CO479	Wall Rail, Seismic	11
	Floor Anchor Bracket, Action Office® Series 1	15
CO485	Floor Anchor Bracket, Action Office® Series 2	16
CO487.	Wall Strip Cover	20
CO508.	Clear Lid	203
CO521.	B-Style Rail-Hanging Shelf	148
CO540	A-Size Drawer Organizer Kit	215
CO541		
CO542		
CO543		

CO544		
CO545		
CO546	B-Size Drawer Organizer Kit	216
CO547		
CO549	Locked Drawer with Lid	196
CO550		
CO551		
CO552	Floor Anchor Adapter Kit	17
CO553	Panel Stiffener, Action Office® Series 1 or Series 2 Top Cap	18
CO554	Frame Stiffener, Ethospace® Top Cap	19
CO555	Counter Top	18:
CO556		
C0561	Locker, with Tambour Door	154
CO562	Locker	159
C0563	Locker on Wheel Base	165
CO564	Locker on Wheel Base, with Tambour Door	160
C0565	Cassette Assembly	208
C0566	Cassette Cover	209
C0567	Cassette Bin	210
C0568	Floor Anchor Bracket, Co/Struc® Lab Module	1/
CO569.	Wall Anchor Bracket	12
CO570	Locker, with Hinged Door	157
CO574	Locker on Wheel Base, with Hinged Door	16
CO600.	Laptop Cart	222
CO601.	Flat Panel Cart	223
CO602.	Coiled Power Cord	22/
CO605.	Utility Basket	225
CO608.	Cable Manager	226
CO692.	Drawer Caddy with Casters	187
CO693.		
CV300.	Endoscopy/Catheter Cabinet	173
CV310.	Endoscope Bracket	175
CV311.	Distal Tube Bracket	176
CV314.	Cabinet Floor	179
CV317.	Shelf Support Kit	177



20-Day or Less/Assigned Lead-Time Order Information — Healthcare

20-Day or Less/Assigned Lead-Time Programs

20-day or less and assigned lead-time products are designed for established Herman Miller dealers with pre-approved credit.

Products are distributed from multiple manufacturing facilities across the United States. All products and options not designated by an Assigned Lead-Time icon $\boxed{\mathbb{A}}$ will ship in 20 business days or less after being acknowledged by Herman Miller. Products and options designated by an $\boxed{\mathbb{A}}$ are on the assigned lead-time program and will ship in 20 days, less than 20 days, or more than 20 days from order acknowledgement.

Ordering Procedure

Please place orders through Order Manager. If this is not available to you, place orders through mail or fax. Verbal purchase orders will not be accepted.

Order Entry fax number: (616) 654 3085.

For more information, contact your Customer Care representative at: $(866)\ 854\ 3048\ \text{EXT}\ 3400$

Shipments and Delivery

Per Herman Miller Terms and Conditions.

Changes and Cancellation

Per Herman Miller Terms and Conditions.

HermanMiller Healthcare April 2014



Packaging/Keyed-Alike Information — Co/Struc® System

Packaging

Packaging varies according to individual product. A surcharge may be added for special packaging requests. Contact Customer Care for specific information.

Keyed-Alike Locks

All components with locks may be specified keyed alike.

Keyed-alike Co/Struc® components with a "CO" or "CV" prefix will be shipped with lock plugs installed. Mobile Technology Support products (CO600, CO601, and CO606) cannot be ordered keyed alike.

Keyed-alike Co/Struc components with a "C" prefix will be shipped with lock chassis installed. Lock plugs will be shipped separately, tagged for product and location as specified. Installation of lock plugs will be on-site according to instructions provided. Herman Miller, Inc., will not be responsible for lock plug installation.

To order keyed-alike locks for "C" prefix components, specify keyed alike. Then place a separate order for plugs and locks:

- List plug part number and key number (232092-XXX); select a key number (or numbers) between 226 and 324.
- 2. List quantity of locks per key number.
- Indicate this is a no charge order.Keyed-alike locks are field installed.

A Medeco master key can be used to open any lock on Co/Struc components with a "CO" or "CV" prefix, except for Mobile Technology Support products CO600, CO601, and CO606. Customers interested in ordering master keys should contact Customer Care with a written request. To remove lock plugs from existing components, loosen the hex nut inside the product and remove the plug. These locks are available in key numbers 01 through 99. To remove lock plugs from locker with electronic lock, remove the cover from the e-lock mounting plate (be careful not to separate battery wires from unit). Next, loosen the hex nut inside the cover and remove the plug.

A Hudson master key can be used to open any lock on Co/Struc components with a "C" prefix. To remove lock plugs from existing components, a removal key or lock change tool is required. Customers interested in ordering master keys, removal keys, or lock change tools should contact Customer Care with a written request. Locks for "C" prefix components are available in key numbers 226 through 324 but are ordered as key numbers 01 through 99. Use the following chart to cross-reference the number selection.

"C" lock	Actual "C" lock
order numbers	plug numbers
1 2 3	226
2	227
3	228
4	229
5	230
6	231
7	232
8	233
9	234
10	235
11	236
12	237
13	238
	239
14	
15	240
16	241
17	242
18	243
19	244
20	245
21 22	246
22	247
23	248
24	249
25	250
26	251
27	252
28	253
29	254
30	255
30 31	256
32	257
33	258
34	259
35	260
36	261
37	262
38	263
39	264
40	265
41	266
42	267
43	268
44	269
45	270
46	271
47	272
48	273
49	274
50	275

"C" lock	Actual "C" lock
order numbers	plug numbers
51	276
52	277
53	278
54	
	279
	280
56	281
57	282
58	283
59	284
60	285
61	286
62	287
63	288
64	
	289
65	290
66	291
67	292
68	293
69	294
70	295
71	296
72	297
73	298
74	299
75	300
76	301
77	302
78	303
79	304
80	305
81	306
82	307
83	308
84	309
85	310
86	311
87	312
88	313
89	314
90	315
91	316
92	317
93	318
94	319
95	
	320
96	321
97	322
98	323
99	324

HermanMiller Healthcare April 2014

For more information about our products and services or to see a list of dealers, please visit us at HermanMiller.com or call (800) 851 1196.

© 2018 Herman Miller, Inc., Zeeland, Michigan

[®] HermanMiller, ♠, Action Office, Aeron, Aside, Avive, Canvas Office Landscape, Caper, Celle, Co/Struc, Eames, Embody, Equa 2, Ergon 3, Ethospace, Exclave, Fine-Tune, FlexFront, Formcoat, Ground Cloth, Herman Miller Healthcare, Layout Studio, Maharam, Meridian, Mirra, Passage, Pellicle, Perspectives, PostureFit, Resolve, Sayl, Setu, Thrive, Tu, Y-Tower, and Z-Axis are among the registered trademarks of Herman Miller, Inc., and its owned subsidiaries.

TM AireWeave, Bubbletack, Cellular Suspension, Compass, Connect, DOT, Durawrap, Everywhere, Flex-Edge, Hopsak 2, Intent, Iota, Keyless Entry, Latitude, Lyris 2, Mora, Motia, Multiscrim, Nelson, Renew, Stackable, Swoop, TriFlex, Twist, and Verus are among the trademarks of Herman Miller, Inc., and its owned subsidiaries.

 ${}^{\rm S}{\rm Chemsurf}^{\otimes}$ is a registered trademark of Wilsonart International.

 ${\sf Corian}^{\circledast}$ is a registered trademark of E.I. du Pont de Nemours and Company Corporation.

FLEXNET™ is a trademark of Milliken & Company.

Floor Saver $\mathsf{Glide}^{\texttt{@}}$ is a registered trademark of Carpin Manufacturing, Inc.

 $\mathsf{Geiger}^{\circledast}$ is a registered trademark of Geiger International.

 ${\rm Microsoft}^{\otimes}$ and ${\rm Natural}^{\otimes}$ Keyboard are registered trademarks of Microsoft Corporation.

Place Textiles $^{\!\scriptscriptstyle{(\!0)}}$ is a registered trademark of Place Textiles.